
<b>Instructions for use</b>	Expiratory valve set for HAMILTON-C2/C3/C6 .....	3
<b>Gebrauchsanweisung</b>	Expirationsventil für HAMILTON-C2/C3/C6 .....	5
<b>Instructions d'utilisation</b>	Système de valves expiratoires pour HAMILTON-C2/C3/C6 .....	8
<b>Instrucciones de uso</b>	Conjunto de válvula espiratoria para HAMILTON-C2/C3/C6 .....	11
<b>Istruzioni per l'uso</b>	Set valvola espiratoria per HAMILTON-C2/C3/C6 .....	14
<b>Instruções de uso</b>	Conjunto de válvula expiratória para HAMILTON-C2/C3/C6 .....	17
<b>使用说明</b>	HAMILTON-C2/C3/C6 呼气阀套件 .....	20
<b>使用説明書</b>	HAMILTON-C2/C3/C6用呼気弁セット .....	22
<b>Инструкции по эксплуатации</b>	Комплект клапана выдоха для аппаратов ИВЛ HAMILTON-C2/C3/C6 .....	24
<b>Kullanım talimatları</b>	HAMILTON-C2/C3/C6 için ekspiratuvar valf seti .....	27


**REF** 160176

**CE** **EAC**

**HAMILTON**  
**MEDICAL**

- 1** Metal plate toward ventilator  
 Metallplättchen weist zum Beatmungsgerät  
 Plaque métallique orientée vers le ventilateur  
 Placa de metal hacia el respirador  
 Piastra metallica verso il ventilatore  
 Placa de metal (na direção do respirador)  
 金属片朝向呼吸机  
 人工呼吸器側 金属プレート  
 Металлическая пластина, направленная  
 к аппарату ИВЛ  
 Ventilátore bakan metal plaka

- 2** Expiratory valve membrane  
 Expirationsventilmembran  
 Membrane de valve expiratoire  
 Membrana de la válvula espiratoria  
 Membrana della valvola espiratoria  
 Membrana da válvula expiratória  
 呼气阀膜  
 呼吸弁ダイヤフラム  
 Мембрана клапана выдоха  
 Ekspiratuvar valf membrani
- 3** Expiratory valve housing  
 Expirationsventilkörper  
 Boîtier de valve expiratoire  
 Carcasa de la válvula espiratoria  
 Calotta della valvola espiratoria  
 Invólucro da válvula expiratória  
 呼气阀外壳  
 呼吸弁ハウジング  
 Корпус клапана выдоха  
 Ekspiratuvar valf yuvasi
- 4** Expiratory valve locking ring  
 Sicherungsring des Expirationsventils  
 Bague de verrouillage de la valve expiratoire  
 Anillo de bloqueo de la válvula espiratoria  
 Anello di bloccaggio della valvola espiratoria  
 Anel de bloqueio da válvula expiratória  
 呼气阀锁环  
 呼吸弁ロックリング  
 Уплотнительное кольцо клапана выдоха  
 Ekspiratuvar valf kilit halkası

## Conventions

### ⚠ WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

### ⚠ CAUTION

*Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.*

### NOTICE

Emphasizes information of particular importance.

## 1 Overview

The expiratory valve set (PN 160176) can be used with HAMILTON- C2/C3/C6 ventilators.

### 1.1 Intended use

The Hamilton Medical single-use expiratory valve set for the HAMILTON-C2/C3/C6 is an accessory for the exhalation gas path. Its main function is to control expiratory gas flow and maintain user-set pressure levels. The Hamilton Medical single-use expiratory valve set is not intended for use with equipment other than Hamilton Medical ventilators.

The single-use expiratory valve set is intended for use by qualified, trained personnel.

## 2 Safety information

### ⚠ WARNING

- Reprocessing of Hamilton Medical single-use products can affect the product properties and may cause injury to the patient. For example, a change to the surface structure during reprocessing

may lead to a change in the tear strength or cause actual cracking. Furthermore, an altered surface structure may result in, for example, a microbial aggregation of spores, allergens, and pyrogens, or cause an increase in the number of particles released as a result of chemical changes in the material properties.

- Hamilton Medical does not assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user. Any liability on the part of Hamilton Medical for any damages resulting from the use of reprocessed or reused single-use parts, or use of any non-Hamilton Medical accessories and consumables, is excluded.
- Follow hospital infection control procedures for all tests and use.
- Remove the protective cap before use.
- Discard the product if it fails two subsequent preoperational checks or if there is any sign of damage.

## ⚠ CAUTION

- (USA only) Federal law restricts this device to sale by or on the order of a physician.
- Visually inspect the expiratory valve housing and membrane. Discard the expiratory valve set if there is any sign of damage.
- Always carry out all preoperational checks before connecting the patient to prevent possible patient injury. Make sure another source of ventilatory support is available.
- Ensure the expiratory valve set is in its intended position in the ventilator.
- Handle used products as contaminated goods according to local laws and regulations or hospital internal procedures.

## NOTICE

- Before use, read the appropriate ventilator *Operator's Manual*.
- Replace the product in accordance with hospital infection control procedures, or depending on the patient's secretions and nebulization of medication. Consider using a filter at the expiratory port during nebulization.
- Any incident with the device leading to serious patient injury, death, or a potential threat to public health must

be reported to the manufacturer and the relevant authorities.


- Ensure the membrane of the expiratory valve set is in its intended position in the ventilator.

### 3 Connecting the expiratory valve set

The expiratory valve set must be mounted between the ventilator and the expiratory limb.

#### To install the expiratory valve

- 1 Remove the safety cover.
- 2 Ensure the membrane is properly aligned with the expiratory valve housing and the metal plate faces up (1).
- 3 Position the expiratory valve set in the expiratory port (2) and twist the locking ring clockwise until it locks into place.

### 4 Calibration and preoperational check

Carry out the tightness test. Follow the instructions provided in the ventilator *Operator's Manual*.

### 5 Delivery

The components are delivered clean and ready for clinical use.

### 6 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

### 7 Technical specifications

#### Environmental requirements and storage conditions

Operating temperature	5°C to 40°C (41°F to 104°F)
Storage temperature & humidity	-20°C to 60°C (-4°F to 160°F) 10% to 95% relative humidity, noncondensing
Conical connectors	OD22/ID15 according to ISO 5356-1:2015

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Specifications are subject to change without notice. For all proprietary trademarks and third-party trademarks used by Hamilton Medical AG, see [www.hamilton-medical.com/trademarks](http://www.hamilton-medical.com/trademarks). © 2021 Hamilton Medical AG. All rights reserved.
Hamilton Medical AG
medin Medical Innovations GmbH