
Instructions for use	Expiratory valve set, single use	3
Gebrauchsanweisung	Expirationsventil, für den Einmalgebrauch	6
Instructions d'utilisation	Système de valves expiratoires à usage unique	9
Instrucciones de uso	Conjunto de válvula espiratoria desechable	12
Istruzioni per l'uso	Set valvola espiratoria, monouso	15
Instruções de uso	Conjunto de válvula expiratória, uso único	18
使用说明	呼气阀套件，一次性使用	21
使用説明書	単回使用呼气弁セット	23
Инструкции по эксплуатации	Комплект клапана выдоха, одноразовый	25
Kullanım talimatları	Ekspiratuvar valf seti, tek kullanımlık	28


REF 950158

CE **EAC**

HAMILTON
MEDICAL

- 1** Metal plate toward ventilator
 Metallplättchen weist zum Beatmungsgerät
 Plaque métallique orientée vers le ventilateur
 Placa de metal hacia el respirador
 Plastra metallica verso il ventilatore
 Placa de metal (na direção do respirador)
 金属片朝向呼吸机
 人工呼吸器側 金属プレート
 Металлическая пластина, направленная к аппарату ИВЛ
 Ventilatore bakan metal plaka
- 2** Expiratory valve membrane (blue)
 Expirationsventilmembran (blau)
 Membrane de valve expiratoire (bleue)
 Membrana de la válvula espiratoria (azul)
 Membrana della valvola espiratoria (blu)
 Membrana da válvula expiratória (azul)
 呼气阀膜 (蓝色)
 呼气弁ダイヤフラム (青)
 Мембрана клапана выдоха (голубая)
 Ekspiratuvar valf membranı (mavi)
- 3** Expiratory valve housing (transparent)
 Expirationsventilkörper (transparent)
 Boîtier de valve expiratoire (transparent)
 Carcasa de la válvula espiratoria (transparente)
 Calotta della valvola espiratoria (trasparente)
 Invólucro da válvula expiratória (transparente)
 呼气阀外壳 (透明)
 呼气弁ハウジング (透明)
 Корпус клапана выдоха (бесцветный)
 Ekspiratuvar valf yuvası (şeffaf)

Conventions

⚠ WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

⚠ CAUTION

Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.

NOTICE

Emphasizes information of particular importance.

1 Overview

The expiratory valve set (PN 950158) can be used with HAMILTON-G5/S1 and GALILEO ventilators.

1.1 Intended use

The Hamilton Medical single-use expiratory valve set for HAMILTON-G5/S1 and GALILEO is an accessory for the exhalation gas path. Its main function is to control expiratory gas flow and maintain user-set pressure levels. The Hamilton Medical single-use expiratory valve set is not intended for use with equipment other than Hamilton Medical ventilators.

The single-use expiratory valve set is intended for use by qualified, trained personnel.

2 Safety information

⚠ WARNING

- Reprocessing of Hamilton Medical single-use products can affect the product properties and may cause injury to the patient. For example, a change to the surface structure during reprocessing may lead to a change in the tear strength or cause actual cracking. Furthermore, an altered surface structure may result in, for example, a microbial aggregation of spores, allergens, and pyrogens, or cause an increase in the

number of particles released as a result of chemical changes in the material properties.

- Hamilton Medical does not assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user. Any liability on the part of Hamilton Medical for any damages resulting from the use of reprocessed or reused single-use parts, or use of any non-Hamilton Medical accessories and consumables, is excluded.
- Follow hospital infection control procedures for all tests and use.
- Remove the protective cap before use.
- Discard the product if it fails two subsequent preoperational checks or if there is any sign of damage.
- Adding attachments or other parts/assemblies to a breathing system can change the pressure gradient across the ventilator, which can adversely affect ventilator performance.

⚠ CAUTION


- (USA only) Federal law restricts this device to sale by or on the order of a physician.
- Visually inspect the expiratory valve housing and membrane. Discard the expiratory valve set if there is any sign of damage.
- Always carry out all preoperational checks before connecting the patient to prevent possible patient injury. Make sure another source of ventilatory support is available.
- Ensure the expiratory valve set is in its intended position in the ventilator.
- Handle used products as contaminated goods according to local laws and regulations or hospital internal procedures.

NOTICE

- Before use, read the appropriate ventilator *Operator's Manual*.
- Replace the product in accordance with hospital infection control procedures, or depending on the patient's secretions and nebulization of medication. Consider using a filter at the expiratory port during nebulization.

- Any incident with the device leading to serious patient injury, death, or a potential threat to public health must be reported to the manufacturer and the relevant authorities.
- Any bacteria filter, HMEF, or additional accessories in the expiratory limb may substantially increase flow resistance and impair ventilation.
- When adding parts to the Hamilton Medical breathing circuit configurations, do not exceed the inspiratory and expiratory resistance values of the ventilator breathing system as specified in the ventilator *Operator's Manual*.
- Ensure the expiratory valve set is in its intended position in the ventilator.
- Ensure the membrane is aligned with the expiratory valve body and the metal faces up.

3 Installing the expiratory valve set

- 1 When holding the expiratory valve set, the metal plate must face up.
- 2 Position the expiratory valve set and turn until it locks into place.

4 Calibration and preoperational check

Carry out the tightness test. Follow the instructions provided in the ventilator *Operator's Manual*.

5 Delivery

The components are delivered clean, and ready for clinical use.

6 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

7 Environmental and storage conditions

Environment	Specification
Operating temperature	10°C to 40°C (50°F to 104°F)
Storage temperature	-10°C to 60°C (14°F to 140°F)
Operating humidity	5% to 85% relative humidity, noncondensing
Conical connectors	OD22/ID15 according to ISO 5356-1:2015

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Specifications are subject to change without notice. For all proprietary trademarks (®) and third-party trademarks (§) used by Hamilton Medical AG, see www.hamilton-medical.com/trademarks. © 2021 Hamilton Medical AG. All rights reserved.
Hamilton Medical AG
medin Medical Innovations GmbH