Volumetric Capnography User Guide

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10107270/00 | 2020-11-30

**ERC**

**HAMILTON**  
**MEDICAL**  
Intelligent Ventilation since 1983

---

English .....	3
Deutsch .....	11
Français .....	19
Español .....	27
Italiano.....	35
Português.....	43
日本語 .....	51
中文 .....	59
Русский .....	67

# User guide

## Volumetric capnography

2020-11-30

10107270/00

---

© 2020 Hamilton Medical AG. All rights reserved. Printed in Switzerland.

No part of this publication may be reproduced, stored in a database or retrieval system, or transmitted in any form or by any means, electronic, mechanical, or by photocopying, recording, or otherwise, without prior written permission of Hamilton Medical AG.

This document may be revised, replaced, or made obsolete by other documents by Hamilton Medical AG at any time and without notice. Ensure that you have the most current applicable version of this document; if in doubt, contact the technical support department of Hamilton Medical AG, Switzerland. While the information set forth herein is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this document shall limit or restrict in any way Hamilton Medical AG's right to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, Hamilton Medical AG has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

The equipment must be operated, serviced, or upgraded only by trained professionals. Hamilton Medical AG's sole responsibility with respect to the equipment and its use is as stated in the limited warranty provided in the device *Operator's Manual*.

Hamilton Medical AG shall not be liable for any loss, cost, expense, inconvenience, or damage that may arise out of misuse of the product, or if non-Hamilton Medical AG parts were used when replacing parts, or if serial numbers were amended, deleted, or removed.

If returning parts to Hamilton Medical AG, be sure to use the standard Hamilton Medical returned goods authorization (RGA) procedure. Disposal of parts shall follow all local, state, and federal regulation with respect to environmental protection.

For all proprietary as well as third-party trademarks used by Hamilton Medical AG, see [www.hamilton-medical.com/trademarks](http://www.hamilton-medical.com/trademarks). Product and/or company names marked with a <sup>§</sup> symbol may be the trademarks and/or registered trademarks of their respective owners, including but not limited to Aerogen<sup>§</sup>, Nihon Kohden<sup>§</sup>, Masimo<sup>§</sup>, Capnostat<sup>§</sup>, Bluetooth<sup>§</sup>, and Wi-Fi<sup>§</sup>.

#### **Manufacturer**

Hamilton Medical AG  
Via Crusch 8, CH-7402 Bonaduz  
Switzerland  
Phone: (+41) 58 610 10 20  
Fax: (+41) 58 610 00 20  
[info@hamilton-medical.com](mailto:info@hamilton-medical.com)  
[www.hamilton-medical.com](http://www.hamilton-medical.com)

1	Introduction .....	6
2	Three phases of the volumetric capnogram .....	6
3	CO <sub>2</sub> elimination (V'CO <sub>2</sub> ).....	7
4	End-tidal CO <sub>2</sub> (PetCO <sub>2</sub> and FetCO <sub>2</sub> ) .....	8
5	Airway dead space (VDaw).....	8
6	Alveolar minute ventilation (V'alv) .....	8
7	Capnogram shape (slopeCO <sub>2</sub> ).....	9
8	Formulas .....	9
9	References .....	10

## 1 Introduction

Volumetric capnography is an important tool to evaluate the quality and quantity of ventilation. CO<sub>2</sub> monitoring data is helpful for the assessment of a patient's airway integrity and to ensure proper endotracheal tube placement, among other applications.

Volumetric capnography measurements that are monitored by the ventilator include:

- CO<sub>2</sub> elimination ( $V'CO_2$ )
- End-tidal partial pressure CO<sub>2</sub> ( $P_{et}CO_2$ )
- End-tidal fraction of CO<sub>2</sub> ( $F_{et}CO_2$ )
- Airway dead space ( $V_{Daw}$ )
- Alveolar minute ventilation ( $V'_{alv}$ )
- Capnogram shape ( $slope_{CO_2}$ )

These parameters are described in further detail in this guide.

Additional information about volumetric capnography can be found on [www.hamilton-medical.com](http://www.hamilton-medical.com) at MyHamilton.

## 2 Three phases of the volumetric capnogram


The alveolar concentration of carbon dioxide (CO<sub>2</sub>) is the result of metabolism, cardiac output, lung perfusion, and ventilation. Change in the concentration of CO<sub>2</sub> reflects perturbations in any or a combination of these factors.

Volumetric capnography provides continuous monitoring of CO<sub>2</sub> production, ventilation/perfusion (V/Q) status, and airway patency, as well as function of the ventilator breathing circuit itself.

Expired gas receives CO<sub>2</sub> from three sequential compartments of the airways, forming three recognizable phases on the expired capnogram. A single breath curve in volumetric capnography exhibits these three characteristic phases of changing gas mixtures – they refer to the airway region in which they originate:

- *Phase I:* Represents CO<sub>2</sub>-free gas from the apparatus and anatomical dead space
- *Phase II:* Represents gas coming from regions that are in the transition between anatomic and alveolar gas compartments; receives CO<sub>2</sub> from units with different perfusion and ventilatory rates
- *Phase III:* (Alveolar plateau) represents the CO<sub>2</sub>-rich gas from the alveoli

Figure 1. Phases of the volumetric capnogram

Using features from each phase, physiologic measurements can be calculated.

### 3 CO<sub>2</sub> elimination (V'CO<sub>2</sub>)

CO<sub>2</sub> elimination (V'CO<sub>2</sub>) is the net exhaled volume of CO<sub>2</sub> in milliliters per minute (ml/min). V'CO<sub>2</sub> measurements permit assessment of metabolic rate (for example, V'CO<sub>2</sub> is high with sepsis and fever) and treatment progress.

V'CO<sub>2</sub> is obtained by adding together VCO<sub>2</sub> measurements over several breaths and dividing the sum by the total time in minutes (Noe 1963). Steady-state conditions are essential to interpret the V'CO<sub>2</sub> values (Brandt 1999). V'CO<sub>2</sub> represents CO<sub>2</sub> elimination but not necessarily CO<sub>2</sub> production. *Normal* values for V'CO<sub>2</sub> are provided in Table 1.

To convert a time-based capnogram into a volumetric capnogram, CO<sub>2</sub> must be combined with flow.

By combining the FetCO<sub>2</sub>/Time graph (Figure 2) with the Flow/Time graph (Figure 3), you can derive the volume of CO<sub>2</sub> exhaled in one breath (Figure 4).

The area under the expiratory waveform minus the area under the inspiratory waveform is the net transfer of CO<sub>2</sub> out of the lungs per breath, or VCO<sub>2</sub>.

Figure 2. Typical capnogram of patient on pressure-controlled ventilation, showing FCO<sub>2</sub>/time<sup>1</sup>


Figure 3. Typical spirogram, patient on pressure-controlled ventilation<sup>2</sup>



Figure 4. Combination of capnogram and spirogram (fractional end-tidal CO<sub>2</sub> concentration/volume)<sup>3</sup>

<sup>1</sup> Inspiration starts at time 0; exhalation starts at approximately 2.75 seconds. Note that inspiratory gas initially contains CO<sub>2</sub> (rebreathing) from artificial deadspace.

<sup>2</sup> The flow to the patient (inspiration) is negative, while the flow from the patient (exhalation) is positive. The expiratory flow waveform is an exponential decay curve. Note that in spontaneously breathing subjects, the flow waveforms may be different shapes.

<sup>3</sup> ViCO<sub>2</sub> is the volume of inspired CO<sub>2</sub>, while VeCO<sub>2</sub> is the volume of exhaled CO<sub>2</sub>. The net elimination of CO<sub>2</sub> is VeCO<sub>2</sub> – ViCO<sub>2</sub>. ViCO<sub>2</sub>, a negative volume indicating rebreathed CO<sub>2</sub>, is normally omitted.

## 4 End-tidal CO<sub>2</sub> (PetCO<sub>2</sub> and FetCO<sub>2</sub>)

End-tidal CO<sub>2</sub> is the maximum partial pressure of CO<sub>2</sub> exhaled during a tidal breath, measured just before the start of inspiration. It represents the final portion of air that was involved in the exchange of gases in the alveolar area, thus providing a reliable index of CO<sub>2</sub> partial pressure in the arterial blood under certain circumstances.

End-tidal CO<sub>2</sub> is either measured as a partial pressure (PetCO<sub>2</sub>) or as a fractional concentration of CO<sub>2</sub> in dry gas (FetCO<sub>2</sub>).

Normal values for PetCO<sub>2</sub> and FetCO<sub>2</sub> are provided in Table 1.


## 5 Airway dead space (VDaw)

The series, anatomical, or airway dead space measurement, VDaw, provides an effective in-vivo measure of volume lost in the conducting airways; this represents the CO<sub>2</sub>-free gas from the apparatus and the anatomical dead space.

A relative increase in dead space points to a rise in respiratory insufficiency and can be regarded as an indicator of the patient's current condition.

Using the Aström method<sup>4</sup>, VDaw can be calculated as the smallest measurable dead space. The calculation<sup>5</sup> comprises a number of computational steps that take the slope of the alveolar plateau into account.

Figure 5. VDaw

Normal values for VDaw are provided in Table 1.

## 6 Alveolar minute ventilation (V'alv)

Alveolar minute ventilation (V'alv) permits assessment of actual alveolar ventilation (as opposed to minute ventilation). V'alv is defined as the difference between minute ventilation and V'Daw (airway dead space).

Not all gas that enters the alveoli participates in gas exchange. Some gas ends up in non- or under-perfused lung spaces. To measure the efficiency of alveolar ventilation, PaCO<sub>2</sub> must be determined from an arterial blood gas sample. The ratio of mixed-to-ideal alveolar partial pressure is a measure of alveolar efficiency (Severinghaus 1957).

Normal values for V'alv are provided in Table 1.

<sup>4</sup> Aström 2000; see Section 9

<sup>5</sup> Wolff 1989, Aström 2000; see Section 9

## 7 Capnogram shape (slopeCO2)

The slope of the alveolar plateau (slopeCO2) is defined by the PetCO2 capnogram shape and can indicate the volume/flow status of the lungs.

Normal values for slopeCO2 are provided in Table 1.

Table 1. Examples of normal or expected values in mechanically ventilated patients<sup>6</sup>

Description	Unit <sup>7</sup>	Normal	Reference
VDaw	ml BTPS	2.2 ml/kg IBW	Radford 1954
slopeCO2	%CO2/l	31324 * Vt - 1.535	Aström 2000
V'CO2 <sup>8</sup>	ml/min STPD	2.6 to 2.9 ml/min/kg	Weissmann 1986, Wolff 1986
FetCO2 <sup>9</sup>	%	5.1 to 6.1%	Wolff 1986
PetCO2	mmHg	32 to 42 mmHg	www.oem.respironics.com
V'alv	l/min	0.052 to 0.070 l/min/kg	Kiiski, Takala 1994 <sup>8</sup>

## 8 Formulas

### Alveolar tidal ventilation (Vtalv)

$$V_{talv} = V_t - V_{Daw}$$

### Alveolar minute ventilation (V'alv)

$$V'_{alv} = f * V_{talv}$$

### Volume of CO2 eliminated in one breath (VCO2)

$$V_{CO2} = V_{eCO2} - V_{iCO2}$$

### Fractional concentration of CO2 in exhaled gas (FetCO2)

$$F_{etCO2} = V'_{CO2} / V_{MinVol}$$

### Partial pressure of CO2 in exhaled gas (PetCO2)

$$P_{etCO2} = F_{etCO2} * (P_b - P_{H2O})$$

<sup>6</sup> These values are for illustration purposes and do not replace physician-directed treatment.

<sup>7</sup> Bulk gas volumes, such as minute ventilation and tidal volume, are usually measured in BTPS. Specific gas volumes are expressed in STPD. Conversion factors can be found in physics textbooks.

<sup>8</sup>  $V'_{CO2} = V'_{alv} * F_{etCO2}$

<sup>9</sup>  $F_{etCO2} = P_{etCO2} / (P_b - P_{H2O})$

## 9 References

- Aström E, Niklason L, Drefeldt B, Bajc M, Jonson B. Partitioning of dead space – a method and reference values in the awake human. *Eur Respir J*. 2000 Oct; 16(4):659-664.1
- Brandi LS, Bertolini R, Santini L, Cavani S. Effects of ventilator resetting on indirect calorimetry measurement in the critically ill surgical patient. *Crit Care Med*. 1999 Mar; 27(3):531-539.
- Kiiski R, Takala J, Eissa NT. Measurement of alveolar ventilation and changes in deadspace by indirect calorimetry during mechanical ventilation: a laboratory and clinical validation. *Crit Care Med*. 1991 Oct; 19(10):1303-1309.
- Noe FE. Computer analysis of curves from an infrared CO<sub>2</sub> analyzer and screen-type airflow meter. *J Appl Physiol* 1963; 18:149-157.
- Radford EP. Ventilation standards for use in artificial respiration. *N Engl J Med* 1954; 251:877-883.
- Severinghaus JW, Stupfel M. Alveolar dead space as an index of distribution of blood flow in pulmonary capillaries. *J Appl Physiol* 1957; 10:335-348.
- Weissman C, Kemper M, Elwyn DH, Askanazi J, Hyman AI, Kinney JM. The energy expenditure of the mechanically ventilated critically ill patient. An analysis. *Chest*. 1986 Feb; 89(2):254-259.
- Wolff G, Brunner JX, Weibel W, et al. Anatomical and series dead space volume: concept and measurement in clinical practice. *Appl Cardiopul Pathophysiol* 1989; 2:299-307.
- Wolff G, Brunner JX, Grädel E. Gas exchange during mechanical ventilation and spontaneous breathing. Intermittent mandatory ventilation after open heart surgery. *Chest*. 1986 Jul; 90(1):11-17.

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$$\text{PetCO2} = \text{FeCO2} * (\text{Pb} - \text{PH2O})$$

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<sup>26</sup>  $\text{V'CO2} = \text{V'alv} * \text{FetCO2}$

<sup>27</sup>  $\text{FetCO2} = \text{PetCO2}/(\text{Pb} - \text{PH2O})$

## 9 Références

- Aström E, Niklason L, Drefeldt B, Bajc M, Jonson B. Partitioning of dead space – a method and reference values in the awake human. *Eur Respir J.* 2000 Oct; 16(4):659-664.1
- Brandi LS, Bertolini R, Santini L, Cavani S. Effects of ventilator resetting on indirect calorimetry measurement in the critically ill surgical patient. *Crit Care Med.* 1999 Mar; 27(3):531-539.
- Kiiski R, Takala J, Eissa NT. Measurement of alveolar ventilation and changes in deadspace by indirect calorimetry during mechanical ventilation: a laboratory and clinical validation. *Crit Care Med.* 1991 Oct; 19(10):1303-1309.
- Noe FE. Computer analysis of curves from an infrared CO<sub>2</sub> analyzer and screen-type airflow meter. *J Appl Physiol* 1963; 18:149-157.
- Radford EP. Ventilation standards for use in artificial respiration. *N Engl J Med* 1954; 251:877-883.
- Severinghaus JW, Stupfel M. Alveolar dead space as an index of distribution of blood flow in pulmonary capillaries. *J Appl Physiol* 1957; 10:335-348.
- Weissman C, Kemper M, Elwyn DH, Askanazi J, Hyman AI, Kinney JM. The energy expenditure of the mechanically ventilated critically ill patient. An analysis. *Chest.* 1986 Feb; 89(2):254-259.
- Wolff G, Brunner JX, Weibel W, et al. Anatomical and series dead space volume: concept and measurement in clinical practice. *Appl Cardiopul Pathophysiol* 1989; 2:299-307.
- Wolff G, Brunner JX, Grädel E. Gas exchange during mechanical ventilation and spontaneous breathing. Intermittent mandatory ventilation after open heart surgery. *Chest.* 1986 Jul; 90(1):11-17.

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<sup>35</sup>  $V'_{CO2} = V'_{alv} * F_{etCO2}$

<sup>36</sup>  $F_{etCO2} = P_{etCO2} / (P_b - P_{H2O})$

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- Aström E, Niklason L, Drefeldt B, Bajc M, Jonson B. Partitioning of dead space – a method and reference values in the awake human. *Eur Respir J.* 2000 Oct; 16(4):659-664.1
- Brandi LS, Bertolini R, Santini L, Cavani S. Effects of ventilator resetting on indirect calorimetry measurement in the critically ill surgical patient. *Crit Care Med.* 1999 Mar; 27(3):531-539.
- Kiiski R, Takala J, Eissa NT. Measurement of alveolar ventilation and changes in deadspace by indirect calorimetry during mechanical ventilation: a laboratory and clinical validation. *Crit Care Med.* 1991 Oct;19(10):1303-1309.
- Noe FE. Computer analysis of curves from an infrared CO<sub>2</sub> analyzer and screen-type airflow meter. *J Appl Physiol* 1963; 18:149-157.
- Radford EP. Ventilation standards for use in artificial respiration. *N Engl J Med* 1954; 251:877-883.
- Severinghaus JW, Stupfel M. Alveolar dead space as an index of distribution of blood flow in pulmonary capillaries. *J Appl Physiol* 1957; 10:335-348.
- Weissman C, Kemper M, Elwyn DH, Askanazi J, Hyman AI, Kinney JM. The energy expenditure of the mechanically ventilated critically ill patient. An analysis. *Chest.* 1986 Feb; 89(2):254-259.
- Wolff G, Brunner JX, Weibel W, et al. Anatomical and series dead space volume: concept and measurement in clinical practice. *Appl Cardiopul Pathophysiol* 1989; 2:299-307.
- Wolff G, Brunner JX, Grädel E. Gas exchange during mechanical ventilation and spontaneous breathing. Intermittent mandatory ventilation after open heart surgery. *Chest.* 1986 Jul; 90(1):11-17.

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<sup>62</sup>  $V' CO2 = Valv/min * FetCO2$

<sup>63</sup>  $FetCO2 = ETCO2 / (Pb - PH2O)$

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- Aström E, Niklason L, Drefeldt B, Bajc M, Jonson B. Partitioning of dead space – a method and reference values in the awake human. *Eur Respir J.* 2000 Oct; 16(4):659-664.1
- Brandi LS, Bertolini R, Santini L, Cavani S. Effects of ventilator resetting on indirect calorimetry measurement in the critically ill surgical patient. *Crit Care Med.* 1999 Mar; 27(3):531-539.
- Kiiski R, Takala J, Eissa NT. Measurement of alveolar ventilation and changes in deadspace by indirect calorimetry during mechanical ventilation: a laboratory and clinical validation. *Critical Care Medicine* 1991 Oct; 19(10):1303-1309.
- Noe FE. Computer analysis of curves from an infrared CO2 analyzer and screen-type airflow meter. *J Appl Physiol* 1963; 18:149-157.
- Radford EP. Ventilation standards for use in artificial respiration. *N Engl J Med* 1954; 251:877-883.
- Severinghaus JW, Stupfel M. Alveolar dead space as an index of distribution of blood flow in pulmonary capillaries. *J Appl Physiol* 1957; 10:335-348.
- Weissman C, Kemper M, Elwyn DH, Askanazi J, Hyman AI, Kinney JM. The energy expenditure of the mechanically ventilated critically ill patient. An analysis. *CHEST Journal* 1986 Feb; 89(2):254-259.
- Wolff G, Brunner JX, Weibel W, et al. Anatomical and series dead space volume: concept and measurement in clinical practice. *Appl Cardiopul Pathophysiol* 1989; 2:299-307.
- Wolff G, Brunner JX, Grädel E. Gas exchange during mechanical ventilation and spontaneous breathing. Intermittent mandatory ventilation after open heart surgery. *CHEST Journal* 1986 Jul; 90(1):11-17.

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