
<b>Instructions for use</b>	Breathing circuit set, dual limb.....	7
<b>Gebrauchsanweisung</b>	Beatmungsschlauchset, mit zwei Schenkeln.....	12
<b>Instructions d'utilisation</b>	Kit de circuit respiratoire à deux branches .....	18
<b>Instrucciones de uso</b>	Equipo respiratorio de dos ramas .....	24
<b>Istruzioni per l'uso</b>	Set circuito paziente, branca doppia .....	29
<b>Instruções de uso</b>	Kit de respiração, alça dupla .....	34
<b>使用说明</b>	呼吸装置, 双肢.....	39
<b>使用説明書</b>	デュアルリム呼吸回路セット.....	43
<b>Инструкции по эксплуатации</b>	Комплект дыхательного контура с двумя патрубками.....	48
<b>Kullanım Talimatları</b>	Solunum devresi seti, çift parçalı .....	54

**REF** 260180, 260182, 260170, 260169

CE 0197

ERC


**HAMILTON**  
**MEDICAL**

<p>1 To patient Zum Patienten Port inspiratoire Hacia el paciente Al paziente Para o paciente 至病人 吸気側接続口 Порт «к пациенту» Hastaya</p>	<p>4 Inspiratory limb Inspirationsschenkel Branche inspiratoire Rama inspiratoria Branca inspiratoria Alça inspiratória 吸気肢 吸気回路 Патрубок вдоха Inspiratuvar parça</p>	<p>7 T-piece or Y-piece, with pressure line T-Stück oder Y-Stück mit Druckmessleitung Pièce en T ou pièce en Y, avec ligne de pression Pieza en T o pieza en Y, con línea de presión Raccordo a T o raccordo a Y con linea di pressione Peça em "T" ou peça em "Y", com linha de pressão 带有耐压管路的 T 形管或 Y 形管 T ピースまたは Y ピース、圧力ライン付き T- или У-образный коннектор с линией контроля давления T parçası veya Y parçası, basınç hattı ile birlikte</p>
<p>2 From patient Vom Patienten Port expiratoire Desde el paciente Dal paziente Do paciente 自病人 呼気側接続口 Порт «от пациента» Hastadan</p>	<p>5 Expiratory limb Expirationsschenkel Branche expiratoire Rama espiratoria Branca espiratoria Alça expiratória 呼気肢 呼気回路 Патрубок выдоха Ekspiratuvar parça</p>	<p>8 Patient interface (mask or nasal prongs)* Patientenanschluss (Maske oder Nasenprongs)* Interface patient (masque ou embouts nasaux)* Interfaz del paciente (máscara o adaptadores nasales)* Interfaccia paziente (maschera o nasal prongs)* Interface do paciente (máscara ou cateteres nasais)* 病人接口 (面罩或鼻塞)* 患者インターフェース (マスクまたは鼻カニューレ)* Интерфейс пациента (маска или наконечники назальной канюли)* Hasta arayüzü (maske veya nasal prong)*</p>
<p>3 Inspiratory filter* Inspirationfilter* Filtre inspiratoire* Filtro inspiratorio* Filtro inspiratorio* Filtro inspiratório* 吸気过滤器* 吸気フィルタ* Фильтр вдоха* Inspiratuvar filtre*</p>	<p>6 Pressure line (blue) Druckmessleitung (blau) Ligne de pression (bleue) Línea de presión (azul) Linea di pressione (azzurro) Linha de pressão (azul) 耐压管路 (蓝色) 圧力ライン (青) Линия контроля давления (голубая) Basınç hattı (mavi)</p>	<p>9 Elbow adapter Ellenbogenadapter Adaptateur coudé Adaptador acodado Adattatore a gomito Adaptador de cotovelo 弯接头 エルボアダプタ Угловой адаптер Dirsek adaptörü</p>


\* Not included | Nicht im Lieferumfang enthalten | Non inclus | No incluido | Non incluso | Não incluído | 不含 | 含まれていない | Не входит в комплект | Dahil değildir



Use with pressure line (nCPAP, nCPAP-PC) (C1/T1/MR1 only) | Verwendung mit Druckmessleitung (Modi nCPAP, nCPAP-PC) (nur C1/T1/MR1) | Utilisation avec ligne de pression (nCPAP, nCPAP-PC) (C1/T1/MR1 uniquement) | Se usa con la línea de presión (nCPAP, nCPAP-PC) (solo C1/T1/MR1) | Utilizzare con la linea di pressione (nCPAP, nCPAP-PC) (solo C1/T1/MR1) | Uso com linha de pressão (nCPAP, nCPAP-PC) (somente C1/T1/MR1) | 使用耐圧管路 (nCPAP、nCPAP-PC) (仅限 C1/T1/MR1) | 圧力ラインと使用 (nCPAP、nCPAP-PC) (C1/T1/MR1のみ) | Использование с линией контроля давления (режимы nCPAP и nCPAP-PC) (только для C1/T1/MR1) | Basınç hattı ile birlikte kullanın (nCPAP, nCPAP-PC) (sadece C1/T1/MR1)

<p><b>1</b> To patient Zum Patienten Port inspiratoire Hacia el paciente Al paziente Para o paciente <b>至病人</b> 吸気側接続口 Порт «к пациенту» Hastaya</p>	<p><b>5</b> Expiratory limb Expirationsschenkel Branche expiratoire Rama espiratoria Branca espiratoria Alça expiratória <b>呼气肢</b> 呼气回路 <b>Патрубок выдоха</b> Ekspiratuvar parça</p>	<p><b>8</b> HMEF/HME* HMEF/HME* FECH/ECH* HMEF/HME* HMEF/HME* HMEF/HME* HMEF/HME* 人工鼻 (HMEF/HME) * <b>Теплооблагодобменник (с фильтром или без него)*</b> HMEF/HME*</p>
<p><b>2</b> From patient Vom Patienten Port expiratoire Desde el paciente Dal paziente Do paciente <b>自病人</b> 呼気側接続口 <b>Порт «от пациента»</b> Hastadan</p>	<p><b>6</b> Swivel Y-piece Schwenkbares Y-Stück Pièce en Y articulée Pieza en Y inclinable Raccordo a Y girevole Peça em "Y" giratória <b>旋转 Y 形管</b> 旋回 Y ピース <b>У-образный коннектор</b> с шарнирным креплением Muylu Y parçası</p>	<p><b>9</b> nCPAP-PS flow sensor position nCPAP-PS: Position des Flow-Sensors Position du capteur de débit nCPAP-AI Posición del sensor de flujo nCPAP-PS Posizione del sensore di flusso per nCPAP-PS Posição do sensor fluxo nCPAP-PS <b>nCPAP-PS 流量传感器位置</b> nCPAP-PS でのフローセンサの位置 <b>Расположение датчика потока nCPAP-PS</b> nCPAP-PS akış sensörü konumu</p>
<p><b>3</b> Inspiratory filter* Inspirationsfilter* Filtre inspiratoire* Filtro inspiratorio* Filtro inspiratorio* Filtro inspiratório* <b>吸気过滤器 *</b> 吸気フィルタ * <b>Фильтр вдоха*</b> Inspiratuvar filtre*</p>	<p><b>7</b> Flow sensor Flow-Sensor Capteur de débit Sensor de flujo Sensore di flusso Sensor fluxo <b>流量传感器</b> フローセンサ <b>Датчик потока</b> Akış sensörü</p>	<p><b>10</b> Elbow adapter Ellenbogenadapter Adaptateur coudé Adaptador acodado Adattatore a gomito Adaptador de cotovelo <b>弯接头</b> エルボアダプタ <b>Угловой адаптер</b> Dirsek adaptörü</p>
<p><b>4</b> Inspiratory limb Inspirationsschenkel Branche inspiratoire Rama inspiratoria Branca inspiratoria</p>	<p>Alça inspiratória <b>吸気肢</b> 吸気回路 <b>Патрубок вдоха</b> Inspiratuvar parça</p>	<p>* Not included   Nicht im Lieferumfang enthalten   Non inclus   No incluido   Non incluso   Não incluído   <b>不含</b>   含まれていない   <b>Не входит в комплект</b>   Dahil değildir</p>

Use with flow sensor (all modes except nCPAP, nCPAP-PC) | Verwendung mit Flow-Sensor (alle Modi außer nCPAP, nCPAP-PC) | Utilisation avec capteur de débit (tous modes excepté nCPAP, nCPAP-PC) | Se usa con sensor de flujo (en todos los modos, excepto en nCPAP y nCPAP-PC) | Utilizzare con il sensore di flusso (tutte le modalità tranne nCPAP, nCPAP-PC) | Uso com sensor fluxo (todos os modos, exceto nCPAP, nCPAP-PC) | 使用流量传感器 (除 nCPAP 和 nCPAP-PC 外的所有其它模式) | フローセンサと使用 (nCPAP、nCPAP-PCを除くすべてのモード) | Использование с датчиком потока (все режимы, кроме nCPAP и nCPAP-PC) | Akış sensörü ile birlikte kullanın (nCPAP, nCPAP-PC dışındaki tüm modlar)

- 1 Expiratory valve membrane  
 Expirationsventilmembran  
 Membrane de valve expiratoire  
 Membrana de la válvula espiratoria  
 Membrana della valvola espiratoria  
 Membrana da válvula expiratória  
 呼气阀膜  
 呼吸弁ダイヤフラム  
 Мембрана клапана выдоха  
 Ekspiratuvar valf membrani
- 2 Metal plate toward ventilator  
 Metallplättchen weist zum Beatmungsgerät  
 Plaque métallique orientée vers le ventilateur  
 Placa de metal hacia el respirador  
 Piastra metallica verso il ventilatore  
 Placa de metal (na direção do respirador)  
 金属片朝向呼吸机  
 人工呼吸器側 金属プレート  
 Металлическая пластина, направленная  
 к аппарату ИВЛ  
 Ventilatore bakan metal plaka
- 3 Expiratory valve housing  
 Expirationsventilkörper  
 Boîtier de valve expiratoire  
 Carcasa de la válvula espiratoria  
 Calotta della valvola espiratoria  
 Invólucro da válvula expiratória  
 呼气阀外壳  
 呼吸弁ハウジング  
 Корпус клапана выдоха  
 Ekspiratuvar valf yuvasi

## REF

260180	Breathing circuit (1.5 m), Y-piece, flow sensor, pressure line
260182	Breathing circuit (3.0 m), Y-piece, flow sensor, pressure line
260170	Breathing circuit (1.5 m), Y-piece, flow sensor, expiratory valve, pressure line
260169	Breathing circuit (3.0 m), Y-piece, flow sensor, expiratory valve, pressure line

## Conventions

### ⚠ WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

### ⚠ CAUTION

*Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.*

### NOTICE

Emphasizes information of particular importance.

## 1 Overview

The dual limb breathing circuit set (PN 260180) can be used with all Hamilton Medical ventilators. For an overview of the breathing circuit set and its components, see pages 2 through 6.

The dual limb breathing circuit set (PN 260182) can be used with HAMILTON-C1/C2/C3/C6/T1/MR1/G5/S1 ventilators.

The dual limb breathing circuit sets (PN 260170, 260169) can be used with HAMILTON-C1/T1/MR1 ventilators.

### 1.1 Intended use

#### PN 260180

The Hamilton Medical neonatal breathing circuit set with PN 260180 includes a 1.5 m breathing tube, Y-piece, flow sensor, pressure line, connectors, T-piece, and elbow adapters.

This breathing circuit set is intended to be used to connect the HAMILTON-C1/C2/C3/C6/T1/MR1/G5/S1 with the patient interface during invasive or noninvasive ventilation of neonatal patients.

#### PN 260182

The Hamilton Medical neonatal breathing circuit set with PN 260182 includes a 3.0 m breathing tube, Y-piece, flow sensor, pressure line, connectors, T-piece, and elbow adapters.

This breathing circuit set is intended to be used to connect the HAMILTON-C1/C2/C3/C6/T1/MR1/G5/S1 with the patient interface during invasive or noninvasive ventilation of neonatal patients.

The Hamilton Medical neonatal flow sensor is a device to measure patient air flow to be used with the neonatal patient group of Hamilton Medical ventilators.

The Hamilton Medical pressure line is intended to be used to measure the airway pressure in the breathing circuit set with the HAMILTON-C1/T1/MR1 ventilator. The pressure line is for use with nCPAP or nCPAP-PC modes and is connected to the airway measurement port of the ventilator. In all other modes, the Hamilton Medical flow sensor is used to measure patient airflow on Hamilton Medical ventilators.

The breathing circuit set is a medical device intended for use by qualified, trained healthcare professionals.

#### PN 260170, 260169

The Hamilton Medical neonatal ventilator breathing circuit sets with PN 260170 and PN 260169, comprise breathing circuits that include a Y-piece, flow sensor, pressure line, connectors, T-piece, elbow adapters, and expiratory valve set.

Each of these breathing sets is intended to be used to connect the HAMILTON-C1, HAMILTON-T1, and HAMILTON-MR1 ventilators to the patient endotracheal tube or respiratory mask during neonatal ventilation.

The Hamilton Medical neonatal flow sensor is a device to measure patient air flow to be used with the neonatal patient group of Hamilton Medical ventilators.

The Hamilton Medical pressure line is intended to be used to measure the airway pressure in the breathing circuit set with the HAMILTON-C1/T1/MR1 ventilator. The pressure line is for use with nCPAP or nCPAP-PC modes and is connected to the airway measurement port of the ventilator. In all other modes, the Hamilton Medical flow sensor is used to measure patient airflow on Hamilton Medical ventilators.

The Hamilton Medical single-use neonatal expiratory valve set for HAMILTON-C1/T1/MR1 is an accessory for the exhalation gas path.

Its main function is to control expiratory gas flow and maintain user-set pressure levels. The Hamilton Medical single-use expiratory valve set is not intended for use with equipment other than Hamilton Medical ventilators.

The breathing set is a medical device intended for use by qualified, trained healthcare professionals.

## 2 Safety information

### WARNING

- Reprocessing of Hamilton Medical single-use products can affect the product properties and may cause injury to the patient. For example, a change to the surface structure during reprocessing may lead to a change in the tear strength or cause actual cracking. Furthermore, an altered surface structure may result in, for example, a microbial aggregation of spores, allergens, and pyrogens, or cause an increase in the number of particles released as a result of chemical changes in the material properties.
  - Hamilton Medical does *not* assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user. Any liability on the part of Hamilton Medical for any damages resulting from the use of reprocessed or reused single-use parts, or use of any non-Hamilton Medical accessories and consumables, is excluded.
- Follow hospital infection control procedures for all tests and use.
  - Remove the safety cover(s) before use.
  - Discard the product if it fails two subsequent preoperational checks or if there is any sign of damage.
  - To prevent patient or ventilator contamination, *always* use a bacteria filter or HMEF between the patient and the inspiratory port. If no bacteria filter is used, the exhaled gas can contaminate the ventilator.
  - Adding attachments or other parts/assemblies to a breathing system can change the pressure gradient across the ventilator, which can adversely affect ventilator performance.
  - Connect only a single breathing circuit at a time. Do *not* connect any circuit as an extension to another one. Doing so can seriously harm or kill the patient.
  - Route the breathing circuit without tension and without any kinks from the ventilator or humidifier to the patient, and protect against unintended influences.

## CAUTION

- (USA only) Federal law restricts this device to sale by or on the order of a physician.
- Handle used products as contaminated goods according to local laws and regulations or hospital internal procedures.
- Always carry out all preoperational checks before connecting the patient to prevent possible patient injury. Make sure another source of ventilatory support is available.

## NOTICE

- Before use, read the appropriate ventilator *Operator's Manual*.
- Replace the product in accordance with hospital infection control procedures, or depending on the patient's secretions and nebulization of medication. Consider using a filter at the expiratory port during nebulization.
- Any incident with the device leading to serious patient injury, death, or a potential threat to public health must be reported to the manufacturer and the relevant authorities.
- Any bacteria filter, HMEF, or additional accessories in the expiratory limb may substantially increase flow resistance and impair ventilation.

- When adding parts to the Hamilton Medical breathing circuit configurations, do *not* exceed the inspiratory and expiratory resistance values of the ventilator breathing system as specified in the ventilator *Operator's Manual*.
- Certain parts of the breathing circuit may be visible on the MR image. Keep the breathing circuit out of the scanned area as much as possible.
- Ensure the membrane of the expiratory valve set is in its intended position in the ventilator.
- Use of the calibration adapter is required for flow sensor calibration with all ventilators (HAMILTON-C1/T1/MR1 ventilators only with software version 2.1.x and higher.).
- (HAMILTON-C1/T1/MR1 only). A neonatal flow sensor is required with breathing circuits used for all ventilation modes except nCPAP and nCPAP-PC.
- (HAMILTON-C1/T1/MR1 only.) Ensure you select the correct expiratory valve set (adult/pediatric or neonatal) for your patient. If the expiratory valve set does not match the selected patient group on the ventilator, the **Wrong expiratory valve alarm** is generated.

- (HAMILTON-C1/T1/MR1 only.) When using the nCPAP or nCPAP-PC modes, remove the flow sensor and use the pressure line with the breathing circuit.

## 3 Connecting and calibrating components

Setting up a neonatal breathing circuit set comprises the following steps:


To ...	See ...
Connect the neonatal expiratory valve set	Section 3.1
Connect the breathing circuit	Section 3.2
Connect the neonatal flow sensor or pressure line	Section 3.3

### 3.1 Connecting the neonatal expiratory valve set (HAMILTON-C1/T1/MR1 only)

The neonatal expiratory valve set must be mounted between the ventilator and the expiratory limb.

#### To install the expiratory valve

- 1 Remove the safety cover.
- 2 Ensure the membrane is properly aligned with the expiratory valve housing and the metal plate faces up (A).
- 3 Position the expiratory valve set in the expiratory port (B) and twist the locking ring clockwise until it locks into place (C).

### 3.2 Connecting the breathing circuit

#### To connect the breathing circuit

- 1 Connect the breathing circuit as shown on page 2 or page 3, as appropriate.
- 2 Connect a flow sensor or a pressure line (nCPAP, nCPAP-PC modes only) as described next.

### 3.3 Connecting the neonatal flow sensor or pressure line

To ...	See ...
Connect a flow sensor	Section 3.3.1
Connect a pressure line	Section 3.3.2

#### 3.3.1 Connecting a flow sensor

For use with all modes other than nCPAP and nCPAP-PC.

- 1 Connect the flow sensor to the breathing circuit. Ensure the end marked PATIENT faces the patient.
- 2 Connect the blue and clear tubes to the same-colored flow sensor connectors on the ventilator. Blue to blue, clear to white/silver.
- 3 Perform the preoperational checks as described in your ventilator *Operator's Manual*.

### 3.3.2 Connecting the pressure line (nCPAP and nCPAP-PC)

Use the pressure line with the breathing circuit when using the nCPAP or nCPAP-PC modes. Do not use a flow sensor.

The pressure is measured at the T-piece adapter in the inspiratory line close to the patient, or (if available) over the optional pressure measuring connection at the Y-piece of the breathing circuit.

#### To connect a pressure line

- 1 Using an adapter, connect the pressure line to the small inlet at the top of the T- or Y-piece, whichever is used.
- 2 Connect the pressure line to the blue flow sensor connector on the ventilator.
- 3 Perform the preoperational checks as described in your ventilator *Operator's Manual*.

### 4 Preoperational check

Perform all of the required preoperational checks before patient use, including a Leak/Tightness test and necessary calibrations. For detailed instructions, see your ventilator *Operator's Manual*.

### 5 Delivery

The components are delivered clean and ready for clinical use.

## 6 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

## 7 Environmental requirements and storage conditions

### Environmental requirements and storage conditions

Operating temperature & humidity	-15°C to 40°C (5°F to 104°F) 5% to 95% relative humidity, noncondensing
----------------------------------	--

Storage temperature & humidity	-20°C to 60°C (-4°F to 140°F) 5% to 95% relative humidity, noncondensing
--------------------------------	---

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Specifications are subject to change without notice. For all proprietary and third-party trademarks used by Hamilton Medical AG, see [www.hamilton-medical.com/trademarks](http://www.hamilton-medical.com/trademarks). © 2021 Hamilton Medical AG. All rights reserved.
Hamilton Medical AG
medin Medical Innovations GmbH