

Bonaduz, May 28, 2021

Product Security Advisory – PSA-2021-02

1. Introduction

During our regular vulnerability monitoring procedure we discovered an issue with multiple operating systems (<https://us-cert.cisa.gov/ics/advisories/icsa-21-119-04>).

Hamilton Medical is aware, that the “BadAlloc” vulnerability was found in VxWorks and ARM CMSIS operating systems (RTOS) which is used by some of our products.

Hamilton Medical has conducted an impact analysis and risk assessment, and has identified no patient risks associated with these vulnerabilities.

2. Affected Products

The following Hamilton Medical products are affected by the identified vulnerabilities:

- HAMILTON-C1/T1/MR1
- HAMILTON-C2
- HAMILTON-C3
- HAMILTON-C6
- HAMILTON-G5/S1

3. Vulnerability Description

The following section refers to the US-CERT advisory ICSA-21-119-04, which is published online at: <https://www.us-cert.gov/ics/advisories>

1. Integer overflow or wraparound CWE-190

ARM CMSIS RTOS2 versions prior to 2.1.3 are vulnerable to integer wrap-around `inosRtxMemoryAlloc` (local `malloc` equivalent) function, which can lead to arbitrary memory allocation, resulting in unexpected behavior such as a crash or injected code execution.

[CVE-2021-27431](#) has been assigned to this vulnerability. A CVSS v3 base score of 7.3 has been calculated; the CVSS vector string is ([AV:N/AC:L/PR:N/UI:N/S:U/C:L/I:L/A:L](#)).

2. Integer overflow or wraparound CWE-190

Wind River VxWorks several versions prior to 7.0 firmware is vulnerable to weaknesses found in the following functions; `calloc(memLib)`, `mmap/mmap64 (mmanLib)`, `cacheDmaMalloc(cacheLib)` and `cacheArchDmaMalloc(cacheArchLib)`. This improper memory assignment can lead to arbitrary memory allocation, resulting in unexpected behavior such as a crash or a remote code injection/execution.

[CVE-2020-35198](#) and [CVE-2020-28895](#) have been assigned to this vulnerability. A CVSS v3 base score of 7.3 has been calculated; the CVSS vector string is ([AV:N/AC:L/PR:N/UI:N/S:U/C:L/I:L/A:L](#)).

4. Impact

An impact analysis and risk assessment has identified no patient risks associated with the vulnerabilities. All vulnerabilities require physical access to the device.

- **Integer overflow or wraparound**

Devices have limited memory available and the application will never alloc a memory amount which leads to a memory corruption. There is no interface that allows users to remotely control memory allocation.

All Hamilton Medical ventilators except HAMILTON-C1/T1 > V3.0.0 are not connected via Ethernet to a network and/or the internet. There is an RS232 connection possibility for outgoing data only (no bi-directional communication).

HAMILTON-C1/T1 > 3.0.0 may be connected via Ethernet or WiFi to a network and/or the internet. The connection is protected by an independent module between the ventilator and the network. This module does not use an affected operating system.

5. Mitigation

No customer interventions are necessary. In general, Hamilton Medical recommends:

- Prevent physical access to the device by unauthorized personnel
- Pay attention to notifications, alarms, and alerts
- Use only software from official Hamilton Medical channels and have it installed by an authorized technician

As part of the continuous product care and maintenance process, Hamilton Medical will provide fixes for vulnerabilities. Due to the risk profile, as well as technical feasibility, Hamilton Medical will not provide an immediate fix for these vulnerabilities. The patches provided by the RTOS vendors will be included in the regular software updates.

6. Contact Information

For further questions regarding the impact or mitigation of the vulnerabilities, please contact your local sales or service manager. Alternatively, you can directly contact Hamilton Medical at:

product-security.med.global@hamilton-medical.com

Revisions

| Revision | Publication date | Updates |
|----------|------------------|-----------------|
| 1.0 | 2021-05-26 | Initial version |

Released by

| | |
|---|--|
| Hollenberg Sebastian Product Management | Annemarie Weideli Team Leader Regulatory Affairs |
| | |