
<b>Instructions for use</b>	HAMILTON-BC8022, Breathing circuit set	4
<b>Gebrauchsanweisung</b>	HAMILTON-BC8022, Beatmungsschlauchset	8
<b>Instructions d'utilisation</b>	Kit de circuit respiratoire HAMILTON-BC8022	12
<b>Instrucciones de uso</b>	del equipo respiratorio HAMILTON-BC8022	16
<b>Istruzioni per l'uso</b>	HAMILTON-BC8022, set circuito paziente	20
<b>Instruções de uso</b>	para HAMILTON-BC8022, circuito de respiração	24
<b>使用说明</b>	HAMILTON-BC8022, 呼吸装置	28
<b>使用説明書</b>	HAMILTON-BC8022、呼吸回路セット	31
<b>Инструкция по эксплуатации</b>	комплекта дыхательного контура HAMILTON-BC8022	34
<b>Kullanım Talimatları</b>	HAMILTON-BC8022, Solunum devresi seti	38

**REF** 260161


CE 0197

ERC

**HAMILTON**  
**MEDICAL**

Breathing circuit diagram | Schematische Darstellung des Beatmungsschlauchsystems | Schéma du circuit respiratoire | Diagrama del circuito respiratorio | Diagramma del circuito paziente | Esquema do circuito de respiração | 呼吸管路图表 | 呼吸回路の図 | Схема дыхательного контура | Solunum devresi şeması

- 1** Ventilator inspiratory limb  
 Inspirationsschenkel für das Beatmungsgerät  
 Branche inspiratoire du ventilateur  
 Rama inspiratoria del respirador  
 Branca inspiratoria del ventilatore  
 Alça inspiratória do respirador  
**呼吸机吸气肢**  
 人工呼吸器側の吸気回路  
**Патрубок вдоха со стороны аппарата ИВЛ**  
 Ventilator inspiratuar parçasi
- 2** Humidifier chamber  
 Befeuchtkammer  
 Réservoir de l'humidificateur  
 Câmara del humidificador  
 Camera dell'umidificatore  
 Câmara do umidificador  
**湿化器水箱**  
 加温加湿器チャンバ  
**Камера увлажнителя**  
 Nemlendirici haznesi
- 3** Patient inspiratory limb  
 Inspirationsschenkel für den Patienten  
 Branche inspiratoire patient  
 Rama inspiratoria del paciente  
 Branca inspiratoria del circuito paziente  
 Alça inspiratória do paciente  
**病人吸气肢**  
 患者側の吸気回路  
**Патрубок вдоха со стороны пациента**  
 Hastaya giden inspiratuar parça
- 4** Y-piece (OD22/OD22 / ID15/OD22)  
 Y-Stück (AD22/AD22 – ID15/AD22)  
 Pièce en Y (DE 22/DE 22/DI 15/DE 22)  
 Pieza en Y (DE22/DE22/DI15/DE22)  
 Raccordo a Y (D.E. 22/D.E. 22; D.I. 15/D.E. 22)  
 Peça em "Y" (DE22/DE22 / DI15/DE22)  
**Y形管 (外径 22/ 外径 22 / 内径 15/ 外径 22)**  
 Yピース (外径 22/ 外径 22、内径 15/ 外径 22)  
**Y-образный коннектор (НД 22/НД 22/ВД 15/ НД 22)**  
 Y parçasi (DÇ22/DÇ22 / İÇ15/DÇ22)
- 5** Water feed line with spike  
 Wasserzufuhrleitung mit Anstechdorn  
 Ligne d'alimentation en eau avec pointe  
 Vía de suministro de agua con punta  
 Linea di alimentazione dell'acqua con perforatore  
 Linha de alimentação da água com grampo  
**带有尖状灌注器的供水管线**  
 スパイク付き水供給管  
**Линия подачи воды со штекерным соединителем**  
 Sivri uçlu su besleme hattı
- 6** Ventilator expiratory limb  
 Expirationsschenkel für das Beatmungsgerät  
 Branche expiratoire du ventilateur  
 Rama espiratoria del respirador  
 Branca espiratoria del ventilatore  
 Alça expiratória do respirador  
**呼吸机呼气肢**  
 人工呼吸器側の呼気回路  
**Патрубок выдоха со стороны аппарата ИВЛ**  
 Ventilator ekspiratuar parçasi
- 7** Elbow adapter, 90° (ID22/OD22)  
 Ellenbogenadapter, 90° (ID22/AD22)  
 Adaptateur coudé, 90° (DI 22/DE 22)  
 Adaptador acodado, 90° (DI22/DE22)  
 Adattatore a gomito, 90° (D.I. 22/D.E. 22)  
 Adaptador de cotovelo, 90° (DI22/DE22)  
**弯接头, 90° (内径 22/ 外径 22)**  
 エルボアダプタ、90° (内径 22/ 外径 22)  
**Угловой адаптер, 90° (ВД 22/НД 22)**  
 Dirsek adaptörü, 90° (İÇ22/DÇ22)

## Conventions

### ⚠ WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

### ⚠ CAUTION

*Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.*

### NOTICE

Emphasizes information of particular importance.

## 1 Overview

The HAMILTON-BC8022 breathing circuit set (PN 260161) can only be used with the HAMILTON-H900 humidifier. For an overview of the breathing circuit set and its components, see page 2.

### 1.1 Intended use

The HAMILTON-BC8022 breathing set is intended to be used together with compatible Hamilton Medical respiratory gas humidifiers during invasive and noninvasive mechanical ventilation of adult and pediatric patients > 10 kg.

## 2 Safety information

### ⚠ WARNING

- Reprocessing of Hamilton Medical single-use products can affect the product properties and may cause injury to the patient. For example, a change to the surface structure during reprocessing may lead to a change in the tear strength or cause actual cracking. Furthermore, an altered surface structure may result in, for example, a microbial

aggregation of spores, allergens, and pyrogens, or cause an increase in the number of particles released as a result of chemical changes in the material properties.

- Hamilton Medical does *not* assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user. Any liability on the part of Hamilton Medical for any damages resulting from the use of reprocessed or reused single-use parts, or use of any non-Hamilton Medical accessories and consumables, is excluded.
- The breathing circuit set must *not* be covered by any objects, such as sheets, towels, and so on.
- Use caution when positioning the breathing circuit. Place the breathing circuit in such a way that condensation cannot reach the patient.
- Follow hospital infection control procedures for all tests and use.

- Discard the product if it fails two subsequent preoperational checks or if there is any sign of damage.
- Adding attachments or other parts/assemblies to a breathing system can change the pressure gradient across the ventilator, which can adversely affect ventilator performance and patient safety.
- Connect only a single breathing circuit at a time. Do *not* connect any circuit as an extension to another one. Doing so can seriously harm or kill the patient.
- Route the breathing circuit without tension and without any kinks from the ventilator or humidifier to the patient, and protect against unintended influences.
- Do *not* use heated breathing circuits without existing gas flow. If the gas supply is interrupted, the humidifier must be turned off or be in **Standby** mode.
- Do *not* touch the hot plate or the bottom of the chamber. The surfaces can reach a temperature of over 85°C. These hot surfaces radiate heat.
- Heated limbs must *not* be placed directly on the patient's skin. Attach breathing circuits or tubing holders appropriately to avoid mechanical forces on the ET tube.

- Only fill the humidifier chamber with sterile, demineralized water that meets the hospital's hygiene requirements.
- Do *not* operate the device if the water level in the humidifier chamber exceeds the marked maximum.
- Check the breathing circuit for condensation regularly and drain if required.
- The humidifier must always be positioned below patient level. Do *not* operate the humidifier at an angle in excess of 10°.
- Do *not* use this breathing circuit set with humidifiers other than the HAMILTON-H900.

**⚠ CAUTION**

- *(USA only) Federal law restricts this device to sale by or on the order of a physician.*
- *Check the breathing set for damage prior to use. Discard the breathing set if there is any sign of damage.*
- *Handle used products as contaminated goods according to local laws and regulations or hospital internal procedures.*
- *The temperature symbol on the temperature probe must always face up and be visible.*

**NOTICE**

- Before use, read the *HAMILTON-H900 Humidifier Instructions for Use*.
- The refill water should *not* be warmer than 37°C.
- Ensure that the water supply to the humidifier chamber is functioning properly.
- Set the appropriate ventilator alarms.
- Do *not* apply pressure to the water reservoir once the water refill has automatically stopped and the water has reached the maximum level.
- Ensure the humidifier is placed in **Standby** or turned off before disconnecting any components.
- Do *not* add any drugs or medication directly to the water in the water chamber. If the HAMILTON-H900 is used in combination with any medical gases or nebulized medications, follow the *Instructions for use* of the supplied application and make sure that it is suitable for use with active humidification.
- Use the stickers provided to indicate the date on which to replace the breathing circuit set.

- Check the stability of all connections, prior to use.
- Incompatible parts can result in degraded performance, which can affect safety.

### 3 About the breathing circuit and breathing circuit sets

The breathing circuit comprises the limbs with integrated heater wire and integrated temperature probe. The breathing circuit set comprises the breathing circuit and the humidifier chamber with water feed line.

### 4 Connecting the breathing circuit

Note that the limb connectors on the humidifier combine electrical connections with breathing circuit connectors. Ensure proper orientation of electrical contacts on the breathing circuit connectors to match the connecting element on the humidifier.

For an overview of the breathing circuit set and its components, see page 2.

### To connect the breathing circuit set to the humidifier:

- 1 Insert the humidifier chamber completely into the humidifier until it clicks into place.  
To remove it, pull the humidifier chamber out of the humidifier.
- 2 Connect the ventilator inspiratory limb first to the humidifier, and then to the inspiratory port on the ventilator.
- 3 Connect the ventilator expiratory limb to the expiratory port on the ventilator and the Y-piece.
- 4 Connect the patient inspiratory limb first to the humidifier chamber, and then to the Y-piece.
- 5 Insert the humidifier chamber water feed spike into the water reservoir.

### 5 Delivery

The components are delivered clean and ready for clinical use.

### 6 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

## 7 Technical specifications

Parameter	Specification <sup>1</sup>
Flow range	4 to 120 l/min
Connectors <sup>2,3</sup> Ventilator: Patient:	ID22 ID15/OD22
Compliance (at 60 cmH <sub>2</sub> O)	2 ml/cmH <sub>2</sub> O
Compressible volume	1600 ml
Smallest internal diameter	19 mm
Water volume (min/max)	80 ml / 200 ml
Flow resistance (at 30 l/min) Inspiratory: Expiratory:	0.01 cmH <sub>2</sub> O/(l/min) 0.007 cmH <sub>2</sub> O/(l/min)
Flow resistance (at 15 l/min) Inspiratory: Expiratory:	0.008 cmH <sub>2</sub> O/(l/min) 0.005 cmH <sub>2</sub> O/(l/min)
Gas leakage (at 200 cmH <sub>2</sub> O)	< 40 ml/min
Maximum operating pressure	20 kPa
Circuit length	1.95 m (approx. 6 ft)

Parameter	Specification <sup>1</sup>
Compatibility	HAMILTON-H900 humidifier
Operating temperature	Recommended: 18°C to 26°C (65°F to 78°F)
Storage temperature	-20°C to 60°C (-4°F to 140°F)
Operating and storage humidity	5% to 95% relative humidity, noncondensing
Maximum operating duration	28 days
Conforms to	IEC 60601-1:2005/ A1:2012 ISO 80601-2- 74:2017 EN ISO 5356-1:2015 EN ISO 5367:2014 IEC 62366-1:2015 ISO 10993-1:2018 ISO 13485:2016

1. All values presented are worst case scenarios.
2. Internal diameter (ID) and outer diameter (OD) in mm.
3. 15 mm and 22 mm conical connectors according to ISO 5356-1.

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Specifications are subject to change without notice. For all proprietary trademarks and third-party trademarks used by Hamilton Medical AG, see [www.hamilton-medical.com/trademarks](http://www.hamilton-medical.com/trademarks). © 2021 Hamilton Medical AG. All rights reserved.
Hamilton Medical AG
medin Medical Innovations GmbH
Hamilton Medical, Inc.