HAMILTON-H900

Instructions for use
Gebrauchsanweisung
Instructions d'utilisation
Instrucciones de uso
Istruzioni per l'uso
Instruções de uso
Инструкции по эксплуатации
Kullanım talimatları
使用说明

REF 950001 (230 V) | 950004 (115 V) | 950008 (100 V)

Software version | Softwareversion | Version logicielle | Versión de software
Versione software | Versão de software | Версия программного обеспечения
Yazılım sürümü | 软件版本 1.10x

624431/10 | 2023-07-14

CE 0197 EAC

HAMILTON
MEDICAL

English	5
Deutsch	45
Français	87
Español	129
Italiano.....	169
Português.....	211
Русский	251
Türkçe.....	295
中文.....	335

Instructions for use

HAMILTON-H900

2023-07-14

624431/10

© 2023 Hamilton Medical AG. All rights reserved. Printed in Switzerland.

No part of this publication may be reproduced, stored in a database or retrieval system, or transmitted in any form or by any means, electronic, mechanical, or by photocopying, recording, or otherwise, without prior written permission of Hamilton Medical AG.

This document may be revised, replaced, or made obsolete by other documents by Hamilton Medical AG at any time and without notice. Ensure that you have the most current applicable version of this document; if in doubt, contact the technical support department of Hamilton Medical AG, Switzerland. While the information set forth herein is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this document shall limit or restrict in any way Hamilton Medical AG's right to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, Hamilton Medical AG has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

The equipment must be operated, serviced, or upgraded only by trained professionals. Hamilton Medical AG's sole responsibility with respect to the equipment and its use is as stated in the limited warranty provided in this manual.

Hamilton Medical AG shall not be liable for any loss, cost, expense, inconvenience, or damage that may arise out of misuse of the product, or if non-Hamilton Medical AG parts were used when replacing parts, or if serial numbers were amended, deleted, or removed.

If returning parts to Hamilton Medical AG, be sure to use the standard Hamilton Medical returned goods authorization (RGA) procedure. Disposal of parts shall follow all local, state, and federal regulation with respect to environmental protection.

For all proprietary as well as third-party trademarks used by Hamilton Medical AG, see www.hamilton-medical.com/trademarks.

Manufacturer

Hamilton Medical AG
Via Crusch 8, CH-7402 Bonaduz,
Switzerland
Phone: (+41) 58 610 10 20
Fax: (+41) 58 610 00 20
info@hamilton-medical.com
www.hamilton-medical.com

Distributor in USA

Hamilton Medical, Inc.
4655 Aircenter Circle, Reno, NV
89502
Phone: (775) 858-3200
Toll-free: (800) 426-6331
Fax: (775) 856-5621
info@hamiltonmedical.com
www.hamilton-medical.com

	Preface.....	10
1	Safety information	11
1.1	Electromagnetic susceptibility	11
1.2	Electrical power	12
1.3	Fire and other hazards	12
1.4	General operation and setup	12
1.5	Breathing circuits and water chamber	13
1.6	Humidifier and breathing circuit positioning.....	14
1.7	Monitoring and alarms	14
1.8	Maintenance	14
2	Overview.....	15
2.1	About the HAMILTON-H900 humidifier	15
2.2	Physical descriptions	16
2.2.1	About the humidifier	16
2.2.2	About the main display.....	16
2.2.3	About the humidifier breathing circuits.....	18
2.2.4	About the status indicators on the humidifier	18
2.3	Navigating windows and controls	19
2.3.1	Using the control buttons	19
2.3.2	Using the sliders	19
3	Preparing the humidifier for use	19
3.1	Connecting to a power source.....	20
3.2	Connecting the breathing circuit to the humidifier.....	20
3.3	Turning the humidifier on and off.....	20
3.4	Testing the alarms	21
3.5	Specifying humidifier settings	21
3.6	Connecting the patient.....	21

4	Working with the HAMILTON-H900 humidifier	21
4.1	Entering/exiting Standby	21
4.2	About the humidifier operating modes	22
4.2.1	Changing modes	22
4.2.2	Auto and manual control settings	23
4.3	Changing humidity using temperature controls	24
4.4	Monitoring humidification	25
4.4.1	Displaying additional monitoring parameters	26
4.5	Day and Night display brightness settings	26
5	Responding to humidifier alarms	26
5.1	Temporarily silencing an alarm	28
5.2	Adjusting alarm loudness	28
5.3	Troubleshooting alarms	29
6	Maintenance	31
6.1	Cleaning, disinfection, and sterilization	31
6.2	Preventive maintenance	32
7	Configuration	33
7.1	Accessing Configuration mode	33
7.2	Changing configuration settings	33
7.3	Restoring the factory default settings	34
8	Parts and accessories	34
9	Specifications	35
9.1	Physical characteristics	35
9.2	Environmental requirements	36
9.3	Electrical specifications	36
9.4	Control settings	37
9.5	Monitored parameters	38
9.6	Configuration	38

9.7	Technical performance data.....	39
9.8	Essential performance.....	40
9.9	Functional description of humidification system	40
9.10	Symbols used on device labels and packaging.....	40
9.11	Standards and approvals.....	41
9.12	Disposal.....	41
9.13	Warranty	42
9.13.1	Miscellaneous.....	43

Preface

Be sure to read the documentation before using the device or accessories.

Conventions used in this guide

In this manual:

- Button names are shown in a **bold** font.
- The graphics shown may not exactly match what you see in your environment.
- Not all features are available in all markets.

Safety messages are displayed as follows:

WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

CAUTION

Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.

NOTICE

Emphasizes information of particular importance.

In tables, safety messages are indicated as follows:

 **WARNING!**

 **CAUTION!**

 **NOTICE!**

To download the latest version of this manual or other documents, visit the MyHamilton website:

<https://www.hamilton-medical.com/MyHamilton>

Hamilton Medical offers the Hamilton Medical College, which provides a variety of learning modules free of charge. To register, go to:
<http://college.hamilton-medical.com/>

CAUTION

(USA only): Federal law restricts this device to sale by or on the order of a physician.

Intended use

The HAMILTON-H900 humidifier is intended for respiratory gas conditioning during invasive and noninvasive mechanical ventilation.

The device is intended for use in the hospital and institutional environment where health care professionals provide patient care.

The HAMILTON-H900 humidifier is a medical device intended for use by qualified, trained personnel under the direction of a physician and within the limits of its stated technical specifications.

1 Safety information

This section provides safety information related to setting up, operating, and servicing the humidifier.

Carefully review all parts of this safety section before setting up and using the humidifier.

Be sure to review the Instructions for Use before using the humidifier.

Be sure to read the Instructions for Use provided with any devices and accessories used with the humidifier before use.

If you have questions about any of the information in this manual, contact your Hamilton Medical representative or technical service personnel.

1.1 Electromagnetic susceptibility

WARNING

- **MR UNSAFE.** Keep away from magnetic resonance imaging (MRI) equipment. The HAMILTON-H900 poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.
- Functioning of the humidifier may be impaired by the operation of high-frequency surgical equipment, microwaves, shortwaves, or strong magnetic fields in close proximity.
- To prevent increased emissions, decreased immunity, or interrupted operation of the HAMILTON-H900 humidifier or any accessories, use only accessories or cables that are expressly stated in this manual or in the Hamilton Medical e-catalog.

- Additional equipment connected to medical electrical equipment must comply with the respective IEC or ISO standards (for example, IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the requirements for medical electrical systems (see IEC 60601-1, clause 16).
- Anybody connecting additional equipment to medical electrical equipment configures a medical system and is, therefore, responsible that the system complies with the requirements for medical electrical systems. Note that local laws take priority over the above-specified requirements. If you have questions about how to proceed, consult your Hamilton Medical representative or technical service department.
- The humidifier is *not* protected against the effects of the discharge of a cardiac defibrillator.

NOTICE

- The HAMILTON-H900 humidifier requires special precautions regarding EMC (Electromagnetic Compatibility), and must be installed and put into service according to the EMC information provided in the *HAMILTON-H900 EMC Declarations*.
- Portable and mobile RF communications equipment can affect the HAMILTON-H900 humidifier.

1.2 Electrical power

WARNING

- To minimize the risk of electric shock, plug the humidifier power cord into an appropriate grounded power receptacle. It is the hospital's responsibility to ensure that the receptacle is properly grounded (earth).
- *Do not* use if the power cord is damaged.
- In case of incorrect voltage, power failure, or disconnection from the power supply, the humidifier emits an audible whistling sound. Turn the device off immediately and check for correct voltage.
- Ensure the correct voltage is used.

1.3 Fire and other hazards

WARNING

Do *not* use the humidifier in hazardous areas or with flammable gases.

1.4 General operation and setup

WARNING

- Before using the humidifier on the patient, verify that the breathing circuit is correctly connected to the ventilator as follows:
 - The blue inspiratory limb is connected to the *To patient* inspiratory port.
 - The white expiratory limb is connected to the *From patient* expiratory port.
- Use only parts and accessories specified in Section 8. Doing so ensures proper humidifier operation and prevents possible patient injury.
- Do not use the humidifier at an altitude above 4,000 meters or outside a temperature range of 10°C to 40°C. Using the humidifier above this altitude or outside of this temperature range can affect the quality of the therapy or injure the patient.
- Incorrect humidity or temperature settings can harm the patient.
- Do *not* use the humidifier during transfer of a ventilated patient outside the hospital.
- Turn ON the humidifier *after* turning on the ventilator.
- Turn OFF the humidifier *before* turning off the ventilator.
- Do not cover the humidifier IR measurement cells.
- Invasive (INV) mode must *only* be used with intubated and tracheostomized patients.
- Take additional precautions in case of allergic reactions.
- Regularly check the breathing circuit for condensation and drain it, if required.
- Device humidification performance may be negatively affected by the simultaneous use of a nebulizer.
- Do *not* touch the hot plate or the bottom of the water chamber. The surfaces can reach a temperature of over 85°C. These hot surfaces radiate heat.

CAUTION

- Position the humidifier where the primary power supply can be easily disconnected.
- Use only parts and accessories specified in Section 8, and in the Hamilton Medical e-catalog, or that are specified as being compatible with the humidifier.

Doing so ensures proper humidification, avoids degraded performance, and keeps your warranty in force.

- If the humidifier is turned Off then On, it will automatically start in Invasive (INV) mode unless configured differently.
- Check the breathing circuit set for damage prior to use. Discard the breathing circuit set if there is any sign of damage.
- The humidifier heating element and heater wires are automatically turned on when the water chamber and breathing circuits are correctly mounted, and the humidifier is turned on or is in Standby.
- If the ambient or gas inlet temperature, and/or flow rate is outside the recommended range, physiological humidity levels may not be achieved.
- The humidifier still may be powered, even if the display is not illuminated. See Table 1.

NOTICE

- Before using a breathing circuit set on a patient, you must perform a tightness test. If this test fails, it may be repeated once. After two consecutive failed tests, the breathing set must be replaced by a new one.

- A fine, breath-dependent condensation (fogging) forming in the limb, flex tube, or flow sensor indicates that the humidity is properly set.
- All actions generate an audible sound.

1.5 Breathing circuits and water chamber

WARNING

- Failure to correctly connect the breathing circuit to the ventilator can injure the patient.
- To prevent accidental disconnection of the breathing circuit during use or transport, use only breathing circuits that comply with ISO 5367 or ISO 80601-2-74.
- Only fill the water chamber with sterile, demineralized water that meets the hospital's hygiene requirements.
- Do not tilt the water chamber.
- Do not fill any drug into the humidifier chamber.
- Ensure the water level does not exceed the maximum fill level.

NOTICE

- If the humidifier does not detect the breathing circuit, replace the components.
- The refill water should be no warmer than 37°C.
- Ensure that the water supply to the chamber is functioning properly.
- Ensure the humidifier is in Standby before disconnecting any components.

1.6 Humidifier and breathing circuit positioning

WARNING

- The humidifier must *always* be positioned below patient level.
- Do *not* operate the humidifier at an angle in excess of 10° relative to the floor.
- Be sure to route the breathing circuit without tension and without any kinks from the humidifier to the patient.
- Heated breathing limbs must *not* be placed directly on the patient's skin.
- If the humidifier is used adjacent to, or stacked on, other medical electrical equipment, verify the humidifier's normal operation in the configuration in which it will be used.
- Place the breathing circuit in such a way that liquid condensate cannot reach the patient.
- Attach breathing circuits or tubing clips appropriately to avoid mechanical forces on the ET tube.

NOTICE

- Check the stability of all connections prior to use.
- All limb connectors on the humidifier combine electrical connections with breathing circuit connectors. Ensure proper orientation of electrical contacts to match the connecting element on the humidifier.

1.7 Monitoring and alarms

NOTICE

- When an alarm is active, the display alternates between showing active alarm symbols and the current chamber exit temperature.
- Upon turning off the humidifier, the alarm loudness returns to the default value of 5.
- All technical alarms, detailed technical information, and maintenance procedures are described in the *HAMILTON-H900 Service Manual*.

1.8 Maintenance

WARNING

- Modifications to the humidifier are *not* permitted.
- *Always* disconnect the humidifier from electrical power before cleaning and disinfection to reduce the risk of electrical shock.

CAUTION

- *Handle used breathing circuits and water chambers as contaminated goods according to local laws and regulations or hospital internal procedures.*
- *To ensure proper servicing and to prevent possible physical injury, only Hamilton Medical authorized service personnel may service the humidifier using information provided in the HAMILTON-H900 Service Manual.*
- *Use only replacement parts supplied by Hamilton Medical.*

- Do not attempt service procedures other than those specified in the *HAMILTON-H900 Service Manual*.
- Use only approved cleaning agents for cleaning and disinfection.

NOTICE

- For specific information on cleaning, disinfecting, and sterilizing auto-clavable (reusable) accessories and components, refer to the appropriate *Reprocessing Guide* and *Instructions for Use* provided with each part.
- (USA only) Use only EPA-registered cleaning agents for cleaning and disinfection.

2 Overview

This guide provides information about the set up and use of the HAMILTON-H900 humidifier.

The humidifier is designed for use with Hamilton Medical, as well as third-party, ventilators.

2.1 About the HAMILTON-H900 humidifier

The HAMILTON-H900 humidifier heats and humidifies respiratory gases for adult, pediatric, and neonatal patients.

It supports the following common respiratory modes:

- Invasive ventilation
- Noninvasive ventilation
- High flow oxygen therapy

The system comprises the humidifier housing, display, controls, water chamber, heater plate, and electrical connections for a heated breathing circuit.

It offers the following main features:

- Adjustable temperature and humidity controls
- Monitoring and alarm capabilities
- Integration of humidifier monitored data and controls with supported Hamilton Medical ventilators¹

¹ Not available in all markets.


2.2 Physical descriptions

This section provides an overview of the humidifier and related breathing circuit sets.

2.2.1 About the humidifier


Figures 1 through 5 provide an overview of the humidifier.

Figure 1. HAMILTON-H900 humidifier, front view

- | | |
|--------------------------------------|-----------------------------|
| 1 Ventilator inspiratory limb (blue) | 5 Heater plate |
| 2 Patient inspiratory limb (blue) | 6 Water chamber limb (blue) |
| 3 Water feed tube | 7 Front panel and display |
| 4 Visual alarm indicator | |

Figure 2. HAMILTON-H900 humidifier, rear view

- | | |
|--------------------|-----------------------------------|
| 1 Mounting bracket | 3 Potential equalization terminal |
| 2 RS-232 port | 4 AC power socket |

2.2.2 About the main display

You can directly access all settings, modes, alarms, and controls from the main display during humidification.

Figures 3 through 5 illustrate the front panel controls, indicators, and display.

Figure 3. HAMILTON-H900 humidifier, front panel controls


- | | |
|--|--------------------------|
| 1 Temperature monitoring window button | 5 Audio pause key |
| 2 Temperature control sliders | 6 Power/Standby key |
| 3 Exp. temp increase button | 7 Auto button |
| 4 INV/NIV/HiFlow mode button | 8 Display (see Figure 5) |

Figure 4. HAMILTON-H900 humidifier, front panel status indicators

- | | |
|-------------------------------------|----------------------------------|
| 1 Limb connection status indicators | 3 Power/Standby status indicator |
| 2 Audio pause indicator | |

Figure 5. HAMILTON-H900 humidifier, display

- | | |
|-------------------------------------|----------------------------------|
| 1 Ventilator connection enabled | 6 Current chamber exit temp |
| 2 Max. temperature gradient setting | 7 Chamber exit temp setting |
| 3 Temperature gradient setting | 8 Min. chamber exit temp setting |
| 4 Min. temperature gradient setting | 9 Active modes/settings |
| 5 Max. chamber exit temp setting | |

All of the elements on the display are shown for clarification purposes only; during operation, they may not all appear at the same time.

2.2.3 About the humidifier breathing circuits

The HAMILTON-H900 humidifier supports a variety of single and dual limb breathing circuits for adult, pediatric, and neonatal patients.

For details, see your breathing circuit *Instructions for use*.

2.2.4 About the status indicators on the humidifier

Indicators on the front of the humidifier show important humidifier status information. See Figure 4.

Table 1. Status indicators

Symbol	Description
Power/Standby key with status indicator	
	<i>Green:</i> The humidifier is in operation.
	<i>Orange:</i> The humidifier is in Standby.
	<i>Blue:</i> The humidifier is turned off and connected to a power source.
	<i>Unlit:</i> The humidifier is turned off and <i>not</i> connected to a power source.
Limb connection status indicator	
	<i>Green:</i> The limb is connected correctly.
	<i>Orange:</i> The humidifier is testing the limb connection.
	<i>Red:</i> The limb connection is either faulty or a limb is not connected.
Audio pause key with indicator	
	<i>Red:</i> Audio pause is active.

2.3 Navigating windows and controls

This section describes how to use the humidifier buttons and sliders.

2.3.1 Using the control buttons

The buttons on the front panel (Figure 3) let you select modes, change settings, and see an overview of key temperature settings.

Table 2. Control buttons (see Figure 3)

Button	Description
	Exp. temperature increase Increases the temperature in the expiratory limb.
	Auto Enables Auto mode. See Section 4.2.2.
	INV/NIV/HiFlow operating modes Switches between INV, NIV, and HiFlow modes. See Section 4.2.
	Temperature monitoring window Displays the current Y-piece temperature, chamber exit temperature, and temperature gradient values. See Section 4.4.

2.3.2 Using the sliders

In manual mode you can change the temperature settings using the control sliders (Figure 3).

Table 3. Temperature control sliders

Slider	Description
Temperature gradient slider	
	Adjusts the temperature difference between the chamber exit and the Y-piece. See Section 4.3.
Chamber exit temp slider	
	Adjusts the chamber exit temperature. See Section 4.3.

3 Preparing the humidifier for use

Before proceeding, review the safety information in Section 1.

Preparing the humidifier for use comprises the following steps:

- A one-time configuration of settings by a medical technician (see Table 4)
- For each new patient, humidifier setup tasks performed by medical caregivers (see Table 5)

Table 4. Configuration and setup by a medical technician

To ...	See ...
<i>These one-time initial configuration tasks are completed by technical personnel.</i>	
Mount and position the humidifier	See the applicable <i>Installation Guide</i> for the ventilator or trolley
Configure the humidifier settings	Section 7
<i>Optional:</i> Connect the humidifier to the RS-232 COM port of any supported Hamilton Medical ventilator	See the applicable <i>Operator's Manual</i> for the ventilator.

Table 5. Configuration and setup by medical caregivers

To ...	See ...
<i>The following tasks are performed by medical personnel caring for patients.</i>	
Connect to a power source	Section 3.1
Connect the breathing circuit	Section 3.2
Turn on the humidifier	Section 3.3
Test the alarms	Section 3.4
Specify humidifier settings	Section 3.5
Connect the patient	Section 3.6

3.1 Connecting to a power source

Always check the reliability of the primary power outlet before plugging in the humidifier and the ventilator.

To connect the humidifier to a power source

- ▶ Connect the humidifier to an outlet that supplies AC power.²

3.2 Connecting the breathing circuit to the humidifier

Hamilton Medical provides a variety of breathing circuits for adult, pediatric, and neonatal patients.

For details, see the specific breathing circuit *Instructions for Use*.

3.3 Turning the humidifier on and off

When the humidifier is turned on, it automatically runs in the configured default mode.

Be sure to turn the ventilator on *before* turning on the humidifier.

To turn on the humidifier

- ▶ Press and hold  (Power/Standby) for approximately 3 seconds.

The humidifier runs a self-test. Once complete, the humidifier automatically starts in invasive (INV) mode unless configured differently.

² If using the humidifier with a supported Hamilton Medical ventilator, connect the humidifier power cable to the dedicated power socket on the ventilator.

To turn off the humidifier

- ▶ Press and hold  for approximately 3 seconds.

Be sure to turn the humidifier off *before* turning off the ventilator.

3.4 Testing the alarms

Before connecting a new patient to the humidifier, verify that the alarms operate correctly.

Ensure that the humidifier is turned on.

To test the humidifier alarms

1. Generate an alarm by removing the water chamber from the humidifier.
2. Verify that:
 - An acoustic alarm is audible
 - The visual alarm indicator is flashing (Figure 4)
 - The corresponding alarm symbol appears on the display (Table 9)

Resolve the alarm by re-inserting the water chamber into the humidifier, and verify that the alarm is reset.

3.5 Specifying humidifier settings

Select the humidification mode for your patient.³ For details, see Section 4.2.

3.6 Connecting the patient

Connect the breathing circuit to the appropriate patient interface.

The humidifier is now set up and ready for use.

4 Working with the HAMILTON-H900 humidifier

Table 6. Operation overview

For details about ...	See ...
Turning the humidifier on/off	Section 3.3
Entering/exiting Standby	Section 4.1
Accessing humidifier controls	Section 2.3
Humidifier operating modes: INV, NIV, HiFlow	Section 4.2
Auto and Manual settings	Section 4.2.2
Changing humidity using temperature controls	Section 4.3
Monitored humidifier parameters	Section 4.4
Humidifier alarms	Section 5

4.1 Entering/exiting Standby

CAUTION

In Standby, all settings are retained and the heat output is reduced.

Standby is a waiting mode that lets you maintain current settings while the humidifier is not in operation.

When in Standby, the Chamber exit temp and Temperature gradient settings are reduced.⁴ For details, see Section 9.7.

³ When used with a Hamilton Medical ventilator that supports integrated operation, the humidifier matches the operating mode of the ventilator. For details, see your ventilator *Operator's Manual*.

⁴ During nebulization, the Temperature gradient is *not* changed during Standby.

To put the humidifier into Standby

- ▶ Briefly press  (Power/Standby).
When in **Standby**, the display shows the elapsed time the humidifier has been in Standby. The Power/Standby indicator light is orange.

To exit Standby and start humidification

- ▶ Press .
The humidifier resumes normal operation; the Power/Standby status indicator light is green.

4.2 About the humidifier operating modes

The HAMILTON-H900 offers the following operating modes: Invasive (INV), Noninvasive (NIV), and High flow oxygen therapy (HiFlow).

The operating mode selection determines the initial temperature settings, both at the water chamber exit (Chamber exit temp) and at the Y-piece (Temperature gradient), as well as the allowed temperature ranges for each of these controls.

These temperature ranges differ by patient group. See Table 18.

4.2.1 Changing modes

The humidifier automatically starts up in Invasive (INV) mode, unless configured differently in the custom default settings.

When changing to Invasive (INV) or Noninvasive (NIV) mode, the humidifier automatically sets the control settings to Auto.

To change modes

- ▶ Press  (Modes) until the desired mode (INV, NIV, HiFlow) is displayed.
The selected mode appears on the display.

Figure 6. INV mode (1) selected


Figure 7. NIV mode (1) selected



Figure 8. HiFlow mode (1) selected

4.2.2 Auto and manual control settings

You can use the humidifier in **Auto** mode, or you can manually adjust the settings. By default, the humidifier starts in **Auto** mode.

The **chamber exit temp** and **Temperature gradient** are set using either of the following methods:

- Loaded from the configured default settings (**Auto** mode)
- Set manually by the user


Automatic settings (Auto)

When set to **Auto**, the humidifier loads the configured default settings for the selected mode (see Section 7) and uses them to control the gas temperature.

At low ambient temperatures, the humidifier may automatically adjust the water chamber temperature.

When **Auto** mode is selected, the text **AUTO** appears on the main display.

Figure 9. Auto mode (1) selected


Manual settings

When you change the **Chamber exit temp** or **Temperature gradient** using the sliders, the humidifier switches to **Manual** mode. However, it still continues to automatically control the temperature to reach the specified settings.

When settings are changed manually, **Auto** no longer appears on the main display; the space is empty.

Figure 10. Manual mode (1) is in effect

4.2.2.1 Enabling Auto mode

After manually changing either of the temperature settings, you can re-enable **Auto** mode.

To re-enable Auto mode

- ▶ Press  (**Auto**).

Changing the operating mode to **INV** or **NIV** also automatically switches the humidifier to **Auto** mode.

4.3 Changing humidity using temperature controls

NOTICE

If the humidifier is in Auto mode, dragging either slider switches humidification to Manual mode.

NOTICE

In HiFlow mode, you cannot change the Temperature gradient using the slider.

You can change the Temperature gradient to a new value in Configuration. See Section 7.

You can adjust the following controls on the humidifier⁵:

Table 7. Adjustable humidifier controls

Control	Description
Chamber exit temp	<p>Temperature at the water chamber exit.</p> <p>The possible range of values for this control depends on the selected humidifier operating mode. For details, see Tables 16 and 18.</p> <p>Higher values result in higher absolute humidity.</p>
Temperature gradient	<p>The difference between the temperature at the water chamber exit and at the Y-piece.</p>
Exp. temp increase	<p>When selected, the humidifier provides additional heat in the expiratory limb.⁶</p>

The maximum allowed temperature at the patient Y-piece is 42°C. The combination of the Chamber exit temp plus the Temperature gradient cannot exceed this limit. For example, if the Temperature gradient is set to 2°C, the highest possible setting for the Chamber exit temp is 40°C.

Furthermore, the Temperature gradient setting takes precedence over the Chamber exit temp setting. This means if a combination of settings would result in a Y-piece temperature greater than 42°C, then the Temperature gradient will remain unchanged, but the Chamber exit temp will be adjusted.

For example, if Chamber exit temp is set to 40°C, you can set Temperature gradient to 3°C even though the combination exceeds 42°C. Once the Temperature gradient setting is accepted, the Chamber exit temp value is automatically set to 39°C.

To set the Chamber exit temperature


- ▶ Touch the Chamber exit temp slider at the current setting:
 - Drag left to decrease the chamber exit temperature
 - Drag right to increase the chamber exit temperature

The changes are confirmed by an audible beep and are applied immediately.

⁵ When used with a Hamilton Medical ventilator that supports integrated humidifier control, you can also change the settings directly on the ventilator. For details, see your ventilator *Operator's Manual*.

⁶ Single use, dual limb heated breathing circuits only.

Figure 11. Adjusting Chamber exit temperature

To set a Chamber exit temperature > 39°C

To increase the Chamber exit temp above 39°C, the slider must be activated twice.

1. Touch the Chamber exit temp slider at the current setting and drag right to reach the limit of 39°C.
2. Touch and drag the slider again to set the chamber exit temperature above 39°C.

To set the temperature gradient

- ▶ Touch the Temperature gradient slider at the current setting:
 - Drag right to decrease the temperature gradient
 - Drag left to increase the temperature gradient

The changes are confirmed by an audible beep and are applied immediately.

Figure 12. Adjusting Temperature gradient

For example, if the target Y-piece temperature is 39°C and the Chamber exit temp is set to 37°C, set the Temperature gradient to 2°C.

4.4 Monitoring humidification

Monitoring data is available in the Temperature monitoring window, which you can access at any time.⁷

The following parameters are monitored:

- Chamber exit temp
- Temperature gradient
- Y-piece temperature

⁷ When used with a Hamilton Medical ventilator that supports integrated humidifier monitoring, humidifier data is shown on the ventilator display. For details, see your ventilator *Operator's Manual*.

4.4.1 Displaying additional monitoring parameters

To display the Temperature monitoring window


- ▶ Press  (Temperature monitoring window).
The window closes after 30 seconds or when you press the button again.

Figure 13. Temperature monitoring window

- | | |
|-------------------------------|--------------------------------|
| 1 Current Y-piece temperature | 3 Current Temperature gradient |
| 2 Current chamber exit temp | |

4.5 Day and Night display brightness settings

The humidifier display backlight supports two different display settings. Use these settings to set the brightness of the display.

To switch between the Day and Night settings

- ▶ Press and hold  for 5 seconds.

To set the default day and night brightness, see Section 7.

5 Responding to humidifier alarms

Humidifier-related alarms notify you of conditions that require your attention. Alarms are indicated as follows:

- Graphically on the humidifier display
- With a sequence of 3 or 5 audible beeps
- The alarm lamp above the water chamber flashes
- When used with a Hamilton Medical ventilator that supports integrated humidifier monitoring, alarm messages are also indicated on the ventilator display.⁸

The alarms are categorized as high priority, medium priority, or technical fault, as described in Table 8.

See Table 9 for a complete list of alarms.

⁸ Not available on all ventilators. See your ventilator *Operator's Manual*.

Table 8. Alarm indicators

Alarm type	Alarm status indicator	Audio	Action required
High priority	Red, flashing Alarm icon is shown on the display	A sequence of 5 beeps, repeated until the alarm is reset	The humidifier needs immediate attention.
Medium priority	Yellow, flashing Alarm icon is shown on the display	A sequence of 3 beeps, repeated until the alarm is reset	The humidifier needs prompt attention.
Technical fault	Red, flashing Technical fault (TF) number is shown on the display	A sequence of 3 beeps, repeated until the alarm is reset	<ul style="list-style-type: none"> Record the technical fault (TF) number⁹. Remove the humidifier from use. Have the humidifier serviced.

⁹ A list of all technical fault (TF) numbers is available in the *HAMILTON-H900 Service Manual*.

5.1 Temporarily silencing an alarm

One component of an alarm is the audible sound. With most alarms you can pause (silence) the alarm sound for two minutes at a time.

To temporarily silence an alarm

- ▶ Press  (Audio pause) on the front of the humidifier.
The audible humidifier alarm is muted for two minutes.
Pressing the key a second time cancels the **Audio pause**.

The **Audio pause** backlight is continuously lit in red while an **Audio pause** is active.

When the time expires and the issue has not yet been resolved, the alarm sounds again.

5.2 Adjusting alarm loudness

WARNING

Be sure to set the auditory alarm loudness above the ambient sound level. Failure to do so can prevent you from hearing and recognizing alarm conditions.

You can set the loudness of the audible alarm.

By default, the loudness is set to 5. Upon turning off the humidifier, the alarm loudness is reset to the default value.

To adjust the alarm loudness

1. If needed, press  (Power/Standby) to put the humidifier into Standby.
2. Press  for 3 seconds.
The display shows the current alarm loudness setting.
3. Adjust the alarm loudness using the **Chamber exit temp** slider (see Figure 3).
The humidifier confirms your selection with a beep at the new loudness setting.
4. Press  to confirm your selection and return to normal operation.

5.3 Troubleshooting alarms

Table 9 lists the humidifier alarms, their graphical representation on the humidifier, a description, and suggested corrective actions.

Table 9. HAMILTON-H900 alarms

Alarm	Alarm icon	Description	Action needed
High priority			
Humidifier tilt		Humidifier is dangerously inclined. The humidifier is at a 10° angle or higher relative to the floor.	<ul style="list-style-type: none"> • Check the mounting of the humidifier. • Check the ventilator trolley. • Operate the humidifier at an angle less than 10° relative to the floor.
Temperature high		The gas temperature at the water chamber or at the Y-piece is above the set value.	<ul style="list-style-type: none"> • Check whether the breathing circuit is covered by the patient's bed covers. • Check whether the breathing circuit or the humidifier chamber is directly exposed to sunlight. • Replace the breathing circuit.
Water level high		The water level in the water chamber is above the maximum level mark.	<ul style="list-style-type: none"> • Empty humidifier chamber to reduce the water level. • Replace the humidifier chamber. • Operate the humidifier at an angle less than 10° relative to the floor.
Technical fault	Technical fault (TF) number is shown on the display	Technical fault due to humidifier malfunction.	<ul style="list-style-type: none"> • Check humidifier operation and all connections. • Replace the humidifier and have it serviced. • If a technical fault number is displayed, make a note of it and provide it when the humidifier is serviced.

Alarm	Alarm icon	Description	Action needed
Medium priority			
Temperature low		The gas temperature at the water chamber exit or at the Y-piece is below the set value.	<ul style="list-style-type: none"> • Wait until the system heats up completely (approx. 30 minutes). • Check whether all settings are correct. • Avoid direct air flow from air conditioning and the like to the humidifier and breathing circuit.
Water level low		The water level in the chamber is below the low level mark.	<ul style="list-style-type: none"> • Check water bottle and refill tubing. • If the water bottle is empty, connect a new water bottle. • Refill or exchange empty humidifier chamber. • Operate the humidifier at an angle less than 10° relative to the floor.
Check water chamber		No chamber or incompatible water chamber is inserted.	Insert a new humidifier chamber and connect the breathing circuit.
Check left or right limb		<p>The display and connection indicators show which limb is faulty.</p> <ul style="list-style-type: none"> • No limb or defective limb connected. • No air flow. • A limb is not properly connected. • The WHITE humidifier expiratory limb is connected to the ventilator <i>To patient</i> inspiratory port. 	<ul style="list-style-type: none"> • Insert or reseat the breathing circuit correctly. • Replace the breathing circuit. • Connect the BLUE humidifier inspiratory limb to the ventilator <i>To patient</i> inspiratory port. <p>When the connection indicator is:</p> <p><i>Green:</i> The limb is inserted correctly.</p> <p><i>Orange:</i> The humidifier is testing the connection.</p> <p><i>Red:</i> The connection is either faulty or nonexistent.</p>

6 Maintenance

Before proceeding, review the safety information in Section 1.

This section provides information about humidifier maintenance procedures and schedule, as well as cleaning and disinfection instructions.

All of the procedures in this section are to be performed by the operator.

For additional maintenance requirements, contact your Hamilton Medical service representative. Any documents referenced in this section are available on the MyHamilton website: <https://www.hamilton-medical.com/MyHamilton>

6.1 Cleaning, disinfection, and sterilization

Humidifier components must be regularly cleaned and disinfected, using the cleaning methods and solutions specific to the individual components.

It is important that you choose the appropriate method and materials when cleaning and disinfecting the humidifier and its components, not only to avoid damaging the equipment, but also to avoid cross-contamination.

Cleaning and disinfection information is presented as follows:

- Table 10 lists the applicable humidifier-related components, and indicates which cleaning and disinfection methods can be used for each one, the frequency with which the component must be cleaned/disinfected, and cleaning and disinfection agents.
- For breathing circuit sets and humidifier water chambers, see the *Instructions for use* or the *Reprocessing Guide*.

When working with the humidifier components, cleaning methods, and cleaning agents, keep the following in mind:

- Do *not* attempt decontamination procedures unless specified by Hamilton Medical.
- While we provide guidelines for agents and concentrations to use, if you have specific questions about the use of a particular cleaning or disinfection agent, contact the manufacturer of the agent.

Table 10. Humidifier cleaning and disinfection methods

Part	Frequency	Cleaning/disinfection method	Cleaning agents
Humidifier exterior including: <ul style="list-style-type: none"> • Housing • Display • Heater plate • Power cables • Mounting systems 	After each patient use or as needed.	Wipe with a damp cloth using a registered and approved cleaning/disinfection solution.	Super Sani-Cloth germicidal wipes CaviWipes

6.2 Preventive maintenance

Preventive maintenance on your humidifier must be performed by Hamilton Medical authorized service personnel according to the instructions in the *HAMILTON-H900 Service Manual*.

7 Configuration

The humidifier is delivered with a default configuration of a variety of settings arranged by patient group.

You can change the following default settings in **Configuration mode**:

- Chamber exit temperature
- Temperature gradient
- Startup operating mode (INV/NIV/HiFlow)
- Expiratory temperature increase On/Off
- Backlight Day/Night brightness settings

7.1 Accessing Configuration mode

You can access **Configuration mode** when the humidifier is in **Standby**.

To access configuration mode

1. If needed, press  to enter Standby.
2. Press and hold both  and . Ensure that you press  first. The **Configuration window** appears after 5 seconds.

The humidifier is now in **Configuration mode**.

7.2 Changing configuration settings

Custom settings are persistent until you change them or restore factory settings.

In **Configuration mode**, the buttons on the humidifier front panel are used to navigate the views and make selections.

Table 11 describes how to page through **Configuration**, and how to change settings.

For an example of how to change a setting, see the procedure below on changing the default settings for **Temperature gradient**.

For a list of configuration settings, see Section 9.6.


Table 11. Configuration mode navigation buttons

Button	Description
	Move forward through the views
	Move backward through the views
	Increase a value
	Decrease a value
	Save changes and exit Configuration mode

Example: Changing Temperature gradient

To change the default Invasive, Adult/Ped, Temperature gradient setting

1. In Configuration mode, press  until you see the view for the Invasive, Adult/Ped, Temperature gradient.

2. Press  to increase the value or press  to decrease the value.
3. When finished, press  to save the changes and exit Configuration mode.

7.3 Restoring the factory default settings

You can restore factory default settings at anytime.

To restore the factory default settings

1. In Configuration mode, repeatedly press  until the Reset all custom defaults window appears.
2. Press .

The factory default settings are restored.

8 Parts and accessories

This section lists the parts and accessories available for the HAMILTON-H900 humidifier.

For additional information, refer to the Hamilton Medical e-catalog on the Hamilton Medical website or contact your Hamilton Medical representative.

Table 12. Humidifier parts and accessories

PN	Description
260161	HAMILTON-BC8022, breathing circuit set, dual limb, heated, with water chamber
260186	HAMILTON-BC4022, breathing circuit set, single limb, heated, with water chamber
260185	HAMILTON-BC8010, breathing circuit set, dual limb, heated, with water chamber and incubator extension
260187	HAMILTON-BC4010, breathing circuit set, single limb, heated, with water chamber and incubator extension
260196	HAMILTON-HC322, water chamber, single use
260197	HAMILTON-HC310, water chamber, single use


9 Specifications

9.1 Physical characteristics

Table 13. Physical characteristics

Dimension	Specifications
Weight	Humidifier (with empty water chamber): 2.5 kg (5.5 lbs)
Dimensions	See Figure 14

Figure 14. HAMILTON-H900 dimensions

9.2 Environmental requirements

Table 14. Environmental requirements

Environment		Specifications
Temperature	Operation:	10°C to 40°C
	Storage:	-20°C to 60°C
		Recommended ambient temperature: 18°C to 26°C
Altitude		Up to 4,000 m (13,123 ft)
Atmospheric pressure	Operation and storage:	61 kPa to 106 kPa
Relative humidity	Operation:	30% to 95%, noncondensing
	Storage:	10% to 95%, noncondensing
Water protection		IP21
Gas input temperature		18°C to 31°C (recommended)

9.3 Electrical specifications

Table 15. Electrical specifications

Part number	Input power	Power consumption
950001 (230 V)	220 – 240 V 50/60 Hz	283 VA
950004 (115 V)	110 – 127 V 50/60 Hz	293 VA
950008 (100 V)	100 V 50/60 Hz	268 VA
Potential equalization	Terminal for the connection of a potential equalization conductor according to DIN 42801	

9.4 Control settings

Table 16. Control settings and ranges

Parameter or setting (unit)	Mode	Range	Default	Resolution
Chamber exit temperature (°C)	INV	35 to 41	37.0	0.5
	NIV	30 to 35	31.0	0.5
	HiFlow	33 to 37	35.0	2
Temperature gradient (°C)	INV	-2 to 3	Adult/Ped: 2 Neo: 3	0.5
	NIV	-2 to 3	Adult/Ped: 2 Neo: 3	0.5
	HiFlow	--	2	--
Resulting max. airway temperature (Y-piece) ¹⁰ (°C)	INV	33 to 42	--	--
	NIV	28 to 38	--	--
	HiFlow	35 to 39	--	--
Expiratory temperature increase	All modes	On, Off	On	--
Brightness setting for display	All modes	1 to 5	Day: 5 Night: 2	1

¹⁰ Airway temperature is limited to 42°C by the humidifier.

9.5 Monitored parameters

Table 17. Monitored parameters, ranges, and accuracy

Parameter (units)	Range	Accuracy
Chamber exit temp (°C)	10 to 60	10 to 30: ± 1 30 to 41: ± 0.5 41 to 60: ± 1
Y-piece temperature (°C)	28 to 43	± 0.5

9.6 Configuration

Table 18. Configuration specification

Parameter	Mode	Configuration range	Default setting
Chamber exit temp (°C)	INV	35 to 39	37.0
	NIV	30 to 35	31.0
	HiFlow	33, 35, 37	35.0
Temperature gradient (°C)	INV	0 to 3	Adult/Ped: 2 Neo: 3
	NIV	0 to 3	Adult/Ped: 2 Neo: 3
	HiFlow	0 to 3	2
Startup mode	--	Invasive, Noninvasive, HiFlow	Invasive
Brightness setting for display	All modes	1 to 5	Day: 5 Night: 2

9.7 Technical performance data

Table 19. Technical performance data

Description	Specification		
Flow rates ¹¹	Invasive:	Up to 60 l/min	
	Noninvasive:	Up to 120 l/min	
	HiFlow:	Up to 120 l/min	
Warm-up time	Less than 30 min		
Humidity	At an ambient temperature of 18°C to 26°C:		
	Invasive:	Temperature setting of 37 to 41°C	Minimum humidity: 33 mg H ₂ O/l
	Noninvasive:	Temperature setting of 31 to 35°C	Minimum humidity: 12 mg H ₂ O/l
	HiFlow:	Flow rate ≤ 60 l/min	Minimum humidity: 33 mg H ₂ O/l
Flow rate > 60 l/min		Minimum humidity: 12 mg H ₂ O/l	
Standby heating specifications	Heated breathing circuits:	Y-piece limited to 30°C.	
Display	3 in / 64 x 128 pixels, inverted dot matrix display (backlit)		
Audio pause	120 seconds		
Alarm loudness ¹²	Range = 1 to 8. Default = 5.		

¹¹ For minimal flow rates, see the breathing circuit *Instructions for use*.

¹² Volume at 1 meter distance from humidifier. A setting of 1 = 50 dB(A), 5 (default) = 60 dB(A), and 8 = 65 dB(A), with an accuracy of ±6 dB(A).

9.8 Essential performance

The temperatures applied to the respiratory gas and breathing circuit must be maintained within the specified settings, and must be monitored.

Should the temperature fall outside the specified limits, the humidifier must detect this fact and inform the user through an alarm.

9.9 Functional description of humidification system

The HAMILTON-H900 humidifier is designed to heat and humidify respiratory gases during mechanical ventilation. The respiratory gas is heated and humidified as it passes through a chamber that is partially filled with heated water.

The humidifier has two sensor-controlled heating systems:

- *Heater plate.* Heats the water in the water chamber.
- *Integrated breathing circuit heating system.* Heats the breathing limbs to prevent condensation.

The respiratory gas temperature is monitored at the chamber exit and in the breathing limb at the patient end.

9.10 Symbols used on device labels and packaging

Symbol	Description
	Follow the <i>Instructions for use</i>
	Reference number
	Serial number
	Quantity
	Manufacturer
	Date of manufacture
	Type BF applied part (classification of medical electrical equipment, type BF, as specified by IEC 60601-1)
	Dispose according to Council Directive 2002/96/EC or WEEE (Waste Electrical and Electronic Equipment)
	The TÜV NRTL mark with the indicators "C" and "US" means that the product complies with Canadian requirements and the requirements of US authorities for safety
	Potential equalization

Symbol	Description
	Warning: Hot surface. The heater plate and the bottom of the chamber can reach a temperature of over 85°C.
CE0197	CE Marking of Conformity, seal of approval guaranteeing that the device is in conformance with the Council Directive 93/42/EEC concerning medical devices
IP21	Protected against dripping water and solid particles larger than 12.5 mm
max	Filling level mark on water chamber
IOIOI	RS-232 serial interface for communication with Hamilton Medical ventilators
	MR unsafe The HAMILTON-H900 humidifier poses unacceptable risks to patient, medical staff, or other persons within the MR environment.
	Chinese RoHS
	Temperature probe Marks the position of the temperature probe on the breathing circuit.

9.11 Standards and approvals

The humidifier system meets relevant parts of the following standards, listed in Table 20.

Table 20. Standards and approvals, valid versions

IEC 60601-1:2005/A1:2012
IEC 60601-1-2:2014
IEC 60601-1-6:2010/A1:2013
IEC 60601-1-8:2006/A1:2012
ISO 80601-2-74:2017
EN ISO 5356-1:2015
EN ISO 5367:2014
IEC 62304:2006/A1:2015
IEC 62366-1:2015
ISO 13485:2016

9.12 Disposal

The device must be disposed of according to your institution's protocols and Directive 2002/96/EC.

Breathing circuits and humidifier chambers must be disposed of according to the *Instructions for use* supplied with the breathing circuit sets.

9.13 Warranty

LIMITED WARRANTY

THE WARRANTY DESCRIBED IN THIS AGREEMENT IS IN LIEU OF ANY AND ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. HOWEVER, IMPLIED WARRANTIES ARE NOT DISCLAIMED DURING THE PERIOD OF THIS LIMITED WARRANTY.

Hamilton Medical guarantees its products to be shipped free from defects in material and workmanship. The warranty does not include disposable items. Disposable items and consumable products are considered to be of single use or of limited use only and must be replaced regularly as required for proper operation of the product following the Operator's Manual.

Hamilton Medical shall have no obligations nor liabilities in connection with the product other than what is specified herein, including without limitation, obligations, and/ or liabilities for alleged negligence, or for strict liability. In no event shall the company be liable for incidental or consequential damages, either direct or contingent.

This Limited Warranty shall be void and not apply:

1. If the product has not been installed and connected by an authorized local representative of Hamilton Medical in accordance with the instructions furnished by Hamilton Medical and by a Hamilton Medical representative.
2. If no evidence is present that the occurrence of damage/repair happened within the certified warranty period.
3. If the serial number has been altered, effaced, or removed, and there is no bill of sale or evidence to verify the product's purchase date.
4. If the defects arise from misuse, negligence, or accidents or from repair, adjustment, modification, or replacement made outside Hamilton Medical's factories or other than an authorized service center or authorized service representative.
5. If the product has been modified, or in any nature altered without prior written authorization from Hamilton Medical.
6. If the product is or has been used in any way that is not specified under "Intended Use".
7. If the product has been used by anyone but properly trained personnel under the supervision of a physician.

Replacements and/or repairs furnished under this Limited Warranty do not carry a new warranty, but carry only the unexpired portion of the original Limited Warranty. The warranty of repaired and/or replaced components does not exceed the Limited Warranty of the device.

To obtain service under this Limited Warranty, claimant must promptly notify the country's sales partner of Hamilton Medical regarding the nature of the problem, serial number, and the date of purchase of the Product.

Except as stated above, Hamilton Medical shall not be liable for any damages, claims, or liabilities including, but not limited to, personal bodily injury, or incidental, consequential, or special damages. Nor will Hamilton Medical be liable for any damages, claims, or liabilities including, but not limited to, personal bodily injury, or incidental, consequential, or special damages resulting from misuse of the device or failure to comply with any of the provisions made in this manual.

9.13.1 Miscellaneous

The general terms and conditions of Hamilton Medical shall be applicable. This agreement shall be governed by, and construed in accordance with, the laws of Switzerland and may be enforced by either party under the jurisdiction of the court of Chur, Switzerland.

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medin Medical Innovations GmbH