REF

160021

English | 10203811/00

Software version 2.0.2

2026-03-31

Store this Addendum together with your ventilator Operator's Manual.

This Addendum provides the following information.

For details about...	See...
Corrections/additions to the HAMILTON-C6 Operator's Manual Applies to software version 2.0.2 and higher	Section 1
In2Flow nasal cannula discontinuation Applies to all software versions	Section 2

1 Corrections/additions to the HAMILTON-C6 Operator's Manual

1.1 Safety message updates

The following safety message regarding defibrillator usage for applied parts has been updated.

 **CAUTION**

All applied parts are NOT protected against reanimation with a defibrillator. Disconnect the CO2 sensor before using a defibrillator on the patient.

The following safety message regarding *distributed alarm system* (DAS)-related alarms has been added.

 **WARNING**

When an Audio pause is active, the HAMILTON-H900 medium- and high-priority alarms still generate an audible alarm.

The following safety message regarding the ventilator exhaust port has been added.

 **WARNING**

Do *not* connect any component or device to the exhaust port of the expiratory valve unless authorized by Hamilton Medical.

1.2 Changes to CPR ventilation controls in software version 2.0.2

During CPR ventilation when (S)CMV mode is selected and Peak flow timing philosophy is active, the CPR Configuration window now displays a TI control instead of Peak flow.

The rows below were removed from the following tables in the *Operator's Manual*:

Table 1. CPR ventilation modes, controls, and default settings (Section 10.9.1, Table 10-6, in the ventilator *Operator's Manual*)

Parameter	Adult/ Ped. default	Neonatal default
Peak flow (l/min) <i>(S)CMV only</i>	60	--

Table 2. CPR default settings (Section 16.9, Table 16-17, in the ventilator *Operator's Manual*)

Parameter	(S)CMV (Adult/ Ped.)	APVcmv (Neo)
Peak flow (l/min)	60	--

1.3 Cybersecurity safety update

The following cybersecurity safety section has been added to the safety information chapter.

Cybersecurity is the practice of protecting computer systems, networks, devices, and data from digital attacks, damage, or unauthorized access.

This section provides cybersecurity-related information for trained personnel using the device. Note that it does not cover all possible device vulnerabilities or address all specific configurations and use cases.

By applying the secure configuration guidelines (Section 1.3.3), device security improves, and the possibility of cybersecurity-related disruptions to essential medical functions are reduced.

The organization using this device is responsible for ensuring its secure deployment, configuration, maintenance, and decommissioning.

Table 3. Cybersecurity overview

For details about ...	See ...
Roles and responsibilities	Section 1.3.1
Summary of risks	Section 1.3.2
Configuration guidelines	Section 1.3.3
Monitoring, logging, and detection	Section 1.3.4
Maintenance	Section 1.3.5
Safety ventilation	Section 1.3.6
Cybersecurity incident response	Section 1.3.7

1.3.1 Roles and responsibilities

Managing cybersecurity involves different tasks based on the different user groups working with the ventilator.

Clinical ventilator users are responsible for:

- Following the secure ventilator operation guidelines (Section 1.3.3)
- Contacting the Hamilton Medical authorized service representative for all maintenance requirements, including software

Hamilton Medical authorized service personnel are responsible for:

- Installing the latest software version available in your region, if requested
- Providing the ventilator network address and other software related information, if requested

The ventilator *Service Manual* (PN 627038) provides this information for service personnel.

1.3.2 Security risk approach

A detailed threat model was developed for the ventilator to perform a thorough risk assessment. Where necessary, appropriate risk controls have been applied.

The HAMILTON-C6 is protected against reasonable cybersecurity threats in a clinical environment. Due to the restricted intended location of the ventilator, and availability of supporting documentation, it is very likely that the ventilator can operate in a secure environment with minimal risk of cyber attacks.

All threats from third-party components are constantly monitored.

1.3.3 Secure configuration guidelines

Contact your Hamilton Medical authorized service representative to ensure secure device installation and configuration.

For details about connecting to external devices, see Section 4.9, *Connecting to external devices* in the ventilator *Operator's Manual*.

To securely operate the ventilator:

- Verify the ventilator device functionality before connecting the device to the patient.
- Only allow authorized personnel to access and operate the ventilator.
- Only connect compatible and Hamilton Medical-approved devices and components to the ventilator.¹
- Do *not* access internal ventilator components. Only Hamilton Medical authorized service personnel may service the ventilator.
- Adhere to the physical security hardening guidelines.

1.3.3.1 Physical security hardening guidelines

Only operate the ventilator in dedicated areas that meet the physical security requirements (for example, regarding access control and monitoring) defined by your facility.

1.3.4 Monitoring, logging, and detection

Once the ventilator is turned on, the device starts collecting logging data about clinically relevant ventilator activities, including ventilator start and shutdown times, alarms, technical notes, setting changes, calibrations, maneuvers, and special functions. For details about the logging data, see Section 10.13, *About the Event log* in the ventilator *Operator's Manual*.

This Operator's Manual provides detailed information about ventilator operation and capabilities.

When conditions generate an alarm, acoustic and visual notification is provided, as described in Chapter 9, *Responding to alarms* in the ventilator *Operator's Manual*.

1.3.5 Maintenance

Maintenance tasks described in this manual clearly indicate whether they are performed by the ventilator user or by a Hamilton Medical approved service technician. For detailed maintenance information, see Chapter 13, *Maintenance* in the ventilator *Operator's Manual*.

Hamilton Medical authorized service personnel perform security-related maintenance tasks. Preventive maintenance is performed as described in Section 13.4, Table 13-4, in the ventilator *Operator's Manual*.

¹ See the Technical Specifications for your device, as well as the Hamilton Medical e-catalog for compatible devices.

1.3.6 Safety ventilation

In the event of damage, error, or technical fault that affects ventilator function, the HAMILTON-C6 automatically modifies operation to ensure essential function.

These special conditions are described in detail in Section 7.9, *Special conditions* in the ventilator *Operator's Manual*. Note that some manual intervention is required from the user.

1.3.7 Cybersecurity incident response

This section describes what to do when a cybersecurity vulnerability occurs.

1.3.7.1 Monitoring and security scanning

Hamilton Medical constantly monitors evolving threats and employs various security scanning methods and technologies. In addition, Hamilton Medical has devices tested independently by third parties.

1.3.7.2 Security advisories

Hamilton Medical releases security advisories to the public in case of potential vulnerabilities. For details, go to:

<https://www.hamilton-medical.com/Services/Cybersecurity.html>

To protect patients from harm, the following information is shared in case of a vulnerability occurrence:

- Information about products that are affected by a vulnerability and how they are affected
- Information about vulnerabilities of components that are used in other products
- Information about IT equipment that may impact the security of medical devices
- Information about attacks, potential attacks, and availability of exploit code
- Confirmation of incidents
- Availability of patches and other security mitigations
- If necessary, additional instructions on how to mitigate the issue as an interim measure

1.3.7.3 Incident reporting and emergency contact information

 WARNING

Do *NOT* use the ventilator if any cybersecurity issue occurs and provide alternative respiratory support until the issue is resolved.

Hamilton Medical employs trained security personnel to react to product security incidents.

To report a security issue with Hamilton Medical, go to:
<https://www.hamilton-medical.com/Services/Cybersecurity.html>

In case of a security-related emergency, contact your country's Hamilton Medical representative.

1.4 Default control changes in software version 2.0.2

The default Flow Pattern for CPR/(S)CMV is now Square.

Table 4. Control settings, ranges, and accuracy (Section 16.6, Table 16-10, in the ventilator *Operator's Manual*)

Parameter or setting (unit)	Default
Flow Pattern	<i>Adult/Ped.</i> 50% decelerating, CPR/(S)CMV: Square <i>Neo:</i> --

2 In2Flow nasal cannula discontinuation

The following In2Flow nasal cannulas have been discontinued:

- In2Flow nasal cannula (Large)
(PN 10076604)
- In2Flow nasal cannula (Medium)
(PN 10076605)
- In2Flow nasal cannula (Small)
(PN 10076606)

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