
Instructions for use	Expiratory valve set	3
Gebrauchsanweisung	Expirationsventil	7
Instructions d'utilisation	Système de valves expiratoires	11
Instrucciones de uso	Conjunto de válvula espiratoria	15
Istruzioni per l'uso	Set valvola espiratoria	19
Instruções de uso	Conjunto de válvula expiratória	23
使用说明	呼气阀套件	27
使用説明書	呼气弁セット	30
Инструкции по эксплуатации	Комплект клапана выдоха	33
Kullanım talimatları	Ekspiratuvar valf seti	37

HAMILTON-C2/C3/C6
(REF 160176)

HAMILTON-C1/T1/MR1
(REF 161189)

HAMILTON-C1/T1/MR1
(REF 161186, 161186X)

REF 160176, 161186, 161186X, 161189


CE


HAMILTON
MEDICAL


HAMILTON-C2/C3/C6
(REF 160176)

HAMILTON-C1/T1/MR1
(REF 161189)

HAMILTON-C1/T1/MR1
(REF 161186, 161186X)

- 1 Protective cap – Remove before use | Schutzkappe – Vor dem Einsatz entfernen | Capuchon de protection – Retirer avant utilisation | Tapa de protección: Retirar antes de usar | Cappuccio protettivo (da rimuovere prima dell'uso) | Tampa de proteção – Remover antes de usar | 保护盖 – 在使用之前取下 | 保護キャップ – 使用前に取り外してください | Защитная крышка (перед использованием снимается) | Korgucucu kapak Kullanmadan önce çıkarın

- Expiratory valve set components (2-4)
Komponenten des Expirationsventils (2-4)
Composants du système de valves expiratoires (2-4)
Componentes del conjunto de válvula espiratoria (2-4)
Componenti del set valvola espiratoria (2-4)
Componentes do conjunto de válvula expiratória (2-4)
呼气阀套件组件 (2-4)
呼吸弁セットコンポーネント (2-4)
Компоненты комплекта клапана выдоха (2-4)
Ekspiratuvar valf seti bileşenleri (2-4)

- 2 Expiratory valve membrane | Expirationsventilmembran | Membrane de valve expiratoire | Membrana de la válvula espiratoria | Membrana della valvola espiratoria | Membrana da válvula expiratória | 呼气阀膜 | 呼吸弁ダイヤフラム | Мембрана клапана выдоха | Ekspiratuvar valf membranı
- 3 Expiratory valve housing | Expirationsventilkörper | Boîtier de valve expiratoire | Carcasa de la válvula espiratoria | Calotta della valvola espiratoria | Invólucro da válvula expiratória | 呼气阀外壳 | 呼吸弁ハウジング | Корпус клапана выдоха | Ekspiratuvar valf yuvası
- 4 Expiratory valve locking ring | Sicherungsventilkörper | Bague de verrouillage de la valve expiratoire | Anillo de bloqueo de la válvula espiratoria | Anello di bloccaggio della valvola espiratoria | Anel de bloqueio da válvula expiratória | 呼气阀锁环 | 呼吸弁ロックリング | Уплотнительное кольцо клапана выдоха | Ekspiratuvar valf kilit halkası

REF	160176	Expiratory valve set HAMILTON-C2/C3/C6
	161186, 161186X	Expiratory valve set HAMILTON-C1/T1/MR1
	161189	Expiratory valve set HAMILTON-C1/T1/MR1

Conventions

WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

CAUTION

Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.

NOTICE

Emphasizes information of particular importance.

1 Overview

The PN 160176 expiratory valve set can be used with the HAMILTON-C2/C3/C6 ventilators.

The PN 161186, 161186X, and 161189 expiratory valve sets can be used with the HAMILTON-C1/T1/MR1 ventilators.

1.1 Intended use

The intended uses for the expiratory valve sets are as follows.

REF 160176

The Hamilton Medical single-use expiratory valve set for the HAMILTON-C2/C3/C6 is an accessory for the exhalation gas path. Its main function is to control expiratory gas flow and maintain user-set pressure levels. The Hamilton Medical single-use expiratory valve set is not intended for use with equipment other than Hamilton Medical ventilators.

The single-use expiratory valve set is intended for use by qualified, trained personnel.

REF 161186, 161186X

The Hamilton Medical single-use adult/pediatric expiratory valve set for HAMILTON-C1/T1/MR1 is an accessory for the exhalation gas path. Its main function is to control expiratory gas flow and maintain user-set pressure levels. The Hamilton Medical single-use expiratory valve set is not intended for use with equipment other than Hamilton Medical ventilators.

The single-use expiratory valve set is intended for use by qualified, trained personnel.

The Hamilton Medical single-use neonatal expiratory valve set for HAMILTON-C1/T1/MR1 is an accessory for the exhalation gas path. Its main function is to control expiratory gas flow and maintain user-set pressure levels. The Hamilton Medical single-use expiratory valve set is not intended for use with equipment other than Hamilton Medical ventilators.

The single-use expiratory valve set is intended for use by qualified, trained personnel.

2 Safety information

WARNING

- Reprocessing of Hamilton Medical single-use products can affect the product properties and may cause injury to the patient. For example, a change to the surface structure during reprocessing may lead to a change in the tear strength or cause actual cracking. Furthermore, an altered surface structure may result in, for example, a microbial aggregation of spores, allergens, and pyrogens, or cause an increase in the number of particles released as a result of chemical changes in the material properties.

- Hamilton Medical does not assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user. Any liability on the part of Hamilton Medical for any damages resulting from the use of reprocessed or reused single-use parts, or use of any non-Hamilton Medical accessories and consumables, is excluded.
- Follow hospital infection control procedures for all tests and use.
- Remove the protective cover(s) before use.
- Discard the product if it fails two subsequent preoperational checks or if there is any sign of damage.

CAUTION

- *(USA only) Federal law restricts this device to sale by or on the order of a physician.*
- *Handle used products as contaminated goods according to local laws and regulations or hospital internal procedures.*
- *Always carry out all preoperational checks before connecting the patient to prevent possible patient injury. Make sure another source of ventilatory support is available.*

NOTICE

- Before use, read the appropriate ventilator *Operator's Manual*.
- Replace the product in accordance with hospital infection control procedures, or depending on the patient's secretions and nebulization of medication. Consider using a filter at the expiratory port during nebulization.
- When adding parts to the Hamilton Medical breathing circuit configurations, do not exceed the inspiratory and expiratory resistance values of the ventilator breathing system as specified in the ventilator *Operator's Manual*.
- Ensure the expiratory valve set is in its intended position in the ventilator.
- Ensure the membrane is aligned with the expiratory valve body, and the metal faces up.
- Any incident with the device leading to serious patient injury, death, or a potential threat to public health must be reported to the manufacturer and the relevant authorities.

3 Installing the expiratory valve set

To install the expiratory valve set

1. Before use, remove the protective cap (1 on page 2).
2. Ensure the membrane is properly aligned with the expiratory valve housing and the metal plate faces up (A in Figures 1 and 2).
3. Position the expiratory valve set in the expiratory port and twist the locking ring clockwise until it locks into place (B in Figures 1 and 2).

Figure 1. Installing the expiratory valve set HAMILTON-C2/C3/C6



Figure 2. Installing the expiratory valve set HAMILTON-C1/T1/MR1

4 Preoperational check

Perform the preoperational check (Leak test (Tightness)) before patient use. Follow the instructions provided in the corresponding ventilator *Operator's Manual*.

5 Delivery

The components are delivered clean and ready for clinical use.

6 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

7 Technical specifications

HAMILTON-C2/C3/C6 (160176)

Maximum operating duration	28 days
Operating temperature	15°C to 35°C (59°F to 95°F)
Transient operating temperature (for ≤ 20 minutes)	5°C to 40°C (41°F to 104°F)
Transport temperature (for ≤ 4 weeks)	-20°C to 60°C (-4°F to 140°F)
Storage temperature	15°C to 35°C (59°F to 95°F)
Storage, transport, and operating, humidity	< 95% relative humidity, noncondensing
Conical connectors	OD22/ID15, according to ISO 5356-1:2015

HAMILTON-C1/T1/MR1
(161186, 161186X, 161189)

Maximum operating duration	28 days
Operating temperature 161186, 161186X, 161189	15°C to 35°C (59°F to 95°F)
Transient operating temperature Adult/ Pediatric (for ≤ 20 minutes) 161186, 161186X	-15°C to 50°C (5°F to 122°F)
Transient operating temperature Neonatal (for ≤ 20 minutes) 161189	-15°C to 40°C (5°F to 104°F)
Transport temperature (for ≤ 4 weeks) 161186, 161186X, 161189	-20°C to 60°C (-4°F to 140°F)
Storage temperature 161186, 161186X, 161189	15°C to 35°C (59°F to 95°F)
Storage, transport, and operating, humidity 161186, 161186X, 161189	< 95% relative humidity, noncondensing
Conical connectors 161186, 161186X, 161189	OD22/ID15, according to ISO 5356-1:2015

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For devices manufactured in Switzerland

Manufacturer

Hamilton Medical AG
Via Crusch 8, 7402 Bonaduz, Switzerland
☎ +41 (0) 58 610 10 20
info@hamilton-medical.com
www.hamilton-medical.com

medin Medical Innovations GmbH
Adam-Geisler-Straße 1
DE-82140 Olching
Germany

For devices manufactured and sold in the USA
Printed in the USA

Manufactured for

Hamilton Medical AG
Via Crusch 8, 7402 Bonaduz, Switzerland
☎ +41 (0) 58 610 10 20
info@hamilton-medical.com
www.hamilton-medical.com

Distributor in USA

Hamilton Medical, Inc.
201 Edison Way, Unit A
Reno, NV 89502-2305, US
☎ (800) 426-6331 (toll free), (775) 858-3200
info@hamiltonmedical.com

10103753/01
2024-05-31

Specifications are subject to change without notice. For all proprietary trademarks and third-party trademarks used by Hamilton Medical AG, see www.hamilton-medical.com/trademarks.

© 2024 Hamilton Medical AG. All rights reserved.