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

30
PetCO₂
mmHg

255
V_{Daw}
ml

8.2
V'_{alv}
l/min

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- Anderson** JT, Owings JT, Goodnight JE. Bedside noninvasive detection of acute pulmonary embolism in critically ill surgical patients. Arch Surg 1999;134(8):869–874; discussion 874–875.
- Aström** E, Niklason L, Drefeldt B, Bajc M, Jonson B. Partitioning of dead space – a method and reference values in the awake human. Eur Respir J. 2000 Oct; 16(4):659-664.
- Blanch** L, Romero PV, Lucangelo U. Volumetric capnography in the mechanically ventilated patient. Minerva Anesthesiol. 2006 Jun;72(6):577-85.
- Erikson**, L, Wollmer, P, Olsson, CG, et al. Diagnosis of pulmonary embolism based upon alveolar dead space analysis. Chest1989;96,357-362.
- Fletcher** R. The single breath test for carbon dioxide [dissertation]. Lund, Sweden: University of Lund, 1980. 2nd edition revised and reprinted, Solna, Sweden: Siemens Elema, 1986.
- Kallet** RH, Daniel BM, Garcia O, Matthay MA. Accuracy of physiologic dead space measurements in patients with acute respiratory distress syndrome using volumetric capnography: comparison with the metabolic monitor method. Respir Care. 2005 Apr;50(4):462-7.
- Kiiski**, Ritva, and Jukka Takala. „Hypermetabolism and efficiency of CO2 removal in acute respiratory failure.“ CHEST Journal 105.4 (1994): 1198-1203.
- Kumar** AY, Bhavani-Shankar K, Moseley HS, Delph Y. Inspiratory valve malfunction in a circle system: pitfalls in capnography. Can J Anaesth 1992;39(9):997–999.
- Nuckton** TJ, Alonso JA, Kallet RH, Daniel BM, Pittet JF, Eisner MD, Matthay MA. Pulmonary dead-space fraction as a risk factor for death in the acute respiratory distress syndrome. N Engl J Med. 2002 Apr 25; 346(17):1281-1286.
- Olsson** K, Jonson B, Olsson CG, Wollmer P. Diagnosis of pulmonary embolism by measurement of alveolar dead space. J Intern Med. 1998 Sep;244(3):199-207.
- Pyles** ST, Berman LS, Modell JH. Expiratory valve dysfunction in a semiclosed circle anesthesia circuit: verification by analysis of carbon dioxide waveform. Anesth Analg 1984;63(5):536–537.
- Radford** EP. Ventilation standards for use in artificial respiration. N Engl J Med 1954; 251:877-883.
- Rodger** MA, Jones G, Rasuli P, Raymond F, Djunaedi H, Bredeson CN, Wells PS. Steady-state end-tidal alveolar dead space fraction and D-dimer: bedside tests to exclude pulmonary embolism. Chest 2001;120(1):115–119.
- Yaron** M, Padyk P, Hutsinpiiler M, Cairns CB. Utility of the expiratory capnogram in the assessment of bronchospasm. Ann Emerg Med. 1996 Oct;28(4):403-7.
- Weissman** C, Kemper M, Elwyn DH, Askanazi J, Hyman AI, Kinney JM. The energy expenditure of the mechanically ventilated critically ill patient. An analysis. Chest. 1986 Feb; 89(2):254-259.
- Wolff** G, Brunner JX, Grädel E. Gas exchange during mechanical ventilation and spontaneous breathing. Intermittent mandatory ventilation after open heart surgery. Chest. 1986 Jul; 90(1):11-17
- Wolff** G, Brunner JX, Weibel W, et al. Anatomical and series dead space volume: concept and measurement in clinical practice. Appl Cardiopul Pathophysiol 1989; 2:299-307.

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作者： Karjaghli Munir
Matthias Himmelstoss

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Handwriting practice lines consisting of four sets of horizontal lines. Each set includes a solid top line, a dashed midline, and a solid bottom line, providing a guide for letter height and placement.

Handwriting practice lines consisting of four sets of horizontal lines. Each set includes a solid top line, a dashed midline, and a solid bottom line, providing a guide for letter height and placement.

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ELO20160411N.02 此处所示产品并非向一般公众出售。此处所提供的信息仅适用于医疗专业人士。务必阅读标签并遵循产品的使用说明。产品规格如有变更、恕不另行通知。某些功能是选配件。并非所有市场均提供所有功能。所有图像仅作图示说明之用，可能未准确展示产品或其用途。有关全部专有商标 (®)和 Hamilton Medical 哈美顿医疗公司使用的第三方商标，请参阅 www.hamilton-medical.com/trademarks。© 2023 Hamilton Medical 哈美顿医疗公司。版权所有。
Hamilton Medical AG