
Instructions for use	HAMILTON-BC8022 Breathing circuit set	4
Gebrauchsanweisung	HAMILTON-BC8022, Beatmungsschlauchset	9
Instructions d'utilisation	Kit de circuit respiratoire HAMILTON-BC8022	15
Instrucciones de uso	HAMILTON-BC8022, equipo respiratorio	20
Istruzioni per l'uso	HAMILTON-BC8022, set circuito paziente	26
Instruções de uso	HAMILTON-BC8022, kit de respiração	31
使用説明書	HAMILTON-BC8022呼吸回路セット	36
使用说明	HAMILTON-BC8022 呼吸装置	41
Инструкция по эксплуатации	комплекта дыхательного контура HAMILTON-BC8022	45
Kullanım Talimatları	HAMILTON-BC8022, Solunum devresi seti	51


REF 260161

CE 0197 **EAC**

HAMILTON
MEDICAL

- | | | |
|--|--|---|
| <p>1 Ventilator inspiratory limb
Inspirationsschenkel für das
Beatmungsgerät Branche
inspiratoire du ventilateur
Rama inspiratoria del
respirador Branca inspiratoria
del ventilatore Alça inspiratória
do respirador 人工呼吸器側の吸
気回路 呼吸机吸気枝
Патрубок вдоха со стороны
аппарата ИВЛ Ventilator
inspiratuvar parçası</p> | <p>4 Y-piece (OD22/OD22 / ID15/OD22) Y-Stück
(AD22/AD22 – ID15/AD22) Pièce en Y (DE 22/
DE 22/DI 15/DE 22) Pieza en Y (DE22/DE22/
DI15/DE22) Raccordo a Y (D.E. 22/D.E. 22;
D.I. 15/D.E. 22) Peça em "Y" (DE22/DE22 /
DI15/DE22) Yピース (外径22/外径22、内径
15/外径22) Y形管 (外径22/外径22 / 内径
15/外径22) Y-образный коннектор
(HД: 22 / HД: 22 / ВД: 15 / HД: 22) Y parçası
(DÇ22/DÇ22 / İÇ15/DÇ22)</p> | <p>7 Elbow adapter, 90° (ID22/OD22)
Ellenbogenadapter, 90° (ID22/AD22)
Adaptateur coudé, 90° (DI 22/DE 22)
Adaptador acodado, 90° (DI22/DE22)
Adattatore a gomito, 90° (D.I. 22/
D.E. 22) Adaptador de cotovelo, 90°
(DI22/DE22) エルボアダプタ、90° (内
径22/外径22) 弯接头, 90° (内径 22/
外径 22) Угловой адаптер, 90°
(ВД: 22 / HД: 22) Dirsek adaptörü, 90°
(İÇ22/DÇ22)</p> |
| <p>2 Humidification chamber
Befeuchterkammer Réservoir
d'humidification Câmara
de agua Camera di
umidificazione Câmara de
umidificação 加湿チャンバ 湿
化器水箱 Камера увлажнения
 Nemlendirme haznesi</p> | <p>5 Patient expiratory limb Expirationsschenkel
für den Patienten Branche expiratoire
patient Rama espiratoria del paciente
Branca espiratoria del circuito paziente
Alça expiratória do paciente 患者側の呼気回
路 病人呼気枝 Патрубок выдоха со стороны
пациента Hasta espiratuvar parçası</p> | |
| <p>3 Patient inspiratory limb
Inspirationsschenkel für den
Patienten Branche inspiratoire
patient Rama inspiratoria del
paciente Branca inspiratoria del
circuito paziente Alça
inspiratória do paciente 患者側
の吸気回路 病人吸気枝
Патрубок вдоха со стороны
пациента Hastaya giden
inspiratuvar parça</p> | <p>6 Water feed line with spike Wasserzufuhr-
leitung mit Anstechdorn Ligne d'alimentation
en eau avec pointe Vía de suministro de agua
con punta Linea di alimentazione dell'acqua
con perforatore Linha de alimentação da água
com grampo スパイク付き給水チューブ 帶
有尖狀灌注器的供水管線 Линия подачи воды
со штекерным соединителем Sivri uçlu su
besleme hattı</p> | |

Conventions

WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

CAUTION

Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.

NOTICE

Emphasizes information of particular importance.

1 Overview

The HAMILTON-BC8022 breathing circuit set (PN 260161) can only be used with the HAMILTON-H900 humidifier. For an overview of the breathing circuit set and its components, see page 2.

1.1 Intended use

The HAMILTON-BC8022 breathing set is intended to be used together with compatible Hamilton Medical respiratory gas humidifiers during invasive and noninvasive mechanical ventilation of adult and pediatric patients > 10 kg.

2 Safety information

WARNING


- Reprocessing of Hamilton Medical single-use products can affect the product properties and may cause injury to the patient. For example, a change to the surface structure during reprocessing may lead to a change in the tear strength or cause

actual cracking. Furthermore, an altered surface structure may result in, for example, a microbial aggregation of spores, allergens, and pyrogens, or cause an increase in the number of particles released as a result of chemical changes in the material properties.

- Hamilton Medical does *not* assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user. Any liability on the part of Hamilton Medical for any damages resulting from the use of reprocessed or reused single-use parts, or use of any non-Hamilton Medical accessories and consumables, is excluded.

- Adding attachments or other parts/assemblies to a breathing system can change the pressure gradient across the ventilator, which can adversely affect ventilator performance and patient safety.
- Discard the product if it fails two subsequent preoperational checks or if there is any sign of damage.
- Follow hospital infection control procedures for all tests and use.
- Connect only a single breathing circuit at a time. Do *not* connect any circuit as an extension to another one. Doing so can seriously harm or kill the patient.
- Before using the humidifier on the patient, verify that the breathing circuit is correctly connected to the ventilator as follows:
 - The BLUE inspiratory limb is connected to the *To patient* inspiratory port.
 - The WHITE expiratory limb is connected to the *From patient* expiratory port.

Failure to connect the limbs to the proper ventilator ports can result in delivery of overheated gas to the patient.

- Use caution when positioning the breathing circuit. Place the breathing circuit in such a way that condensation cannot reach the patient.
- Heated limbs must *not* be placed directly on the patient's skin. Attach breathing circuits or tubing holders appropriately to avoid mechanical forces on the ET tube.
- Route the breathing circuit without tension and without any kinks from the ventilator or humidifier to the patient, and protect against unintended influences.
- The humidifier must always be positioned below patient level. Do *not* operate the humidifier at an angle in excess of 10°.
- The breathing circuit set must *not* be covered by any objects, such as sheets, towels, and so on.

- Do *not* touch the hot plate or the bottom of the chamber. The surfaces can reach a temperature of over 85°C. These hot surfaces radiate heat.
- Do *not* use heated breathing circuits without existing gas flow. If the gas supply is interrupted, the humidifier must be turned off or be in Standby mode.
- Only fill the humidification chamber with sterile, demineralized water that meets the hospital's hygiene requirements.
- Do *not* operate the device if the water level in the humidification chamber exceeds the marked maximum.
- Regularly check the breathing circuit for condensation, and drain if required.
- Do *not* use this breathing circuit set with humidifiers other than the HAMILTON-H900.

⚠ CAUTION

- (USA only) *Federal law restricts this device to sale by or on the order of a physician.*
- *Check the breathing set for damage prior to use. Discard the breathing set if there is any sign of damage.*
- *Handle used products as contaminated goods according to local laws and regulations or hospital internal procedures.*
- *The temperature symbol on the temperature probe must always face up and be visible.*

NOTICE

- Before use, read the *HAMILTON-H900 Humidifier Instructions for Use*.
- Any incident with the product leading to serious patient injury, death, or a potential threat to public health must be reported to the manufacturer and the relevant authorities.
- The refill water should *not* be warmer than 37°C.
- Ensure that the water supply to the humidification chamber is functioning properly.

- Do *not* apply pressure to the water reservoir once the water refill has automatically stopped and the water has reached the maximum level.
- Do *not* add any drugs or medication directly to the water in the humidification chamber. If the HAMILTON-H900 is used in combination with any medical gases or nebulized medications, follow the *Instructions for use* of the supplied application and make sure that it is suitable for use with active humidification.
- Check the stability of all connections, prior to use.
- Set the appropriate ventilator alarms.
- Ensure the humidifier is placed in Standby or turned off before disconnecting any components.
- Use the stickers provided to indicate the date on which to replace the breathing circuit set.
- Incompatible parts can result in degraded performance, which can affect safety.

3 About the breathing circuit and breathing circuit sets

The breathing circuit comprises the limbs with integrated heater wire and integrated temperature probe. The breathing circuit set comprises the breathing circuit and the humidification chamber with water feed line.

4 Connecting the breathing circuit

Note that the limb connectors on the humidifier combine electrical connections with breathing circuit connectors. Ensure proper orientation of electrical contacts on the breathing circuit connectors to match the connecting element on the humidifier.

For an overview of the breathing circuit set and its components, see page 2.

To connect the breathing circuit set to the humidifier and ventilator

Illustrations for connecting the breathing circuit are provided after the instructions.

1. Hang the patient end of the breathing circuit on the tubing support arm.
2. Holding the humidification chamber in one hand and the breathing limbs to the ventilator in the other, insert the humidification chamber into the humidifier until it clicks into place (Figure 1).

3. Connect the BLUE ventilator-inspiratory limb to the *To patient* inspiratory port on the ventilator (Figure 2).

Ensure that you connect the BLUE inspiratory limb to the *To patient* inspiratory port on the ventilator, and the WHITE expiratory limb to the *From patient* expiratory port on the ventilator.

Connecting the limbs to the incorrect ports can cause a buildup of heat in the humidification chamber, with the risk of delivering overheated gas to the patient.

4. Connect the WHITE expiratory limb to the *From patient* expiratory port on the ventilator (Figure 3).
5. Insert the water spike into the water supply and open the air release cap on the spike (Figure 4).
6. Continue setting up the breathing circuit components, such as flow sensor and patient interface, as required for your patient.

Figure 1. Insert the humidification chamber


Figure 2. Connect BLUE inspiratory limb to To patient inspiratory port on the ventilator


Figure 3. Connect WHITE expiratory limb to From patient expiratory port on the ventilator



Figure 4. Insert water spike into water supply and open the air release cap on the spike

5 Delivery

The components are delivered clean and ready for clinical use.

6 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

7 Technical specifications

Parameter	Specification ¹
Flow range	4 to 120 l/min
Connectors ^{2,3} Ventilator: Patient:	ID22 ID15/OD22
Compliance (at 60 cmH ₂ O)	2 ml/cmH ₂ O
Compressible volume	1600 ml
Smallest internal diameter	19 mm
Water volume (min/max)	80 ml / 200 ml

Parameter	Specification ¹
Flow resistance (at 30 l/min) Inspiratory: Expiratory:	0.01 cmH ₂ O/(l/min) 0.007 cmH ₂ O/(l/min)
Flow resistance (at 15 l/min) Inspiratory: Expiratory:	0.008 cmH ₂ O/(l/min) 0.005 cmH ₂ O/(l/min)
Gas leakage (at 200 cmH ₂ O)	< 40 ml/min
Maximum oper- ating pressure	20 kPa
Circuit length	1.95 m (approx. 6 ft)
Compatibility	HAMILTON-H900 humidifier
Operating temperature	Recommended: 18°C to 26°C (65°F to 78°F)
Storage tempera- ture	-20°C to 60°C (-4°F to 140°F)
Operating and storage humidity	5% to 95% relative humidity, noncon- densing
Maximum oper- ating duration	28 days

Parameter	Specification ¹
Conforms to	IEC 60601-1:2005/ A1:2012 ISO 80601-2- 74:2017 EN ISO 5356-1:2015 EN ISO 5367:2014 IEC 62366-1:2015 ISO 10993-1:2018 ISO 13485:2016

1. All values presented are worst case scenarios.

2. Internal diameter (ID) and outer diameter (OD) in mm.

3. 15 mm and 22 mm conical connectors according to ISO 5356-1.

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For devices manufactured in Switzerland

Manufacturer

Hamilton Medical AG
Via Crusch 8, 7402 Bonaduz, Switzerland
☎ +41 (0) 58 610 10 20
info@hamilton-medical.com
www.hamilton-medical.com

medin Medical Innovations GmbH
Adam-Geisler-Straße 1
DE-82140 Olching
Germany

For devices manufactured and sold in the USA
Printed in the USA

Manufactured for

Hamilton Medical AG
Via Crusch 8, 7402 Bonaduz, Switzerland
☎ +41 (0) 58 610 10 20
info@hamilton-medical.com
www.hamilton-medical.com

Distributor in USA

Hamilton Medical, Inc.
201 Edison Way, Unit A Reno, NV 89502, US
☎ (775) 858-3200
Toll-free: (800) 426-6331
info@hamiltonmedical.com

624481/08
2024-02-07

Specifications are subject to change without notice. For all proprietary trademarks and third-party trademarks used by Hamilton Medical AG, see www.hamilton-medical.com/trademarks.
© 2024 Hamilton Medical AG. All rights reserved.