




Patient-ventilator asynchrony reference card

Asynchrony	Description	On the waveform	Waveform example	Common possible causes
Trigger asynchronies - during the beginning of inspiration				
Delayed triggering	The time interval between the patient's inspiratory effort and the delivery of a mechanical breath is increased	Flow waveform: Look for a longer-than-normal time interval between the positive deflection in flow ① and the delivery of ventilatory support ②		<ul style="list-style-type: none"> Trigger threshold set too high Ventilator pneumatics Presence of AutoPEEP Low respiratory drive Weak inspiratory effort
Ineffective effort	The patient's inspiratory effort fails to trigger the delivery of a mechanical breath	Flow waveform: Look for an abrupt change in the steepness of the waveform ① (decrease in expiratory flow or increase in inspiratory flow) that is not followed by ventilatory support ②		<ul style="list-style-type: none"> Trigger threshold set too high Pressure support too high Set frequency and/or inspiratory time too high (in controlled modes) Tidal volume set too high Presence of AutoPEEP Low respiratory drive Weak inspiratory effort Sedation
Auto triggering	A mechanical breath delivered without an inspiratory effort	Pressure waveform: Look for a delivered mechanical breath showing no drop in airway pressure ① at the beginning of the inspiratory phase		<ul style="list-style-type: none"> Trigger threshold set too low Air leaks in the endotracheal tube cuff, ventilator circuit, or chest tube Flow oscillations (water or secretion in the circuit, cardiac oscillations)
Flow asynchronies - during the gas delivery				
Flow asynchrony	The delivered flow does not meet the patient's inspiratory flow demands	Pressure waveform: Look for an upward concavity ① preceding the end of the mechanical breath		<ul style="list-style-type: none"> Inappropriate selection of ventilation mode (more frequent in volume-controlled modes) High inspiratory effort In volume-controlled modes: Inappropriate flow settings In pressure-controlled modes: Inappropriate P-ramp settings

..... "Correct" waveform, in case of good patient-ventilator synchrony

Patient factors | Ventilator-related factors | Patient-ventilator interface

Asynchrony	Description	On the waveform	Waveform example	Common possible causes
Termination asynchronies - during the end of inspiration				
Double triggering	Two (or more) mechanical breaths are delivered during one single inspiratory effort	Flow waveform: Look for two assisted breaths without expiration between them or with an expiration interval of less than half of the mean inspiratory time (often visually displayed as a waveform with two inspiratory peaks)		<ul style="list-style-type: none"> Cycling criteria (ETS) set too high Pressure support too low P-ramp too short Flow starvation High respiratory drive Time constant too short <p>Double triggering can be an effect of and/or promoted by reverse triggering or early cycling</p>
Early cycling	The duration of the mechanical breath is shorter than the duration of the patient's inspiratory effort	Flow waveform: Look for a small bump ¹ at the beginning of expiration (after peak expiratory flow) followed by an abrupt initial reversal in the expiratory flow ²		<ul style="list-style-type: none"> In pressure support ventilation: <ul style="list-style-type: none"> Cycling criteria (ETS) set too high Low levels of ventilator pressure support Time constant too short In time-cycled ventilation: <ul style="list-style-type: none"> Short inspiratory time
Delayed cycling	The duration of the mechanical breath is longer than the duration of the patient's inspiratory effort	Flow waveform: Look for a change in the slope of the inspiratory flow: a fast decrease ¹ followed by an exponential (less steep) decline ²		<ul style="list-style-type: none"> In pressure support ventilation: <ul style="list-style-type: none"> Cycling criteria (ETS) set too low Pressure support too high P-ramp too long In pressure control ventilation: <ul style="list-style-type: none"> Cycling criteria (ETS) set too low Inspiratory time too long In volume control ventilation: <ul style="list-style-type: none"> Low flow Long inspiratory time High tidal volume

