BiTrac[®] Select[™]

BiTrac MaxShield Select[™]

Pediatric Mask with Anti-Asphyxia Leak Elbow 22mm Male
Interchangeable Elbow Compatible

bg	Инструкции за употреба	lv	Lietošanas norādījumi
cz	Návod k použití	nl	Gebruiksaanwijzing
da	Brugsanvisning	no	Bruksanvisning
de	Gebrauchsanweisung	pl	Sposób użycia
el	Οδηγίες χρήσης	pt	Indicações de Utilização
es	Instrucciones de uso	ro	Instrucțiuni de utilizare
et	Kasutusjuhend	ru	Инструкция по применению
fi	Käyttöohjeet	sk	Návod na použitie
fr	Mode d'emploi	sl	Navodila za uporabo
hr	Upute za uporabu	sr	Uputstvo za upotrebu
hu	Használati útmutató	sv	Bruksanvisning
it	Istruzioni per l'uso	tr	Kullanma Talimatları
lt	Naudojimo instrukcija		

Distributor in the USA:
Hamilton Medical, Inc.
201 Edison Way, Unit A
Reno, NV 89502-2305, USA
Distributed By:
Hamilton Medical AG,
Via Crusch 8, 7402
Bonaduz, Switzerland

Intersurgical UAB
Tel: + 370 387 66 611
e-mail: info@intersurgical.lt
Amnionių g. 60, LT-18170, Pabradė, Lithuania

EC REP QNET BV
Kantstraat 19
NL-5076 NP Haaren
The Netherlands

Pulmonary Inc.
2055 Executive Drive
Indianapolis, IN 46241 USA
www.Pulmonary.com

Intended Use


The BiTrac MaxShield Select™ mask is an accessory for use with CPAP or Bilevel systems for ventilatory support. A minimum of >3.0cm of H₂O at the mask is required.

Contraindications

- 1. This device should not be used on uncooperative or unresponsive patients.
- 2. This device may not be suitable for use on patients taking medication which may cause vomiting, have nocturnal vomiting, hiatal hernia or impaired cardiac sphincter function.
- 3. This mask may not be suitable for use on patients with the following conditions: glaucoma, recent eye surgery or dry eyes.

Warnings

- 1. This mask has an integral leak port that provides a continuous leak path in the patient circuit. It is recommended that the circuit have an additional exhalation device. You may need to adjust the pressure level and oxygen therapy to compensate for the additional leak of the exhalation device.
- 2. This device is not suitable for providing life support ventilation.
- 3. Any unusual chest discomfort, shortness of breath, stomach distension, belching, severe headache, drying eyes, eye pain, eye infections or blurred vision when receiving positive airway pressure or upon awakening should be reported to the clinician.
- 4. To minimize risks of vomiting during sleep, avoid eating or drinking for three (3) hours before use.
- 5. Make sure the mask and headgear are the correct size for the patient's face. If excessive discomfort is experienced, adjust the mask and head gear. If the discomfort continues, consider an alternative mask.
- 6. This CPAP mask should only be used with CPAP systems recommended by your physician or respiratory therapist. A mask should not be used unless the CPAP system is turned ON and operating properly. The vent hole or holes associated with the mask should never be blocked.
Explanation of the Warning: CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation. This warning applies to most models of CPAP systems.

Operating Specifications: 5°C to 40°C
at humidity range of 15% to 95%

Storage Specifications: -20°C to 60°C at relative humidity up to 95%
non-condensing

Before Use

Verify entrainment valve function on anti-asphyxia elbow: With the unit turned OFF, the entrainment valve flap should be open, allowing user to breathe room air.

- Turn CPAP or Bilevel unit ON. The flap valve should close, allowing the user to receive airflow from the unit. If the valve does not function in this manner, replace the mask.
- Do not block the opening on the entrainment valve.
- Check valve before each use for malfunction or blockage caused by patient secretions.

Achieving the Proper Fit:

- 1. Use the size guide to measure the size of the patients face.
- 2. The inside cushion should rest above the eyebrow and directly below the bottom lip. Refer to Figure 1.
- 3. When the correct size is achieved place the proper size mask on the patient.
- 4. Place mask on patients face, and pull the bonnet overhead.
- 5. Clip the mask into place and adjust the straps. Refer to Figure 2.
- 6. Continue with CPAP therapy per your protocols.

Use of Interchangeable Elbow Feature:

- 1. Remove the existing elbow from the mask by holding the adapter ring and twisting the elbow counter clockwise to disengage the bayonet lock (Refer to Figure B.)
- 2. Gently pull the elbow from the mask. (Refer to Figure C.)
- 3. Gently press the desired Interchangeable Compatible Elbow into the mask adapter ring and turn clockwise until fully seated. Do not use excessive force.

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