REF 159001, 159002, 159003, 159005

2021-07-20

English | 10098449/01
Software version 2.8x**Add this Addendum to the front of your ventilator Operator's Manual.**

This Addendum summarizes changes and corrections to the HAMILTON-G5 and HAMILTON-S1 user manuals. Changes apply to software version 2.8x.

1 Corrections to the HAMILTON-G5 and HAMILTON-S1 Operator's Manuals

Preface: Training

Training

Hamilton Medical AG provides online training for users via the Hamilton Medical College. See <https://www.hamilton-medical.com/E-learning-and-Education/College.html>.

Section 1.4.3 Gas supply

 WARNING

Do *not* connect nitric oxide to the oxygen inlet; it is *not* permitted to use the ventilator with nitric oxide or mixtures of nitric oxide.

Section 1.4.2 Electrical: power and batteries

 WARNING

Anybody connecting additional medical equipment to the power sockets on the ventilator configures a medical system and is responsible for ensuring that the system complies with the requirements for medical electrical systems.

Section 1.9.3 O2 sensor

NOTICE

Keep the oxygen sampling site free of any other gases to avoid affecting oxygen measurements.

Section 5.5.2 About the trigger types

The following notice has been added:

NOTICE

IntelliSync+¹ is designed for use with all adult and pediatric patients weighing 10 kg or more.

¹ Not available in all markets.

Section 6.6 O2 enrichment for neonates

The applied oxygen concentration during the enrichment maneuver is increased by 10%.

When adjustable O2 enrichment is available¹, the applied oxygen concentration can be set in the **System > O2** enrichment window.

Section 9.1 Responding to alarms

The following safety messages have been added:

NOTICE

- Alarm conditions, including technical faults/events, that are *not* directly related to a physiological sensor (CO2, SpO2) do *not* affect the function of any attached physiological sensor, including the values of any associated CO2, SpO2, and pulse-rate measurements. Real-time waveforms on the ventilator provide a method for assessing the displayed numeric values.
- The factory default alarm limit settings are set in line with the selected patient group, allowing for unattended monitoring. These settings, however, can *never* replace individual review of the patient and adjustment of alarm limits based on their condition.

Section 9.4 Troubleshooting alarms

The CO2 calibration alarm was omitted from the alarms troubleshooting table. It is listed below.

Alarm	Definition
CO2 calibration needed	<i>Low priority.</i> A previous sensor zero calibration failed.

Action needed

Perform the following checks, repeating the calibration after each one, until calibration is successful:

- Clean or replace airway adapter.
- Perform a zero calibration of the sensor, making sure there is no source of CO2 near the airway adapter.
- Replace the airway adapter.
- Replace the CO2 sensor.
- If the problem persists, have the ventilator serviced.

Section 10.4.2 Suctioning maneuver

The following section provides information about performing closed-suctioning maneuvers.

NOTICE

When performing a closed suctioning maneuver, follow your institution's protocols.

Verify alarm limit settings and consider whether O₂ enrichment should be used prior to performing a closed-suctioning maneuver.

If the Suctioning tool is enabled on your device, ensure O₂ enrichment is not active when performing the closed-suctioning maneuver.

When performing a closed-suctioning maneuver, ventilation continues and the current settings do not need to be adjusted.

You can perform a closed-suctioning maneuver with the following pressure-controlled ventilation modes: APVcmv, APVsimv, P-CMV, P-SIMV, DuoPAP, APRV, SPONT, or ASV.

Section 15.1 Parts and accessories

A lead-free O₂ sensor (PN 10110239) is now available.

The new sensor is for use with HAMILTON-G5/S1 ventilators running software version 2.8x and higher.

The sensor should be replaced every 2 years or when depleted, whichever comes first.

Section 16.2 Environmental requirements

Water Protection

For the HAMILTON-G5/S1 ventilator, the degree of protection (solid particle and liquid ingress) has been updated to a rating of IP22.

Section 16.4 Electrical specifications

The input power specification for the HAMILTON-G5/S1 has been updated as follows:

Element	Specification
Input power	100 to 240 VAC, 50/60 Hz

Section 16.5 Corrections to control settings specifications

The following table lists the corrected values for some control settings.

Table 1. Control settings, ranges, and accuracy

Parameter or setting (unit)	Range: Adult/Pediatric	Range: Neonatal	Default: Adult/Pediatric	Default: Neonatal
End PEEP (cmH ₂ O)	0 to 35 ²	0 to 35 ²	Startup setting = PEEP	Startup setting = PEEP
Nebulizer Duration (min)	5 to 40	5 to 40	30	30
Nebulizer Synchronization	Inspiration, Exhalation, Insp. and Exh.	Inspiration, Exhalation, Insp. and Exh.	Inspiration	Inspiration
PEEP/CPAP (cmH ₂ O)	0 to 50	0 to 25	5	5
Pstart (cmH ₂ O)	0 to 35 ²	0 to 35 ²	Startup setting = PEEP	Startup setting = PEEP

² In some markets, the maximum is 20 cmH₂O.

Section 16.11 Technical performance data

The alarm volume (loudness) is measured at 1 meter distance from the ventilator. A setting of 1 = 57 dB(A), and 10 = 80 dB(A), with an accuracy of ± 6 dB(A).

Section 16.12 Functional description of ventilator system

The following changes have been made to this section:

- The functional description states that the ventilator system contains an integrated air compressing system. This is incorrect, and should instead state that there is a reservoir for compressed air.
- The gas supply and delivery system does not use room air, but instead high-pressure air.
- The description refers to the coordinated operations of the blower and expiratory valve. It should refer to the coordinated operations of the inspiratory and expiratory valves.

Section 16.13 Symbols used on device labels and packaging

The water protection rating of the HAMILTON-G5/S1 ventilator has changed, and the device label has been updated accordingly.

The following table describes the HAMILTON-G5/S1 ingress protection rating.

Symbols used on device labels and packaging

Symbol	Definition
IP22	Protected against dripping water when the device is tilted to a maximum of 15 degrees, and from solid particles larger than 12.5 mm.

2 Corrections and changes to the INTELLiVENT-ASV Operator's Manual for HAMILTON-G5/S1

Section 1.4.11.2 Automatic recruitment maneuvers

Note that automatic recruitment is not available in all markets.

Section 1.4.11.3 Target shift

For software versions 2.81 and higher, the following changes have been made to the PetCO₂ and SpO₂ target shift limits:

Table 2. PetCO₂ target shift limits

PetCO ₂ target shift limits	
All patient conditions	-20 mmHg to 20 mmHg ³

Table 3. SpO₂ target shift limits based on patient condition

SpO ₂ target shift limits	
Brain Injury	-5% ⁴ to +2%
<i>Target shift limits based on all other patient conditions remain unchanged.</i>	

Section 3.2 Technical data

The table below lists the corrected response time performance specifications for the **Oxygen** controller.

Table 4. Performance specifications, Oxygenation controller

	Oxygen	PEEP
Response time (90% of steady state)	N/A, only target range for SpO ₂ specified	6 minutes

³ In some markets, the upper target shift limit is 10 mmHg.

⁴ In some markets, the lower target shift limit is -2%.

3 Corrections and changes to the Pulse oximetry Instructions for use for HAMILTON-G5/S1

SpO₂-related updates

When a supported pulse oximeter is connected to the device, the HAMILTON-G5/S1 ventilators provide integrated monitoring and data display of functional oxygen saturation of arterial hemoglobin (SpO₂) and related pulse oximetry data.

Section 4.3 Nihon Kohden Technical specifications

Ingress protection

For Nihon Kohden pulse oximeters, the degree of protection (solid particle and liquid ingress) has been updated to a rating of IPX2.

Wavelength and maximum light intensity specifications

Nihon Koden pulse oximeters have two wavelengths with peaks in the range of 650 and 950 nm. The maximum light intensity is less than 5.5 mW/sr. The wavelength range can provide especially useful information to clinicians.

Section 5.3 Masimo Technical specifications

Ingress protection

For Masimo pulse oximeters, the degree of protection (solid particle and liquid ingress) has been updated to a rating of IP22.

Wavelength specifications

The wavelength range, provided in the *Pulse Oximetry Instructions for Use*, can provide especially useful information to clinicians.

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