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1. Borges, João et al. "Reversibility of lung collapse and hypoxemia in early acute respiratory distress syndrome (早期急性呼吸窘迫综合征病人肺塌陷及低氧血症的可逆性)". Am J Respir Crit Care Med 174.3 (2006): 268-278.
2. Gattinoni, Luciano et al. "Lung recruitment in patients with the acute respiratory distress syndrome (急性呼吸窘迫综合征病人肺复张)". New England Journal of Medicine 354.17 (2006): 1775-1786.
3. Maggiore, Salvatore et al. "Alveolar Derecruitment At Decremental Positive End-Expiratory Pressure Levels In Acute Lung Injury (急性肺创伤病人递减呼气末正压水平的肺泡去复张)". Am J Respir Crit Care Med 164.5 (2001): 795-801. 网络
4. Grasso, Salvatore et al. "Effects Of High Versus Low Positive End-Expiratory Pressures In Acute Respiratory Distress Syndrome (高/低呼气末正压对急性呼吸窘迫综合征的影响)". Am J Respir Crit Care Med 171.9 (2005): 1002-1008.
5. Demory, Didier et al. "Recruitability Of The Lung Estimated By The Pressure Volume Curve Hysteresis In ARDS Patients (ARDS 病人压力容量曲线滞后估计的肺复张性)". Intensive Care Med 34.11 (2008): 2019-2025.
6. Lu, Qin et al. "A Simple Automated Method For Measuring Pressure - Volume Curves During Mechanical Ventilation (一种用于在机械通气期间测量压力 - 容量曲线的简单自动化方法)". Am J Respir Crit Care Med 159.1 (1999): 275-282. 网络
7. Gattinoni, Luciano et al. "Pressure-Volume Curve Of Total Respiratory System In Acute Respiratory Failure: Computed Tomographic Scan Study (急性呼吸衰竭中总呼吸系统的压力 - 容量曲线: 计算机断层扫描研究)". Am Rev Respir Dis 136.3 (1987): 730-736. 网络
8. Kacmarek, Villar et al. "Open Lung Approach for the Acute Respiratory Distress Syndrome: A pilot, randomized controlled trial (急性呼吸窘迫综合征肺开放治疗方法: 一项初步随机对照试验)". Crit Care Med. 2016 Jan;44(1):32-42.

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