



HAMILTON-C3

The compact high-end ventilator

HAMILTON
MEDICAL



Our passion. Intelligent Ventilation solutions

Ever since we were founded in 1983, our focus has been on supporting the frontline heroes of critical care - with ventilation technologies that are safe, effective, and lung-protective. And we want to lessen the load for those who make extraordinary efforts every day, helping seriously ill patients fight their way back to health.

That is why we are committed to helping medical teams deliver the best respiratory care - to anyone, anywhere. That commitment is present in everything we do.

A handwritten signature in blue ink, appearing to read 'B. Hamilton'.

Bob Hamilton
CEO

Meet the HAMILTON-C3

The HAMILTON-C3 ventilator is a modular high-end ventilation solution for all patient groups. Offering a number of unique features, the HAMILTON-C3 is one of the first ventilators featuring "the bedside assistant" INTELLiVENT®-ASV®. The HAMILTON-C3's compact design and independence from compressed air allow maximum mobility throughout the hospital.

- ✓ Automated control of the patient's ventilation and oxygenation with INTELLiVENT-ASV
- ✓ P/V Tool®: Protective Ventilation Tool for lung assessment and recruitment
- ✓ High flow oxygen therapy
- ✓ High-performance NIV ventilation
- ✓ Compact and independent from compressed air
- ✓ Adult, pediatric, and neonatal ventilation



Slender, flexible, convenient

For intrahospital transport

The high-performance turbine enables the HAMILTON-C3 mechanical ventilator to be completely independent from compressed air, and its integrated high-capacity battery allows you to ventilate your patients during intrahospital transport without the need for an external power source. The compact design of the HAMILTON-C3 mechanical ventilator makes handling much easier. The HAMILTON-C3 can be attached to a trolley with ease.





Ease of use

In close cooperation with users and ventilation experts, our engineers have designed a user interface that is particularly intuitive. Switching between the HAMILTON-C3 and all other Hamilton Medical ventilators is easy, because they are all operated according to the same principles.

The Ventilation Cockpit on the HAMILTON-C3 consolidates the monitoring data and displays it as intuitive graphics. These provide a quick overview of the patient's current ventilation status and provide a reliable basis for therapy decisions.

“

The HAMILTON-C3 ventilators were selected for our NICU because they give us all of the modes that we need, the price met our budget needs, and we liked the ease of use, the battery life of the ventilator, and the ability to use them for transport to and from procedures.

Robert Lopez, Director of Respiratory Care
University Medical Center, Lubbock (TX), USA



The Ventilation Cockpit

1 Main monitoring parameters

All of the main monitoring parameters and alarm limits at a glance. The large characters allow you to see them even from a distance.

2 Dynamic Lung

One quick look shows you tidal volume, lung compliance, patient triggering, and resistance in real-time. The lungs expand and contract in synchrony with the actual breaths.

3 Vent Status

The Vent Status panel displays six parameters related to the patient's dependence on the ventilator. When all values are in the weaning zone, the panel is framed in green, indicating that spontaneous breathing trials or extubation can be considered.

4 Direct access to main controls

Access and adjust the most important controls for the current mode directly on the main display.



Individualized, lung-protective ventilation

The features available on the HAMILTON-C3 help you to individualize your patient's ventilation and to implement a lung-protective ventilation strategy.

Adaptive, lung-protective ventilation with ASV

- ✓ Supports the earliest possible spontaneous breathing by the patient^{1, 2}
- ✓ Shortens the ventilation time in various patient groups^{1, 2}

Your bedside assistant INTELLiVENT-ASV

- ✓ Requires fewer manual adjustments than conventional ventilation, consequently reducing the workload for the healthcare team³
- ✓ Follows the latest recommendations for lung-protective ventilation in terms of tidal volumes, driving pressure, and mechanical power^{4, 5, 6}

Lung assessment and recruitment with the P/V Tool

- ✓ Hysteresis of the pressure/volume curve can be used for assessing the recruitability of the lung at the bedside⁷
- ✓ May reduce the need for assessing recruitability with a CT scan, when using the PV loop in early on-set ARDS⁸

¹ Kirakli C. *Eur Respir J*. 2011 Oct;38(4):774-80

² Chen CW. *Respir Care*. 2011 Jul;56(7):976-83

³ Bialals, E., et al., *Minerva Anesthesiol*, 2016. 82(6): p. 657-68


⁴ Arnal JM. *Intensive Care Med* Exp 2016, 4(Suppl 1):A602

⁵ Arnal, J.-M., M. Saoli, and A. Garnero, *Heart & Lung: The Journal of Cardiopulmonary and Acute Care*. 2019 Nov

⁶ Buiteman-Kruizinga LA. *Crit Care Explor*. 2021 Feb 15;3(2):e0335

⁷ Demory D. *Intensive Care Med*. 2008 Nov;34(11):2019-25

⁸ Chiumello D. *Crit Care Med*. 2020 Oct;48(10):1494-1502



[Adaptive Support Ventilation \(ASV\)](#)

adjusts respiratory rate, tidal volume, and inspiratory pressure continuously, depending on the patient's lung mechanics and effort. ASV adapts ventilation breath-by-breath, 24 hours a day, and from intubation to extubation.



[INTELLiVENT-ASV, your bedside assistant](#)

is an advanced ventilation mode based on ASV. The clinician defines the clinical goal for PetCO₂ and SpO₂. INTELLiVENT-ASV then adjusts CO₂ elimination and oxygenation, and keeps the patient within the predefined ranges. Quick Wean supports the clinician in weaning patients from mechanical ventilation.




[P/V Tool for lung assessment and recruitment](#)

helps you assess recruitability and set PEEP based on respiratory mechanics. It also provides a repeatable method for quickly performing recruitment maneuvers.

Features and options



State-of-the-art ventilation modes



Serial interface for connection to PDMS or patient monitors



Integrated high flow oxygen therapy



On-screen help for alarm troubleshooting



Integrated pneumatic and optional Aerogen nebulizer



High-performance noninvasive ventilation (NIV)




Pulse oximetry (SpO2 and pulse measurement)



Adult, pediatric, and neonatal ventilation



Mainstream (volumetric) and sidestream capnography



High-performance turbine

From the ventilation specialist

Academy

Our Academy provides free, open-access educational content about mechanical ventilation and ventilators. Join here:

www.hamilton-medical.com/Academy.

Universal ventilator consumables

Our accessories and consumables are specially developed for the highest possible patient safety and ease of use. Choose between reusable and disposable parts according to your institutional policies.

Peripheral devices

Our ventilation portfolio includes an active humidifier, the HAMILTON-H900, as well as our automatic cuff pressure controller, IntelliCuff. Both devices may be used with all kinds of mechanical ventilators.





More information and free simulation software:
www.hamilton-C3.com



Manufacturer:

Hamilton Medical AG

Via Crusch 8, 7402 Bonaduz, Switzerland

+41 58 610 10 20

info@hamilton-medical.com

www.hamilton-medical.com

ELO202503035.00

The products shown here are not available for purchase by the general public. The information provided here is intended for healthcare professionals only. Always read the labels and follow the product's instructions for use. Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. All images are used for illustrative purposes only and may not accurately represent the product or its use. For all proprietary trademarks (®) and third-party trademarks used by Hamilton Medical AG see www.hamilton-medical.com/trademarks. ©2025 Hamilton Medical AG. All rights reserved.

HAMILTON-C3