INTELLiVENT-ASV

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1	Closed-loop ventilation mode (IntelliVent-ASV) in intensive care unit: a randomized trial of ventilation delivered (重症监护室闭环通气模式 (IntelliVent-ASV): 一项对输送通气量的随机试验)	3
2	Evaluation of fully automated ventilation: a randomized controlled study in post-cardiac surgery patients (全自动通气评估: 心脏手术后病人的随机对照研究)	4
3	A prospective comparison of the efficacy and safety of fully closed-loop control ventilation (Intellivent-ASV) with conventional ASV and SIMV modes (全闭环控制通气 (Intellivent-ASV) 与常规 ASV 及 SIMV 模式疗效和安全性的前瞻性比较)	5
4	Fully automated closed-loop ventilation is safe and effective in post-cardiac surgery patients (全自动闭环通气是安全、有效的心脏手术后病人的通气模式)	6
5	Prospective randomized crossover study of a new closed-loop control system versus pressure support during weaning from mechanical ventilation (机械通气撤机过程中新闭环控制系统与压力支持的前瞻性随机交叉研究)	7
6	Feasibility study on full closed-loop control ventilation (INTELLiVENT-ASV) in ICU patients with acute respiratory failure: a prospective observational comparative study (急性呼吸衰竭 ICU 病人全闭环控制通气 (INTELLiVENT-ASV) 可行性研究: 前瞻性观察比较研究)	8
7	Safety and efficacy of a fully closed-loop control ventilation (INTELLiVENT-ASV) in sedated ICU patients with acute respiratory failure: a prospective randomized crossover study (全闭环控制通气 (INTELLiVENT-ASV) 对急性呼吸衰竭 ICU 镇静治疗病人的安全性和疗效: 前瞻性随机交叉研究)	9
8	A rational framework for selecting modes of ventilation (用于通气模式选择的一种合理构架)	10
9	Adaptive support ventilation with and without end-tidal CO ₂ closed loop control versus conventional ventilation (适应性支持通气 (有/无呼气末二氧化碳闭环控制) 与常规通气的比较)	10
10	A pilot prospective study on closed loop controlled ventilation and oxygenation in ventilated children during the weaning phase. (对撤机阶段通气儿童闭环控制通气和氧合状态的初步前瞻性研究)	11
	其它文献	12
11	Ventilation with lower tidal volumes as compared with traditional tidal volumes for acute lung injury and acute respiratory distress syndrome (急性肺损伤和急性呼吸窘迫综合征病人的低潮气量通气与常规潮气量通气的比较)	12
12	Higher versus lower positive end-expiratory pressures in patients with acute respiratory distress syndrome (急性呼吸窘迫综合征病人的高/低呼气末正压通气比较)	13

13	Conservative oxygen therapy in mechanically ventilated patients: a pilot before-and-after trial (机械通气病人的保守氧气疗法: 初步前后对照试验)	13
14	Optimal duration of a sustained inflation recruitment maneuver in ARDS patients (对 ARDS 病人进行持续充气肺复张操作的最佳时长)	14
15	Higher vs lower positive end-expiratory pressure in patients with acute lung injury and acute respiratory distress syndrome: systematic review and meta-analysis (急性肺损伤和急性呼吸窘迫综合征病人的高/低呼气末正压通气比较: 系统论述和 meta 分析)	14
16	BTS guideline for emergency oxygen use in adult patients (用于成年病人紧急用氧的 BTS 指南)	15
17	Low mortality rate in adult respiratory distress syndrome using low-volume, pressure-limited ventilation with permissive hypercapnia: a prospective study (通过降低潮气量并限制压力的通气方式及允许性高碳酸血症降低成人呼吸窘迫综合征病人的死亡率: 前瞻性研究)	15
18	Bench-to-bedside review: hypercapnic acidosis in lung injury: from 'permissive' to 'therapeutic' (实验室到床旁的转化医学论述: 肺损伤中的高碳酸血性酸中毒: 从“允许性”到“治疗性”高碳酸血性酸中毒)	16
19	The work of breathing (呼吸做功)	17

Closed-loop ventilation mode (IntelliVent-ASV) in intensive care unit: a randomized trial of ventilation delivered (重症监护室闭环通气模式 (IntelliVent-ASV): 一项对输送通气量的随机试验)

Bialais E, Wittebole X, Vignaux L, Roeseler J, Wysocki M, Meyer J, Reychler G, Novotni D, Sottiaux T, Laterre PF, Hantson P

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Evaluation of fully automated ventilation: a randomized controlled study in post-cardiac surgery patients (全自动通气评估: 心脏手术后病人的随机对照研究)

Lellouche F, Bouchard PA, Simard S, L'Her E, Wysocki M
Intensive Care Med. 2013 Mar;39(3):463-71
PMID 23338569, <http://www.ncbi.nlm.nih.gov/pubmed/23338569>

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A prospective comparison of the efficacy and safety of fully closed-loop control ventilation (Intellivent-ASV) with conventional ASV and SIMV modes (全闭环控制通气 (Intellivent-ASV) 与常规 ASV 及 SIMV 模式疗效和安全性的前瞻性比较)

Abutbul A, Svirin S, Zbedat W, Linton DM, van Heerden PV
S Afr J Crit Care. 2014 Aug;30(1):28-32

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Fully automated closed-loop ventilation is safe and effective in post-cardiac surgery patients (全自动闭环通气是安全、有效的心脏手术后病人的通气模式)

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
Prospective randomized crossover study of a new closed-loop control system versus pressure support during weaning from mechanical ventilation (机械通气撤机过程中新闭环控制系统与压力支持的前瞻性随机交叉研究)

Clavieras N, Wysocki M, Coisel Y, Galia F, Conseil M, Chanques G, Jung B, Arnal JM, Matecki S, Molinari N, Jaber S

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Adaptive support ventilation with and without end-tidal CO2 closed loop control versus conventional ventilation (适应性支持通气 (有/无呼气末二氧化碳闭环控制) 与常规通气的比较)

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A pilot prospective study on closed loop controlled ventilation and oxygenation in ventilated children during the weaning phase. (对撤机阶段通气儿童闭环控制通气和氧合状态的初步前瞻性研究)

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Ventilation with lower tidal volumes as compared with traditional tidal volumes for acute lung injury and acute respiratory distress syndrome (急性肺损伤和急性呼吸窘迫综合征病人的低潮气量通气与常规潮气量通气的比较)

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Higher versus lower positive end-expiratory pressures in patients with acute respiratory distress syndrome (急性呼吸窘迫综合征病人的高/低呼气末正压通气比较)

Brower RG, Lanken PN, MacIntyre N, Matthay MA, Morris A, Ancukiewicz M, Schoenfeld D, Thompson BT; National Heart, Lung, and Blood Institute ARDS Clinical Trials Network
N Engl J Med. 2004 Jul 22;351(4):327-36
PMID 15269312, <http://www.ncbi.nlm.nih.gov/pubmed/15269312>

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Conservative oxygen therapy in mechanically ventilated patients: a pilot before-and-after trial (机械通气病人的保守氧气疗法：初步前后对照试验)

Suzuki S, Eastwood GM, Glassford NJ, Peck L, Young H, Garcia-Alvarez M, Schneider AG, Bellomo R
Crit Care Med. 2014 Jun;42(6):1414-22
PMID 24561566, <http://www.ncbi.nlm.nih.gov/pubmed/24561566>

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Optimal duration of a sustained inflation recruitment maneuver in ARDS patients (对 ARDS 病人进行持续充气肺复张操作的最佳时长)

Arnal JM, Paquet J, Wysocki M, Demory D, Donati S, Granier I, Corno G, Durand-Gasselín J

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Higher vs lower positive end-expiratory pressure in patients with acute lung injury and acute respiratory distress syndrome: systematic review and meta-analysis (急性肺损伤和急性呼吸窘迫综合征病人的高/低呼气末正压通气比较: 系统论述和 meta 分析)

Briel M, Meade M, Mercat A, Brower RG, Talmor D, Walter SD, Slutsky AS, Pullenayegum E, Zhou Q,

Cook D, Brochard L, Richard JC, Lamontagne F, Bhatnagar N, Stewart TE, Guyatt G

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O'Driscoll BR, Howard LS, Davison AG; British Thoracic Society
Thorax. 2008 Oct;63 Suppl 6:vi1-68
PMID 18838559, <http://www.ncbi.nlm.nih.gov/pubmed/18838559>

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OTIS AB