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1	A randomized controlled trial comparing the ventilation duration between Adaptive Support Ventilation and Pressure Assist/Control Ventilation in medical patients in the ICU (ICU 内科病人适应性支持通气时间与压力辅助/控制通气时间比较的随机对照试验) ..	4
2	Adaptive support ventilation for faster weaning in COPD: a randomised controlled trial (适应性支持通气用于 COPD 病人更快撤机: 随机对照试验)	5
3	Adaptive support ventilation for fast tracheal extubation after cardiac surgery: a randomized controlled study (适应性支持通气用于心脏手术后快速气管拔管: 随机对照研究)	6
4	A randomized controlled trial of adaptive support ventilation mode to wean patients after fast-track cardiac valvular surgery (快通道心脏瓣膜手术后病人撤机过程中适应性支持通气模式的随机对照试验)	7
5	Adaptive Support Ventilation versus Synchronized Intermittent Mandatory Ventilation with Pressure Support in weaning patients after orthotopic liver transplantation (原位肝移植术后病人撤机过程中适应性支持通气与同步间歇指令通气辅以压力支持的比较)	7
6	Randomized controlled trial comparing adaptive-support ventilation with pressure-regulated volume-controlled ventilation with automode in weaning patients after cardiac surgery (心脏手术后病人撤机过程中适应性支持通气与带自动模式的压力调节容量控制通气比较的随机对照试验)	8
7	Automatic "respirator/weaning" with adaptive support ventilation: the effect on duration of endotracheal intubation and patient management (采用适应性支持通气的自动“呼吸机/撤机”: 对气管内插管时间和病人管理的影响)	8
8	Effects of implementing adaptive support ventilation in a medical intensive care unit (内科重症监护室中实施适应性支持通气的效用)	9
9	Adaptive support ventilation for complete ventilatory support in ARDS: a pilot randomized controlled trial (适应性支持通气用于 ARDS 完全通气支持: 初步随机对照试验)	10
10	Human versus Computer Controlled Selection of Ventilator Settings: An Evaluation of Adaptive Support Ventilation and Mid-Frequency Ventilation (呼吸机设置的人工选择与计算机控制选择比较: 适应性支持通气和中频通气的评估)	11
11	Adaptive support ventilation versus conventional ventilation for total ventilatory support in acute respiratory failure (急性呼吸衰竭中适应性支持通气与常规通气(完全通气支持)的比较)	12
12	Automatic selection of breathing pattern using adaptive support ventilation (利用适应性支持通气自动选择呼吸形式)	13
13	A randomized controlled trial of 2 protocols for weaning cardiac surgical patients receiving adaptive support ventilation (接受适应性支持通气的的心脏手术病人 2 个撤机方案的随机对照试验)	14

14	Comparing the effects of adaptive support ventilation and synchronized intermittent mandatory ventilation on intubation duration and hospital stay after coronary artery bypass graft surgery (比较适应性支持通气与同步间歇指令通气对冠状动脉搭桥术后病人插管时间及住院时间的影响)	14
15	Intelligent ventilation in the intensive care unit (重症监护室的智能通气)	15
16	Adaptive support and pressure support ventilation behavior in response to increased ventilatory demand (通气需求增加时适应性支持和压力支持通气的响应特性)	15
17	Determinants of tidal volumes with adaptive support ventilation: a multicentre observational study (采用适应性支持通气时潮气量的决定因素: 多中心观察研究)	16
18	Evaluation of adaptive support ventilation in paralysed patients and in lung model (评估瘫痪病人及肺模型中的适应性支持通气)	16
19	Clinical experience with adaptive support ventilation for fast-track cardiac surgery (适应性支持通气用于快通道心脏手术的临床经验)	17
20	Patient-ventilator interactions during partial ventilatory support: a preliminary study comparing the effects of adaptive support ventilation with synchronized intermittent mandatory ventilation plus inspiratory pressure support (部分通气支持过程中病人与呼吸机的人机交互: 适应性支持通气与同步间歇指令通气(辅以吸气压力支持)效用比较的初步研究)	17
21	Automatic weaning from mechanical ventilation using an adaptive lung ventilation controller (利用适应性肺通气控制器实施机械通气自动撤机)	18
22	Effects of adaptive support ventilation and synchronized intermittent mandatory ventilation on peripheral circulation and blood gas markers of COPD patients with respiratory failure (适应性支持通气和同步间歇指令通气对呼吸衰竭 COPD 病人外周循环及血气指标的影响)	18
23	Comparison of 3 modes of automated weaning from mechanical ventilation: a bench study (机械通气三种自动撤机模式比较: 实验室研究)	19
24	Correlation between the %MinVol setting and work of breathing during adaptive support ventilation in patients with respiratory failure (呼吸衰竭病人适应性支持通气过程中 %MinVol 设置与呼吸做功的相关性)	19
25	Adaptive support ventilation prevents ventilator-induced diaphragmatic dysfunction in piglet: an in vivo and in vitro study (适应性支持通气预防小猪呼吸机相关性膈肌功能障碍: 体内外研究)	20
26	Adaptive support ventilation: an appropriate mechanical ventilation strategy for acute respiratory distress syndrome? (适应性支持通气: 一种适用于急性呼吸窘迫综合征的机械通气策略?)	20

27	A comparison of adaptive Support Ventilation (ASV) and Conventional Volume-Controlled Ventilation on Respiratory Mechanics in Acute Lung Injury/ARDS (急性肺损伤/ARDS中适应性支持通气 (ASV) 与常规容量控制通气对呼吸力学的影响比较)	21
28	Adaptive support ventilation for gynaecological laparoscopic surgery in Trendelenburg position: bringing ICU modes of mechanical ventilation to the operating room (适应性支持通气用于 Trendelenburg 截石位体位妇科腹腔镜手术: 将 ICU 机械通气模式引入手术室) .	21
29	Adaptive Support Ventilation as the sole mode of ventilatory support in chronically ventilated patients (适应性支持通气作为长期通气病人的唯一通气支持模式)	22
30	Adaptive lung ventilation (ALV) during anesthesia for pulmonary surgery: automatic response to transitions to and from one-lung ventilation (适应性肺通气 (ALV) 用于肺手术麻醉: 自动响应单肺通气的转换)	22
31	Continuous use of an adaptive lung ventilation controller in critically ill patients in a multi-disciplinary intensive care unit (适应性肺通气控制器在跨学科重症监护室危重病人中的持续使用)	23
32	Automatic selection of tidal volume, respiratory frequency and minute ventilation in intubated ICU patients as start up procedure for closed-loop controlled ventilation (将自动选择 ICU 插管病人潮气量、呼吸频率和分钟通气量作为闭环控制通气的启动程序)	23
	其它文献	24
33	Adaptive support ventilation (适应性支持通气)	24
34	The work of breathing (呼吸做功)	24
35	Automated versus non-automated weaning for reducing the duration of mechanical ventilation for critically ill adults and children (自动撤机与非自动撤机在减少危重成人及儿童病人机械通气时间方面的效用比较)	24
36	Adaptive support ventilation: State of the art review (适应性支持通气: 最新论述)	25
37	Closed loop mechanical ventilation (闭环机械通气)	26

A randomized controlled trial comparing the ventilation duration between Adaptive Support Ventilation and Pressure Assist/Control Ventilation in medical patients in the ICU (ICU 内科病人适应性支持通气时间与压力辅助/控制通气时间比较的随机对照试验)

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Adaptive support ventilation for faster weaning in COPD: a randomised controlled trial (适应性支持通气用于 COPD 病人更快撤机：随机对照试验)

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A randomized controlled trial of adaptive support ventilation mode to wean patients after fast-track cardiac valvular surgery (快通道心脏瓣膜手术后病人撤机过程中适应性支持通气模式的随机对照试验)

Zhu F, Gomersall CD, Ng SK, Underwood MJ, Lee A

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Adaptive Support Ventilation versus Synchronized Intermittent Mandatory Ventilation with Pressure Support in weaning patients after orthotopic liver transplantation (原位肝移植术后病人撤机过程中适应性支持通气与同步间歇指令通气辅以压力支持的比较)

Celli P, Privato E, Ianni S, Babetto C, D' Arena C, Guglielmo N, Maldarelli F, Paglialonga G, Rossi M, Berloco PB, Ruberto F, Pugliese F

Transplant Proc. 2014 Aug 20 [Epub ahead of print]

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Randomized controlled trial comparing adaptive-support ventilation with pressure-regulated volume-controlled ventilation with automode in weaning patients after cardiac surgery (心脏手术后病人撤机过程中适应性支持通气与带自动模式的压力调节容量控制通气比较的随机对照试验)

Gruber PC, Gomersall CD, Leung P, Joynt GM, Ng SK, Ho KM, Underwood MJ
Anesthesiology. 2008 Jul;109(1):81-7
PMID 18580176, <http://www.ncbi.nlm.nih.gov/pubmed/18580176>

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Effects of implementing adaptive support ventilation in a medical intensive care unit (内科重症监护室中实施适应性支持通气的效用)

Chen CW, Wu CP, Dai YL, Perng WC, Chian CF, Su WL, Huang YC

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Adaptive support ventilation for complete ventilatory support in ARDS: a pilot randomized controlled trial (适应性支持通气用于 ARDS 完全通气支持: 初步随机对照试验)

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Human versus Computer Controlled Selection of Ventilator Settings: An Evaluation of Adaptive Support Ventilation and Mid-Frequency Ventilation (呼吸机设置的人工选择与计算机控制选择比较：适应性支持通气和中频通气的评估)

Mireles-Cabodevila E, Diaz-Guzman E, Arroliga AC, Chatburn RL
Crit Care Res Pract. 2012 Sep;2012:204314
PMID 23119152, <http://www.ncbi.nlm.nih.gov/pubmed/23119152>

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Adaptive support ventilation versus conventional ventilation for total ventilatory support in acute respiratory failure (急性呼吸衰竭中适应性支持通气与常规通气(完全通气支持)的比较)

Iotti GA, Polito A, Belliato M, Pasero D, Beduneau G, Wysocki M, Brunner JX, Braschi A, Brochard L, Mancebo J, Ranieri VM, Richard JC, Slutsky AS

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
A randomized controlled trial of 2 protocols for weaning cardiac surgical patients receiving adaptive support ventilation (接受适应性支持通气的心脏手术病人 2 个撤机方案的随机对照试验)

Tam MK, Wong WT, Gomersall CD, Tian Q, Ng SK, Leung CC, Underwood MJ
J Crit Care. 2016 Jun;33:163-8
PMID 27006266, <http://www.ncbi.nlm.nih.gov/pubmed/27006266>

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Comparing the effects of adaptive support ventilation and synchronized intermittent mandatory ventilation on intubation duration and hospital stay after coronary artery bypass graft surgery (比较适应性支持通气与同步间歇指令通气对冠状动脉搭桥术后病人插管时间及住院时间的影响)

Yazdannik A, Zarei H, Massoumi G
Iran J Nurs Midwifery Res. 2016 Mar-Apr;21(2):207-12.
PMID 27095997, <http://www.ncbi.nlm.nih.gov/pubmed/27095997>

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Intelligent ventilation in the intensive care unit (重症监护室的智能通气)

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Adaptive support and pressure support ventilation behavior in response to increased ventilatory demand (通气需求增加时适应性支持和压力支持通气的响应特性)

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Evaluation of adaptive support ventilation in paralysed patients and in lung model (评估瘫痪病人及肺模型中的适应性支持通气)

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Clinical experience with adaptive support ventilation for fast-track cardiac surgery (适应性支持通气用于快通道心脏手术的临床经验)

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Effects of adaptive support ventilation and synchronized intermittent mandatory ventilation on peripheral circulation and blood gas markers of COPD patients with respiratory failure (适应性支持通气和同步间歇指令通气对呼吸衰竭 COPD 病人外周循环及血气指标的影响)

Han L, Wang Y, Gan Y, Xu L
Cell Biochem Biophys. 2014 Apr;70(1):481-4
PMID 24748176, <http://www.ncbi.nlm.nih.gov/pubmed/24748176>

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Comparison of 3 modes of automated weaning from mechanical ventilation: a bench study (机械通气三种自动撤机模式比较: 实验室研究)

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Correlation between the %MinVol setting and work of breathing during adaptive support ventilation in patients with respiratory failure (呼吸衰竭病人适应性支持通气过程中 %MinVol 设置与呼吸做功的相关性)

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A comparison of adaptive Support Ventilation (ASV) and Conventional Volume-Controlled Ventilation on Respiratory Mechanics in Acute Lung Injury/ARDS (急性肺损伤/ARDS中适应性支持通气 (ASV) 与常规容量控制通气对呼吸力学的影响比较)

Choi I, Choi J, Hong S, Lim C, Koh Y
Kor J crit Care Med. 2009 Aug; 24(2): 59-63

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Adaptive Support Ventilation as the sole mode of ventilatory support in chronically ventilated patients (适应性支持通气作为长期通气病人的唯一通气支持模式)

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Adaptive lung ventilation (ALV) during anesthesia for pulmonary surgery: automatic response to transitions to and from one-lung ventilation (适应性肺通气 (ALV) 用于肺手术麻醉: 自动响应单肺通气的转换)

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Continuous use of an adaptive lung ventilation controller in critically ill patients in a multi-disciplinary intensive care unit (适应性肺通气控制器在跨学科重症监护室危重病人中的持续使用)

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Automatic selection of tidal volume, respiratory frequency and minute ventilation in intubated ICU patients as start up procedure for closed-loop controlled ventilation (将自动选择 ICU 插管病人潮气量、呼吸频率和分钟通气量作为闭环控制通气的启动程序)

Laubscher TP, Frutiger A, Fanconi S, Jutzi H, Brunner JX
Int J Clin Monit Comput. 1994 Feb;11(1):19-30
PMID 8195655, <http://www.ncbi.nlm.nih.gov/pubmed/8195655>

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Automated versus non-automated weaning for reducing the duration of mechanical ventilation for critically ill adults and children (自动撤机与非自动撤机在减少危重成人及儿童病人机械通气时间方面的效用比较)

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Adaptive support ventilation: State of the art review (适应性支持通气：最新论述)

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Closed loop mechanical ventilation (闭环机械通气)

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Otis AB