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1	Mortality and pulmonary mechanics in relation to respiratory system and transpulmonary driving pressures in ARDS (ARDS 病人死亡率和肺力学与呼吸系统和跨肺驱动压有关)	4
2	Mechanical ventilation guided by esophageal pressure in acute lung injury (食道压引导的急性肺损伤病人机械通气)	5
3	Lung stress and strain during mechanical ventilation for acute respiratory distress syndrome (急性呼吸窘迫综合征病人机械通气期间的肺应力和应变)	6
4	Esophageal and transpulmonary pressures in acute respiratory failure (急性呼吸衰竭病人的食道压和跨肺压)	7
5	Recruitment maneuvers: using transpulmonary pressure to help Goldilocks (肺复张术: 使用跨肺压来获得恰到好处的效果)	8
6	Recruitment maneuvers and positive end-expiratory pressure titration in morbidly obese ICU patients. (使用肺复张操作和呼气末正压滴定治疗病态肥胖 ICU 病人。)	9
7	Volume delivered during recruitment maneuver predicts lung stress in acute respiratory distress syndrome (在肺复张操作过程中输送的容量可预测急性呼吸窘迫综合征病人的肺应力)	9
8	Effect of body mass index in acute respiratory distress syndrome (急性呼吸窘迫综合征病人体重指数的影响)	10
9	Transpulmonary pressure and gas exchange during decremental PEEP titration in pulmonary ARDS patients (肺源性 ARDS 病人递减 PEEP 滴定期间的跨肺压和气体交换)	11
10	Acute respiratory distress syndrome caused by pulmonary and extrapulmonary disease. Different syndromes? (肺源性及非肺源性疾病引起的急性呼吸窘迫综合征是不同的疾病?)	12
11	Alterations of lung and chest wall mechanics in patients with acute lung injury: effects of positive end-expiratory pressure (急性肺损伤病人肺及胸壁力学的变化: 呼气末正压的影响)	13
12	The occlusion tests and end-expiratory esophageal pressure: measurements and comparison in controlled and assisted ventilation (阻断试验和呼气末食道压: 控制和辅助通气的比较与测量)	14
13	In vivo calibration of esophageal pressure in the mechanically ventilated patient makes measurements reliable (机械通气病人食道压体内校准使测量可靠)	15
14	Do spontaneous and mechanical breathing have similar effects on average transpulmonary and alveolar pressure? A clinical crossover study (自主和机械呼吸是否对平均跨肺压和肺泡压力有类似影响? 一项临床交叉研究)	15

15	Non-invasive assessment of lung elastance in patients with acute respiratory distress syndrome (急性呼吸窘迫综合征病人肺弹性的无创评估)	16
16	Positive end expiratory pressure titrated by transpulmonary pressure improved oxygenation and respiratory mechanics in acute respiratory distress syndrome patients with intra-abdominal hypertension (通过跨肺压滴定的呼气末正压可促进腹内高压的急性呼吸窘迫综合征病人的氧合状态和呼吸力学)	16
17	Comparison of 2 correction methods for absolute values of esophageal pressure in subjects with acute hypoxemic respiratory failure, mechanically ventilated in the ICU (比较 2 个在 ICU 进行机械通气后纠正急性低氧性呼吸衰竭受试者食道压绝对值的方法) ..	17
18	ECMO criteria for influenza A (H1N1)-associated ARDS: role of transpulmonary pressure (甲型 H1N1 流感相关性 ARDS 的 ECMO 标准: 跨肺压的作用) ..	18
19	Esophageal pressures in acute lung injury: do they represent artifact or useful information about transpulmonary pressure, chest wall mechanics, and lung stress? (急性肺损伤病人的食道压是跨肺压、胸壁力学及肺应力的假象还是有用信息?)	19
20	Influence of lung and chest wall compliances on transmission of airway pressure to the pleural space in critically ill patients (肺及胸壁顺应性对气道压传送至危重症病人胸膜腔的影响)	19
21	A simple method for assessing the validity of the esophageal balloon technique (一个评估食道气囊技术有效性的简单方法)	20
22	Pulmonary, chest wall, and lung-thorax elastances in acute respiratory failure (急性呼吸衰竭病人的肺、胸壁及胸肺弹性)	20
23	Topography of esophageal pressure as a function of posture in man (人的体位与食道压分布形态的函数关系)	21
24	Value and limitations of transpulmonary pressure calculations during intra-abdominal hypertension (在腹内高压期间跨肺压计算的数值和限制)	21
25	Pleural pressure and optimal positive end-expiratory pressure based on esophageal pressure versus chest wall elastance: incompatible results (基于食道压与胸壁弹性的胸腔压力和最佳呼气末正压: 不一致的结果)	22
26	Volume-related and volume-independent effects of posture on esophageal and transpulmonary pressures in healthy subjects (体位对健康受试者食道压及跨肺压的容量相关性及其容量独立性影响)	22
27	Recruitment and derecruitment during acute respiratory failure: an experimental study (急性呼吸衰竭期间的肺复张和塌陷: 一项试验研究)	23
28	Validation of esophageal pressure occlusion test after paralysis (对麻痹后食道压阻断试验的验证)	23

29	Lung mechanics in sitting and horizontal body positions (在坐位及平躺位的肺力学)	24
	Additional files (其它文献)	25
30	Esophageal and transpulmonary pressure in the clinical setting: meaning, usefulness and perspectives (临床环境下的食道压和跨肺压: 含义、有效性和前景)	25
31	The application of esophageal pressure measurement in patients with respiratory failure (食道压测量对呼吸衰竭病人的应用)	25
32	Transpulmonary pressure: the importance of precise definitions and limiting assumptions (跨肺压: 准确定义和限制性假设的重要性)	26
33	The promises and problems of transpulmonary pressure measurements in acute respiratory distress syndrome (急性呼吸窘迫综合征病人跨肺压测量的预示作用和问题)	26
34	Measurement of esophageal pressure at bedside: pros and cons (食道压临床测量的利弊) .	26
35	Two steps forward in bedside monitoring of lung mechanics: transpulmonary pressure and lung volume (在肺力学临床监测方面向前迈进了两步: 跨肺压和肺容量)	27
36	Goal-directed mechanical ventilation: are we aiming at the right goals? A proposal for an alternative approach aiming at optimal lung compliance, guided by esophageal pressure in acute respiratory failure (目标导向的机械通气: 我们是否瞄准正确的目标? 关于通过食道压引导瞄准急性呼吸衰竭病人最佳肺顺应性的可选择方法的一项建议)	27
37	Esophageal pressure: benefit and limitations (食道压: 优点和不足)	27
38	Targeting transpulmonary pressure to prevent ventilator induced lung injury (确定跨肺压目标以预防呼吸机诱发性肺损伤)	28
39	Esophageal and gastric pressure measurements (食道压和胃压测量)	28
40	Respiratory mechanics in mechanically ventilated patients (机械通气病人的呼吸力学) ..	29

Mortality and pulmonary mechanics in relation to respiratory system and transpulmonary driving pressures in ARDS (ARDS 病人死亡率和肺力学与呼吸系统和跨肺驱动压有关)

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Mechanical ventilation guided by esophageal pressure in acute lung injury (食道压引导的急性肺损伤病人机械通气)

Talmor D, Sarge T, Malhotra A, O'Donnell CR, Ritz R, Lisbon A, Novack V, Loring SH

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Lung stress and strain during mechanical ventilation for acute respiratory distress syndrome (急性呼吸窘迫综合征病人机械通气期间的肺应力和应变)

Chiumello D, Carlesso E, Cadringer P, Caironi P, Valenza F, Polli F, Tallarini F, Cozzi P, Cressoni M, Colombo A, Marini JJ, Gattinoni L

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Esophageal and transpulmonary pressures in acute respiratory failure (急性呼吸衰竭病人的食道压和跨肺压)

Talmor D, Sarge T, O'Donnell CR, Ritz R, Malhotra A, Lisbon A, Loring SH
Crit Care Med. 2006 May;34(5):1389-94
PMID 16540960, <http://www.ncbi.nlm.nih.gov/pubmed/16540960>

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Volume delivered during recruitment maneuver predicts lung stress in acute respiratory distress syndrome (在肺复张操作过程中输送的容量可预测急性呼吸窘迫综合征病人的肺应力)

Beitler JR, Majumdar R, Hubmayr RD, Malhotra A, Thompson BT, Owens RL, Loring SH, Talmor D
Crit Care Med. 2016 Jan;44(1):91-9
PMID 26474111, <http://www.ncbi.nlm.nih.gov/pubmed/26474111>

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Transpulmonary pressure and gas exchange during decremental PEEP titration in pulmonary ARDS patients (肺源性 ARDS 病人递减 PEEP 滴定期间的跨肺压和气体交换)

Rodriguez PO, Bonelli I, Setten M, Attie S, Madorno M, Maskin LP, Valentini R
Respir Care. 2013 May;58(5):754-63
PMID 23051849, <http://www.ncbi.nlm.nih.gov/pubmed/23051849>

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Acute respiratory distress syndrome caused by pulmonary and extrapulmonary disease. Different syndromes? (肺源性及非肺源性疾病引起的急性呼吸窘迫综合征是不同的疾病?)

Gattinoni L, Pelosi P, Suter PM, Pedoto A, Vercesi P, Lissoni A
 Am J Respir Crit Care Med. 1998 Jul;158(1):3-11
 PMID 9655699, <http://www.ncbi.nlm.nih.gov/pubmed/9655699>

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Alterations of lung and chest wall mechanics in patients with acute lung injury: effects of positive end-expiratory pressure (急性肺损伤病人肺及胸壁力学的变化：呼气末正压的影响)

Pelosi P, Cereda M, Foti G, Giacomini M, Pesenti A
Am J Respir Crit Care Med. 1995 Aug;152(2):531-7
PMID 7633703, <http://www.ncbi.nlm.nih.gov/pubmed/7633703>

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The occlusion tests and end-expiratory esophageal pressure: measurements and comparison in controlled and assisted ventilation (阻断试验和呼气末食道压: 控制和辅助通气的比较与测量)

Chiumello D, Consonni D, Coppola S, Froio S, Crimella F, Colombo A
Ann Intensive Care. 2016 Dec;6(1):13
PMID 26868503, <http://www.ncbi.nlm.nih.gov/pubmed/26868503>

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In vivo calibration of esophageal pressure in the mechanically ventilated patient makes measurements reliable (机械通气病人食道压体内校准使测量可靠)

Mojoli F, Iotti GA, Torriglia F, Pozzi M, Volta CA, Bianzina S, Braschi A, Brochard L
Crit Care. 2016 Apr 11;20:98
PMID 27063290, <http://www.ncbi.nlm.nih.gov/pubmed/27063290>

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Do spontaneous and mechanical breathing have similar effects on average transpulmonary and alveolar pressure? A clinical crossover study (自主和机械呼吸是否对平均跨肺压和肺泡压力有类似影响? 一项临床交叉研究)

Bellani G, Grasselli G, Teggie-Droghi M, Mauri T, Coppadoro A, Brochard L, Pesenti A
Crit Care. 2016 Apr 28;20(1):142
PMID 27160458, <http://www.ncbi.nlm.nih.gov/pubmed/27160458>

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Non-invasive assessment of lung elastance in patients with acute respiratory distress syndrome (急性呼吸窘迫综合征病人肺弹性的无创评估)

Garnero A, Tuxen D, Ducros L, Demory D, Donati SY, Durand-Gasselín J, Cooper J, Hodgson C, Arnal JM

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Positive end expiratory pressure titrated by transpulmonary pressure improved oxygenation and respiratory mechanics in acute respiratory distress syndrome patients with intra-abdominal hypertension (通过跨肺压滴定的呼气末正压可促进腹内高压的急性呼吸窘迫综合征病人的氧合状态和呼吸力学)

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Comparison of 2 correction methods for absolute values of esophageal pressure in subjects with acute hypoxemic respiratory failure, mechanically ventilated in the ICU (比较 2 个在 ICU 进行机械通气后纠正急性低氧性呼吸衰竭受试者食道压绝对值的方法)

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Esophageal pressures in acute lung injury: do they represent artifact or useful information about transpulmonary pressure, chest wall mechanics, and lung stress? (急性肺损伤病人的食道压是跨肺压、胸壁力学及肺应力的假象还是有用信息?)

Loring SH, O'Donnell CR, Behazin N, Malhotra A, Sarge T, Ritz R, Novack V, Talmor D
J Appl Physiol. 2010 Mar;108(3):515-22
PMID 20019160, <http://www.ncbi.nlm.nih.gov/pubmed/20019160>

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Influence of lung and chest wall compliances on transmission of airway pressure to the pleural space in critically ill patients (肺及胸壁顺应性对气道压传送至危重症病人胸膜腔的影响)

Jardin F, Genevray B, Brun-Ney D, Bourdarias JP
Chest. 1985 Nov;88(5):653-8
PMID 3902386, <http://www.ncbi.nlm.nih.gov/pubmed/3902386>

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A simple method for assessing the validity of the esophageal balloon technique (一个评估食道气囊技术有效性的简单方法)

Baydur A, Behrakis PK, Zin WA, Jaeger M, Milic-Emili J. Am Rev Respir Dis 1982 Nov;126(5):788-91
PMID 7149443, <http://www.ncbi.nlm.nih.gov/pubmed/7149443>

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Pulmonary, chest wall, and lung-thorax elastances in acute respiratory failure (急性呼吸衰竭病人的肺、胸壁及胸肺弹性)

Katz JA, Zinn SE, Ozanne GM, Fairley HB
Chest. 1981 Sep;80(3):304-11
PMID 6944170, <http://www.ncbi.nlm.nih.gov/pubmed/6944170>

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Topography of esophageal pressure as a function of posture in man (人的体位与食道压分布形态的函数关系)

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Value and limitations of transpulmonary pressure calculations during intra-abdominal hypertension (在腹内高压期间跨肺压计算的数值和限制)

Cortes-Puentes GA, Gard KE, Adams AB, Faltsek KA, Anderson CP, Dries DJ, Marini JJ

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Pleural pressure and optimal positive end-expiratory pressure based on esophageal pressure versus chest wall elastance: incompatible results (基于食道压与胸壁弹性的胸腔压力和最佳呼气末正压: 不一致的结果)

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Volume-related and volume-independent effects of posture on esophageal and transpulmonary pressures in healthy subjects (体位对健康受试者食道压及跨肺压的容量相关性及容量独立性影响)

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Recruitment and derecruitment during acute respiratory failure: an experimental study (急性呼吸衰竭期间的肺复张和塌陷: 一项试验研究)

Pelosi P, Goldner M, McKibben A, Adams A, Eccher G, Caironi P, Losappio S, Gattinoni L, Marini JJ

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Validation of esophageal pressure occlusion test after paralysis (对麻痹后食道压阻断试验的验证)

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Lung mechanics in sitting and horizontal body positions (在坐位及平躺位的肺力学)

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Additional files (其它文献)

Esophageal and transpulmonary pressure in the clinical setting: meaning, usefulness and perspectives (临床环境下的食道压和跨肺压: 含义、有效性和前景)

Mauri T, Yoshida T, Bellani G, Goligher EC, Carteaux G, Rittayamai N, Mojoli F, Chiumello D, Piquilloud L, Grasso S, Jubran A, Laghi F, Magder S, Pesenti A, Loring S, Gattinoni L, Talmor D, Blanch L, Amato M, Chen L, Brochard L, Mancebo J; 平台压工作组 (PLUG - 欧洲危重病医学会急性呼吸分会)

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The application of esophageal pressure measurement in patients with respiratory failure (食道压测量对呼吸衰竭病人的应用)

Akoumianaki E(1), Maggiore SM, Valenza F, Bellani G, Jubran A, Loring SH, Pelosi P, Talmor D, Grasso S, Chiumello D, Guérin C, Patroniti N, Ranieri VM, Gattinoni L, Nava S, Terragni PP, Pesenti A, Tobin M, Mancebo J, Brochard L.

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Transpulmonary pressure: the importance of precise definitions and limiting assumptions (跨肺压: 准确定义和限制性假设的重要性)

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The promises and problems of transpulmonary pressure measurements in acute respiratory distress syndrome (急性呼吸窘迫综合征病人跨肺压测量的预示作用和问题)

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Measurement of esophageal pressure at bedside: pros and cons (食道压临床测量的利弊)

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Two steps forward in bedside monitoring of lung mechanics: transpulmonary pressure and lung volume (在肺力学临床监测方面向前迈进了两步: 跨肺压和肺容量)

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Goal-directed mechanical ventilation: are we aiming at the right goals? A proposal for an alternative approach aiming at optimal lung compliance, guided by esophageal pressure in acute respiratory failure (目标导向的机械通气: 我们是否瞄准正确的目标? 关于通过食道压引导瞄准急性呼吸衰竭病人最佳肺顺应性的可选择方法的一项建议)

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Esophageal pressure: benefit and limitations (食道压: 优点和不足)

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Targeting transpulmonary pressure to prevent ventilator induced lung injury (确定跨肺压目标以预防呼吸机诱发性肺损伤)

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
Esophageal and gastric pressure measurements (食道压和胃压测量)

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Respiratory mechanics in mechanically ventilated patients (机械通气病人的呼吸力学)

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