Instructions for use	ProxyCo CO2 sensor set.....	6
Gebrauchsanweisung	ProxyCo CO2-Sensor-Set.....	22
Instructions d'utilisation	Ensemble de capteur de CO2 ProxyCo.....	38
Instrucciones de uso	Conjunto del sensor de CO2 ProxyCo.....	54
Istruzioni per l'uso	Set sensore di CO2 ProxyCo.....	69
Instruções de uso	Conjunto de sensores CO2 ProxyCo.....	84
使用説明書	ProxyCo CO2センサセット.....	99
使用说明	ProxyCo CO2 传感器套装.....	113
Инструкции по эксплуатации	Комплект CO2-датчика ProxyCo.....	125
Kullanım Talimatları	ProxyCo CO2 Sensör seti.....	141

REF 10169182, 10169159, 10083174, 10100344

CE 0197


HAMILTON
MEDICAL

Instructions for Use





ProxyCo CO2 Sensor Set

2026-02-26

10117068/00

- 1 ProxyCo CO2 sensor | ProxyCo CO2-Sensor | Capteur de CO2 ProxyCo | Sensor de CO2 ProxyCo | Sensore di CO2 ProxyCo | Sensor CO2 ProxyCo | ProxyCo CO2 センサ | ProxyCo CO2 传感器 | CO2-датчик ProxyCo | ProxyCo CO2 sensörü
- 2 ProxyCo CO2 sensor connectors (straight or angled) | Stecker für den ProxyCo CO2-Sensor (gerade oder abgewinkelt) | Connecteurs de capteur de CO2 ProxyCo (droits ou coudés) | Conectores del sensor de CO2 ProxyCo (rectos o en ángulo) | Connettori (dritto o angolato) per sensore di CO2 ProxyCo | Conectores do sensor CO2 ProxyCo (retos ou angulados) | ProxyCo CO2 センサコネクタ (ストレートまたは角度 90°) | ProxyCo CO2 传感器接头 (直套管或直角接头) | Прямые или угловые разъемы для подключения CO2-датчика ProxyCo | ProxyCo CO2 sensör konnektörleri (düz veya açılı)
- 3 ProxyCo CO2 airway adapter, adult/pediatric | ProxyCo CO2-Atemwegsadapter, Erwachsene/pädiatrische Patienten | Adaptateur de CO2 ProxyCo, adulte/enfant | Adaptador de vía aérea de CO2 ProxyCo, pacientes adultos/pediátricos | Adattatore vie aeree per CO2 ProxyCo, pazienti adulti/pediatrici | Adaptador de vías aéreas CO2 ProxyCo, adulto/pediatrico | ProxyCo CO2 エアウェイアダプタ、成人 / 小児用 | ProxyCo CO2 气道接口，成人 / 儿童 | Адаптер воздуховода для CO2-датчика ProxyCo для взрослых/детей | ProxyCo CO2 hava yolu adaptörü, yetişkin/pediyatrik
- 4 ProxyCo CO2 airway adapter, neonatal | ProxyCo CO2-Atemwegsadapter, Neonaten | Adaptateur de CO2 ProxyCo, néonatal | Adaptador de vía aérea de CO2 ProxyCo, pacientes neonatos | Adattatore vie aeree per CO2 ProxyCo, pazienti neonatali | Adaptador de vías aéreas CO2 ProxyCo, neonatal | ProxyCo CO2 エアウェイアダプタ、新生児用 | ProxyCo CO2 气道接口，新生児 | Адаптер воздуховода для CO2-датчика ProxyCo для младенцев | ProxyCo CO2 hava yolu adaptörü, yenidoğan
- 5 Sensor holder | Sensorhalter | Support de capteur | Soporte del sensor | Supporto per sensore | Suporte do sensor | センサホルダー | 传感器支架 | Держатель датчика | Sensör tutucu

10169182 ProxyCo CO2 sensor set

10169159 ProxyCo CO2 sensor set, angled connector

10083174 ProxyCo CO2 airway adapter, adult/pediatric

10100344 ProxyCo CO2 airway adapter, neonatal

Conventions	7	Preventive maintenance	16
Sensor set overview	7	Technical specifications	16
Intended use	7	Embedded firmware	16
Potential adverse effects on stated performance	8	Compatibility	16
Clinical benefits	8	Physical characteristics	17
Residual Risks	9	Environmental and storage conditions	18
Safety information	9	Technical performance data	18
Setting up and connecting the sensor and airway adapter	12	Technical description	20
Connecting the sensor and airway adapter to the		Disposal	20
breathing circuit.	12	Glossary and abbreviations	21
Performing a zero calibration of the airway adapter	13		
Disconnecting the sensor from the ventilator	14		
Storing the sensor when not in use	14		
Monitoring	14		
Alarms	14		
Cleaning and disinfection	14		
Approved cleaning/disinfection agents	15		
ProxyCo CO2 sensor, sensor holder, and cable cleaning			
and disinfection	15		

1 Conventions

WARNING

A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient, or in damage to the medical device or other property.

NOTICE

A NOTICE provides additional information intended to avoid inconvenience during operation.

2 Sensor set overview

The ProxyCo CO₂ sensor is a mainstream sensor that uses infrared (IR) absorption spectrometry to continuously measure CO₂ in the breathing gas, and to use this data to calculate pCO₂ and etCO₂.

The ProxyCo CO₂ sensor is connected to an airway adapter and the ventilator. After performing a zero calibration, the system is ready for use.

The CO₂ measurement is based on the absorption of a specific wavelength of nondispersive infrared (NDIR) light by CO₂. The sensor's response to a known concentration stored in the sensor's memory during factory calibration allows it to measure the CO₂ concentration in the airway adapter.

2.1 Intended use

Category / Part number	Description
Intended purpose / use	
10169182 10169159 10083174 10100344	The ProxyCo CO ₂ sensor and ProxyCo CO ₂ airway adapters, are intended to provide carbon dioxide monitoring to a host monitoring system.
Intended patient / target group	
10169182 10169159	Adult, pediatric, and neonatal patients
10083174	Adult and pediatric patients
10100344	Neonatal patients
Intended user / user group	
10169182 10169159 10083174 10100344	Qualified, trained healthcare professional under the direction of a physician.
Intended use environment	
10169182 10169159 10083174 10100344	The ProxyCo CO ₂ sensor and ProxyCo CO ₂ airway adapters are intended to be used only together with a compatible mechanical ventilator (with the corresponding software option installed); therefore, they will be used in the same environment as the ventilator.

Category / Part number	Description
Indications	
10169182	Continuous CO2 monitoring allowing time based and volumetric capnography assessment during artificial ventilation with a compatible mechanical ventilator (with corresponding software option installed).
10169159	
10083174	
10100344	
Contraindications	
10169182	None.
10169159	
10083174	
10100344	
Limitations	
10169182	The ProxyCo CO2 sensor should not be used in the presence of flammable anesthetics. Patients with mixed ventilation, perfusion, or metabolism problems are difficult to interpret and might require additional measurements.
10169159	

2.2 Potential adverse effects on stated performance

Performance of the ProxyCo CO2 sensory system may be adversely affected by the following:

- Gas sample condensate
- Leaks of sampled gas may result in lower pCO2 readings
- Environmental vibrations may increase the noise level of the ProxyCo CO2 sensor pCO2 reading

Note that neither humidity (noncondensing) nor cyclical pressures have any effect on the stated accuracy of the device.

2.3 Clinical benefits

Continuous CO2 monitoring during invasive mechanical ventilation provides several critical benefits that significantly improve patient care and outcomes:

- Early detection of complications
- Reduced invasive procedures
- Enhanced accuracy compared to other methods
- Support for optimal ventilator management
- Detection of hidden problems due to continuous assessment

2.4 Residual Risks

This product has been developed in accordance with Regulation (EU) 2017/745 and has undergone a detailed evaluation concerning its safety and effectiveness. All known risks have been mitigated with appropriate measures in accordance with ISO 14791.

It is recommended to periodically use blood gas analysis (BGA) to verify plausibility of the measured values and support clinical decision-making.

3 Safety information

Review all safety information in this section before setting up the ProxyCo CO2 sensor and airway adapter.

If you have questions about any of the information in this guide, contact your Hamilton Medical technical representative.

WARNING

- Read your ventilator's *Operator's Manual* and the *CO2 monitoring Instructions for use* before using the ProxyCo CO2 sensor and airway adapters.
- **MR UNSAFE.** Keep the ProxyCo CO2 sensor set away from magnetic resonance imaging (MRI) equipment. The device poses unacceptable risks to the patient, medical staff, and other persons within the MRI environment.
- The ProxyCo CO2 sensor set must *not* be used in an oxygen-rich environment, with flammable anesthetics, or in insufficiently ventilated areas. **Danger of fire.**
- Use only Hamilton Medical-approved parts and accessories that are specified in this document and in the product e-catalog. Modifications to any components of the ProxyCo CO2 sensor and airway adapter are *not* permitted.
- Use of accessories, transducers, and cables other than those specified by Hamilton Medical can result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment, and may result in improper operation.
- Visually inspect the ProxyCo CO2 sensor and airway adapter for damage or contamination before use. Do *not* use any components that are damaged or fail preoperational checks—replace them immediately.
- If an airway adapter is dropped or has been in contact with contaminated surfaces, it must be discarded; use a new airway adapter.
- Handle used products as contaminated, and dispose of all parts removed from the device according to your institution's protocols. Comply with all local, state, and federal regulations with respect to environmental protection, especially when disposing of the electronic parts.
- Do *not* reprocess, reuse, or resterilize ProxyCo CO2 airway adapters. Reprocessing may compromise the integrity and performance, and pose a risk of cross-contamination or infection. Hamilton Medical does not assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user.

- Follow the cleaning and disinfecting procedures for the ProxyCo CO2 sensor as described in this guide, and the cleaning agent manufacturer's *Instructions for Use*.
For cleaning and disinfecting the components, only use the approved cleaning and disinfection agents listed in Section 9.1.
- Do *not* sterilize the ProxyCo CO2 sensor.
- Portable radio frequency (RF) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm to any part of the breathing gas sensor, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Correct function of the device may be impaired by the operation of high-frequency surgical equipment, microwaves, shortwaves, or strong magnetic fields in close proximity.

CAUTION

- *The ProxyCo CO2 sensor is NOT protected against reanimation with a defibrillator. Disconnect the CO2 sensor before using a defibrillator on the patient.*
- *Adding attachments or other parts/assemblies to the device not listed in this guide can adversely affect device performance and may lead to patient injury.*
- *Always connect all components securely and check for leaks in accordance with standard clinical procedures.*

- *The ProxyCo CO2 sensor is NOT patient isolated. Use of the sensor does not require direct patient contact. If isolation is desired or required, it is the responsibility of the ventilator to provide the necessary isolation.*
- *It is the responsibility of the ventilator to provide appropriate current limiting on the ProxyCo CO2 sensor power supply rails.*
- *Keep all fluids away from any electrical connections.*
- *Do NOT place the ProxyCo CO2 sensor directly on the patient's skin; ensure a cloth barrier is between the sensor and skin at all times. The sensor can reach temperatures up to and above 43°C (109°F) depending on ambient temperature.*
- *Always use the correct airway adapter for the patient group.*
- *Position airway adapters with windows in a vertical position to help prevent patient secretions from pooling on the windows. If pooling occurs, replace the airway adapter.*
- *Positioning of cables:*
 - *Do NOT position the cables in any manner that may cause patient entanglement or strangulation.*
 - *Support the cables to avoid stress on the ET tube.*
 - *Do NOT apply excessive tension to any cables.*
- *To ensure valid measurements:*
 - *Always perform a zero calibration when changing the airway adapter.*
 - *Always perform a zero calibration with the ProxyCo CO2 sensor connected to the airway adapter.*

- Calibrate the flow sensor after connecting the airway adapter to the breathing circuit.
- If there are any questions about the CO₂ values from the sensor, review the CO₂ values from a blood gas analysis (BGA), if applicable.
- Nebulization can affect CO₂ measurements.
- Check for breath sounds over the stomach if the pCO₂ curve shows (very low) values that do not correlate with the patient's breaths after tracheal intubation.
- Do NOT use with anesthetic agents, as this may result in incorrect CO₂ readings. The influence of isoflurane, sevoflurane, and norflurane on measurements is provided in Section 11.5.1.
- Periodically check the sensor, airway adapter, and tubing for excessive moisture or secretion build-up. Replace or clean the sensor/tubing if needed, and, if required, replace the airway adapter. Excessive moisture can affect measurements.
- To reduce the risk of fire and electric shock, before cleaning and disinfecting components, always disconnect the ProxyCo CO₂ sensor from electrical power and the patient, and ensure it has cooled down to room temperature.
- If the sensor is dropped, do the following, in order:
 - Clean the sensor (Section 9.2).
 - Visually check the sensor for damage.
 - Perform a zero calibration to ensure it is still functioning properly (Section 4.2).

NOTICE

- Any incident with the device leading to serious patient injury, death, or a potential threat to public health must be reported to the manufacturer and the relevant authorities.
- We recommend additional independent monitoring devices, including pulse oximeters measuring SpO₂, to be used during therapy. The operator assumes full responsibility for proper therapy and patient safety in all situations.
- Use of an oxygen concentrator may result in an increased baseline (up to a reading of 0.2% CO₂ during inspiration).
- The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A).

If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment. The ProxyCo CO₂ sensor complies with the IEC 60601-1-2 EMC (Electromagnetic Compatibility) Collateral Standard.
- (USA only) Federal law restricts this device to sale by or on the order of a physician.

4 Setting up and connecting the sensor and airway adapter

Before proceeding, review the safety information in Section 3.

Setting up the CO₂ sensor for use comprises the following steps:

1. Connecting the sensor and airway adapter to the breathing circuit (Section 4.1)
2. Performing a zero calibration of the airway adapter (Section 4.2)

For information on enabling CO₂ monitoring on the ventilator, see the *CO₂ monitoring Instructions for use* for your ventilator.

4.1 Connecting the sensor and airway adapter to the breathing circuit

To set up CO₂ monitoring

See *Figures A, B, and E on pages 4 and 5.*

1. Connect the ProxyCo CO₂ sensor connector (**A.2**) to the CO₂ port (**E**) on the ventilator.
2. Attach the ProxyCo CO₂ sensor to the appropriate airway adapter for the patient (**A.3, A.4**).
Align the narrow end of the airway adapter with the narrow end of the sensor, and press the components together until they click into place.
3. Perform a zero calibration of the airway adapter (Section 4.2).

4. Connect the airway adapter to the breathing circuit proximal to the patient (**B**), with the airway adapter windows in a vertical position.

Do *not* place the airway adapter between the ET tube and any connected adapter, as this may allow patient secretions to accumulate in the adapter.

5. Ensure the sensor cable faces *away* from the patient, and secure the cable safely out of the way.
6. Perform the rest of the preoperational checks on the ventilator (flow sensor calibration and Leak test). For details, see your ventilator *Operator's Manual*.
7. Start ventilation.
8. Verify the quality of the connection, as described next. Note that the CO₂ sensor has a 30 second warm-up time.

To verify the quality of the connection

- Check the capnogram (CO₂ waveform) on the ventilator display.

If the values are not as expected, check the sensor/airway adapter connections, and if needed, re-calibrate the airway adapter (Section 4.2).

4.2 Performing a zero calibration of the airway adapter

Perform a zero calibration:

- When changing airway adapters
- When a calibration-related alarm is generated on the ventilator

To perform zero calibration

See Figure D on page 5.

1. Disconnect the airway adapter from the breathing circuit, if necessary.
2. Ensure the airway adapter is attached to the CO₂ sensor (D).
3. Ensure the sensor has been connected to the ventilator for at least two (2) minutes.
4. Access the CO₂ sensor zero calibration function on your ventilator and start the calibration.

For details, see the *CO₂ monitoring Instructions for use* for your ventilator.

5. During calibration:
 - Keep the sensor away from all sources of CO₂, including the patient's and your own exhaled breath, as well as the ventilator exhaust port.
 - Do *not* move the sensor.
6. After the ventilator indicates a successful zero calibration, reconnect the airway adapter to the breathing circuit as described in Section 4.1.

In case of zero calibration failure

Perform the following checks, repeating zero calibration after each listed step until zero calibration is successful.

1. Check the airway adapter, and replace if necessary.
2. Check the sensor windows for potential pollution.
3. Ensure there is *no* source of CO₂ near the airway adapter.
4. Ensure there is no movement or vibration during zero calibration.
5. Connect a new airway adapter.
6. Connect a new ProxyCo CO₂ sensor.

If the problem persists, have the ventilator serviced.

5 Disconnecting the sensor from the ventilator

To disconnect the sensor

1. Disable the sensor on the ventilator.
See the *CO2 monitoring Instructions for use* for your ventilator for details.
2. Gently pull the connector straight back and disengage from the CO2 port on the ventilator.

6 Storing the sensor when not in use

Ensure the sensor is cleaned and disinfected as described in Section 9.2.

- Store the sensor according to the sensor lifetime (Section 11.5) on the sensor holder when it is disconnected and not in use. See Figure C on page 5.

7 Monitoring

Sensor measurements are displayed on the ventilator using:

- Capnograms, both time based and volume based (single-breath CO2)
- CO2-related monitoring parameters

See the *CO2 monitoring Instructions for use* for your ventilator for a full list of CO2-related monitoring parameters.

8 Alarms

Operator-adjustable and nonadjustable alarms, together with a visual alarm indicator, help ensure your patient's safety. These alarms notify you of conditions that require your attention. They are categorized as high, medium, or low priority.

For details, see the *CO2 monitoring Instructions for use* for your ventilator.

To verify CO2-related alarms are operational

- Check for the low-priority sensor warm-up alarm during the first 30 seconds after connecting the sensor. (The sensor has a 30-second warm-up time.)
If the alarm does not appear, disconnect and reconnect the sensor.

9 Cleaning and disinfection

The ProxyCo CO2 sensor set must be regularly cleaned and disinfected using the cleaning methods and solutions specific to the individual component.

Cleaning/disinfection agents are tested as a specific composition of different chemicals by our certified laboratory to validate their materials compatibility and cleaning/disinfection efficacy. See Section 9.1.

Other surface cleaning/disinfection agents are used at your own risk. Hamilton Medical is *not* responsible for damages, incidents, or adverse events due to the use of any cleaning/disinfection agents that are *not* listed in Section 9.1.

9.1 Approved cleaning/disinfection agents

NOTICE

Use of cleaning and disinfectant wipes may cause marginal color change in some materials. Marginal color change on its own is *not* an indicator of product malfunction.

Cleaning/disinfection wipe	Contact time (min)
EPA-registered cleaning/disinfection wipe	
CaviWipes 2.0	2
Approved cleaning/disinfection wipe	
Incidin OxyWipes S	2
CaviWipes 2.0	2
mikrozid® universal wipes green line	1
Environmental Cross V-Lock	3

9.2 ProxyCo CO2 sensor, sensor holder, and cable cleaning and disinfection

To surface clean a component

Surface clean each component before and after each patient use, or according to your institution's protocols.

1. Unfold a clean disposable cleaning and disinfectant wipe (from Section 9.1), and wipe off visible soil from the component surface.
2. Wet all surfaces, working from top to bottom, from clean to dirtier areas, in an S-shaped fashion.
Should your wipe dry out or become soiled, replace it with a new wipe to complete the area.
Ensure the surface is visually free of soil.
Do *not* overuse the wipe.
3. Dispose of the wipe.

To surface disinfect a component

After cleaning the surface as instructed previously, do the following.

1. Wipe cleaned surfaces again, until visibly wet, with a cleaning and disinfectant wipe listed in Section 9.1.
2. Dispose of the wipe.
3. Wait for the contact time of your chosen wipe.
4. Wait until the surfaces are dry.

5. Check the surfaces for visible damage.

If there is any damage, contact your Hamilton Medical technical representative.

6. After use of the wipes, ensure their packaging lid is firmly sealed.

9.2.1 Single-use (SU) ProxyCo CO2 airway adapters

Do *not* clean, soak, or wash the single-use airway adapter.

Replace the product in accordance with your hospital infection control procedures and/or depending on the patient's secretions.

A single-use airway adapter can be used continuously on a single patient for a maximum of seven (7) days.

10 Preventive maintenance

The ProxyCoCO2 sensor is end-user maintenance free.

11 Technical specifications

The following sections provide the technical specifications for the components.

11.1 Embedded firmware

The ProxyCo CO2 sensor contains electronic programmable elements. The sensor does not include user-configurable software and does not communicate directly with IT networks.

All data processing and all user interactions are performed exclusively via the ventilator to which the sensor is connected.

11.2 Compatibility

The ProxyCo CO2 sensor set is compatible with, and has been verified for use with, the HAMILTON-C6 ventilator. The operator is responsible for ensuring compatibility of the ventilator with the ProxyCo CO2 sensor set.

11.3 Physical characteristics

Materials used

Component	Part number	Materials
ProxyCo CO2 sensor	10169182 10169159	PA, TPE
ProxyCo CO2 airway adapters	10083174 10100344	PC
Sensor holder	10106089	PC

Dimensions

Dimensions	Specifications
ProxyCo CO2 sensor PN 10169182, 10169159	
Weight	< 100 g
Dimensions, sensor head (L x H x W)	24 x 32 x 48 mm
Cable length	3.2 m
ProxyCo CO2 airway adapter, adult/pediatric PN 10083174	
Weight	7.7 g
Dimensions (L x H x W)	57 x 23 x 23 mm
Conical connectors	<i>Proximal:</i> ID15/OD22 <i>Distal:</i> ID15

Dimensions	Specifications
ProxyCo CO2 airway adapter, neonatal PN 10100344	
Weight	7.2 g
Dimensions (L x H x W)	52 x 22 x 23 mm
Conical connectors	<i>Proximal:</i> ID15 <i>Distal:</i> OD15
Sensor holder PN 10106089	
Dimensions (L x H x W)	133 x 40 x 37 mm Fits on standard rail

11.4 Environmental and storage conditions

Environment		Specifications
ProxyCo CO2 sensor and airway adapters PN 10169182, 10169159, 10083174, 10100344		
Temperature	Operating (continuous)	-10°C to 40°C (14°F to 104°F)
	Operating (transient)	Maximum 20 minutes at -20°C to 50°C (-4°F to 122°F)
	Storage	15°C to 30°C (59°F to 86°F)
	Shipping (max. one (1) week)	-20°C to 60°C (-4°F to 140°F)
Relative humidity	Operating	15% to 95%, noncondensing
	Storage and shipping	10% to 95%, noncondensing
Atmospheric pressure		376 to 1100 hPa
Ingress protection of sensor		IP57

11.5 Technical performance data

The data in this table applies to the ProxyCo CO2 sensor set unless otherwise noted.

Description	Specifications
Sensor type	Mainstream CO2 sensor
Principle of operation	Nondispersive infrared (NDIR) technology
Type of operation	Continuous
Classification, according to IEC 60601-1	Transportable Class II ME equipment
Insulation	Applied part of type BF
CO2 measurement range	0 to 220 hPa (0 to 165 mmHg) partial pressure
CO2 accuracy (in volume percent)	± (0.43% + 8% of gas level)
CO2 resolution	0.1 hPa (0.07 mmHg)
Initialization time (time from turning on device to operating performance)	Full specifications within two (2) minutes; capnogram displayed within 30 seconds
Response time (sensor)	< 80 ms
Sampling rate	PCO2: 100 Hz etCO2: 1 Hz
Flow rate	0 to 260 l/min
Respiration rate range	0 to 150 b/min

Description	Specifications
Dead space, airway adapter Adult/pediatric	≤ 5.0 ml
Dead space, airway adapter Neonatal with flow sensor (PN 155500)	3.8 ml
Compensation (using values provided by ventilator)	• Ambient pressure compensation • FiO2 compensation • Heliox compensation
Maximum housing temperature	43°C (109°F) with ambient temperature up to 35°C (95°F)
Data output	Measurement data and sensor status using a proprietary protocol
Data transmission protocol	Proprietary protocol over a serial RS-232 connection. The nominal data rate is 38,400 baud.
Voltage requirements	+ 5.00 VDC ± 5%
Power rating, total consumption	0.7 W typical, steady state
Power loss	After a power loss, the sensor restarts operation with the previous calibration values. The warm-up period is reset.
Usage time (airway adapter) PN 10083174, 10100344	Seven (7) days maximum continuous use on a single patient

Description	Specifications
Lifetime (sensor) PN 10169182, 10169159	5 years (44,000 hours)
Shelf life (airway adapter) PN 10083174, 10100344	5 years

11.5.1 Influence of interfering gases

Influence on the measurement of 5% exhaled CO2 with a listed interfering gas present.

Interfering gas and gas level	Relative influence on measured value
Isoflurane 5%	< 15% increase
Sevoflurane 5%	< 2% decrease
Norflurane 5% (as MDI propellant only)	< 2% decrease
Ethanol 0.1%	< 1% decrease
Isopropanol 0.1%	< 1% decrease
Acetone 0.1%	< 1% decrease
Methane 1%	< 1% decrease

11.6 Technical description

The ProxyCo CO₂ sensor is a mainstream sensor using infrared (IR) absorption spectrometry to continuously measure CO₂ in the breathing gas, and to use this data to calculate pCO₂ and etCO₂¹.

The ProxyCo CO₂ sensor is connected to the ventilator and airway adapter, and after performing a zero calibration, the system is ready for use.

The CO₂ measurement is based on the absorption of a specific wavelength of nondispersive infrared (NDIR) light by CO₂. The sensor's response to a known concentration stored in the sensor's memory during factory calibration allows it to measure the CO₂ concentration in the airway adapter.

The etCO₂ value is calculated from the maximum of a moving average of the partial pressure of CO₂ (pCO₂) over the exhalation phase.

The exhalation phase is determined as follows:

The *beginning* of the exhalation phase is detected when the pCO₂ measurement is above a defined threshold. The *end* of the phase is identified where there is a significant drop in the pCO₂ measurement, indicating the change to inhalation.

The calculated etCO₂ value and the exhalation phase interval used are checked to ensure the values are plausible given the specifications the sensor is using.

To qualify as a valid measurement, the exhalation phase length must be longer than a minimum interval allowed by the supported respiration rate.

For short expiration phases (< 250 ms), the etCO₂ value may be too low.

The etCO₂ value is updated every one (1) second.

The newest value in this time interval is transmitted to the host, regardless of the number of breaths. A quality index is generated on the Hamilton Medical ventilator for etCO₂¹, displayed next to the value in the Monitoring window and when using INTELLiVENT-ASV.

12 Disposal

Used components must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts. This includes adhering to Council Directive 2002/96/EC or WEEE Directive (Waste Electrical and Electronic Equipment).

1. Shown on Hamilton Medical ventilators in the PetCO₂ parameter.

13 Glossary and abbreviations

The table below lists definitions of terms and abbreviations used throughout this document.

Term/abbreviation	Definition
CO ₂	Carbon dioxide
etCO ₂	End tidal CO ₂
ET tube	Endotracheal tube
MRI	Magnetic resonance imaging
pCO ₂	The partial pressure of carbon dioxide
REF	Reference number; also referred to as a part number
SpO ₂	Measured oxygen saturation.
SU	Single use
VDC	Volts direct current

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Specifications are subject to change without notice. For all proprietary trademarks and third-party trademarks used by Hamilton Medical AG, see www.hamilton-medical.com/trademarks.
© 2026 Hamilton Medical AG. All rights reserved.
Hamilton Medical AG
medin Medical Innovations GmbH
Hamilton Medical, Inc.