




Nuflow[®] nasal cannula set

for oxygen therapy in pediatrics

- ✓ High degree of patient comfort
- ✓ Low noise level and system stagnation pressure
- ✓ Individual fixation, gentle to the skin
- ✓ Anatomically adapted shape of the prongs

The Nuflow nasal cannula is made out of soft and skin compatible silicone which helps to perfectly balance therapy and wearing comfort. Even the smallest cannula size provides flow rates up to 8 l/min.



Available cannula sizes

Size	small	medium	large	xlarge
Patient group¹	Premature	Newborns	Infants	Children
Cannula Prod. No*	10072354	10072355	10072356	10072357
Flow limits²	10 l/min	14 l/min	23 l/min	27 l/min
Outer diameter prong	2,3 mm	2,7 mm	3,2 mm	3,7 mm
Septum distance	3,7 mm	4,3 mm	4,8 mm	5,3 mm
Packing unit	10 pieces	10 pieces	10 pieces	10 pieces
Fixation patch size	small	small	large	large
Fixation patch Prod. No	10072358	10072358	10072359	10072359
Packing unit	10 pieces	10 pieces	10 pieces	10 pieces

1 - Depends on size of patient nose

2 - See user manual for details



* The scope of delivery includes:

- ✓ Nasal Cannula in selected size
- ✓ Matching fixation patch
- ✓ Instruction for use
- ✓ Circuit adapter

Additional fixation patches can be ordered separately.

Distributor:

Hamilton Medical, Inc.

4655 Aircenter Circle, Reno, NV 89502, USA

☎ +1 800 426 6331

info@hamiltonmedical.com

www.hamilton-medical.com

ELO202000103.00

Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. For all proprietary trademarks (®) and third-party trademarks (®) used by Hamilton Medical AG see www.hamilton-medical.com/trademarks. © 2020 Hamilton Medical AG. All rights reserved.