


Nuflow[®]

Nasal cannula set



Available sizes:

	1 Nasal cannula set small	2 Nasal cannula set medium	3 Nasal cannula set large	4 Nasal cannula set xlarge
REF	1320	1330	1340	1350
low limits ¹	10,0 l/min	14 l/min	23 l/min	27 l/min
Outer diameter prong	2,3 mm	2,7 mm	3,2 mm	3,7 mm
Septum distance	3,7 mm	4,3 mm	4,8 mm	5,3 mm
Patient group ²	Premature	Premature, Newborns	Infants	Children
Packing Unit	10 pieces	10 pieces	10 pieces	10 pieces
Fixation patches	1380 small	1380 small	1381 large	1381 large

Connectors

Tube connection: OD 10 mm or OD 15 mm available

1 - See user manual for details

2 - Depends on size of patient nose



Material

Cannula Silicone, Zylar, TPE, PP,
Fixation patch Hydrocolloid
- free of PVC, DEHP, natural rubber latex

Consumables

REF1380 Fixation patch, small (for use with Nuflow size REF1320 or REF1330)
REF1381 Fixation patch, large (for use with Nuflow size REF1340 or REF1350)

Combinable devices

medinSINDI®, medinCNO®mini, medin-NC3®, medinBLENDER
Ventilators with high-flow mode and overpressure protection

Regulations

Classification	EU IIa
Country of origin	Germany
Intended use	single patient use, non invasive, max. 7 days
Manufacturing surrounding	clean room
Disposal	The Nuflow Nasal Cannula and its packaging can be discarded with regular household waste.

Manufacturer:

medin Medical Innovations GmbH

Adam-Geisler-Straße 1

☎ +49 (0)8142 448 46 -0

info@medin-medical.com

www.medin-medical.com

REV01
09.10.2018

Specifications are subject to change without notice. Some features are options. Not all features are available in all markets.
For all proprietary trademarks (®) and third-party trademarks (®) used by medin Medical Innovations GmbH see
www.medin-medical.com/trademarks. © 2018 medin Medical Innovations GmbH. All rights reserved.

Nuflow

