


Transport and storage solutions for the HAMILTON-C2/C3

Docking trolley

- ✓ Developed to facilitate intrahospital transport
- ✓ Connects directly to the patient bed, allowing a single person to move both bed and ventilator simultaneously
- ✓ Designed to accommodate a HAMILTON-C2/C3 ventilator plus standard equipment:

- 1 Patient monitor holder
- 2 Infusion bag holder
- 3 Bed docking arms
- 4 Mount for a septic fluid aspirator
- 5 Holder for two gas cylinders
- 6 Mount for a HAMILTON-H900 humidifier
- 7 Basket
- 8 Handle
- 9 Holder for infusion pumps



Trolley for the HAMILTON-C2/C3

- ✓ For intrahospital transport
- ✓ Mount for a HAMILTON-H900 humidifier
- ✓ Holder for two gas cylinders
- ✓ Mount for a standard rail
- ✓ 4 high-quality independent wheels with brakes
- ✓ Tubing support arm with quick-lock

Universal bed mount kit

- ✓ Used to mount the HAMILTON-C2/C3 to a patient bed
- ✓ Fits onto round and edged rails of various dimensions



Bed mount with cylinder holder

- ✓ Used to mount the HAMILTON-C2/C3 to a patient bed
- ✓ Space for one oxygen cylinder

Product

Part number

Docking trolley

160092

Trolley for the HAMILTON-C2/C3

160170

Universal bed mount kit

160148

Bed mount with cylinder holder

160144

Manufacturer:

Hamilton Medical AG

Via Crusch 8, 7402 Bonaduz, Switzerland

☎ +41 (0)58 610 10 20

info@hamilton-medical.net

www.hamilton-medical.com