



NIV masks and interchangeable elbow portfolio

## BiTrac NIV masks

The BiTrac NIV full face and nasal masks are single patient multiple use masks that provide compatible ventilation for spontaneously breathing patients requiring noninvasive ventilation (NIV).

### Benefits for you and your patients



The mask provides a double-lip cushion.



The OmniClip is an adjustable piece on the forehead pad.



The headgear has adjustable straps with clips.



# BiTrac Select NIV masks interchangeable elbow portfolio

The new interchangeable elbow concept for full face and MaxShield masks



Interchangeable elbow	PN	QTY
Standard - Used with a coaxial or dual limb circuit on a Hamilton Medical ventilator.	10161785	10
Standard Nebulizer - Used to administer nebulized medication on a Hamilton ventilator with a dual limb circuit. Compatible with Aerogen® nebulizers.	10161781	10
Standard Bronchoscopy - Used for bronchoscopy on a Hamilton ventilator.	10161783	10
Anti-Asphyxia Nebulizer* - Used to administer nebulized medication on a BiPAP device. Exhalation port in the circuit. Compatible with Aerogen® nebulizers.	10161780	10
Anti-Asphyxia with leak* - Used for BiPAP devices without an exhalation port in the circuit. Compatible with Aerogen nebulizers.	10161779	10
Anti-Asphyxia* - Used for BiPAP devices with an exhalation port in the circuit.	10161784	10
Anti-Asphyxia Bronchoscopy* - Used for bronchoscopy on a BiPAP device.	10161782	10

\*Not to be used with ventilators

# BiTrac NIV masks

With standard elbow for all Hamilton Medical ventilators, single patient multiple use



## BiTrac® Select MaxShield

- Interchangeable standard elbow 22-mm ID
- Bonnet style headgear
- Recommended for patients 1 year or older and > 7 kg

Size	Part number	QTY
XXS	10161830	10
XS	10161831	10



## BiTrac NIV nasal mask

- Fixed standard elbow 22-mm ID
- OmniClip with silicone forehead pad and head strap
- Adult patients

Size	Part number	QTY
S	10161579	10
M	10161578	10
L	10161577	10
XL	10161576	10



## BiTrac® Select NIV full face mask

- Interchangeable standard elbow 22-mm ID
- OmniClip with silicone forehead pad and head strap
- Adult patients

Size	Part number	QTY
S	10161769	10
M	10161770	10
L	10161771	10
XL	10161772	10



## BiTrac® Select MaxShield

- Interchangeable standard elbow 22-mm ID and head strap
- Adult patients

Size	Part number	QTY
S	10161778	10
L	10161777	10
XL	10161776	10

# BiTrac NIV masks and circuit

With anti-asphyxia vented elbow for BiLEVEL and CPAP application, single patient multiple use



## BiTrac® Select MaxShield\*

- Interchangeable anti-asphyxia elbow 22-mm OD
- Bonnet style headgear
- Recommended for patients 1 year or older and > 7 kg

Size	Part number	QTY
XXS	10161832	10
XS	10161833	10



## BiTrac NIV full face mask\*

- Interchangeable anti-asphyxia elbow 22-mm ID
- OmniClip with silicone forehead pad and head strap
- Adult patients

Size	Part number	QTY
S	10161765	10
M	10161766	10
L	10161767	10
XL	10161768	10



## BiTrac® Select MaxShield\*

- Interchangeable anti-asphyxia elbow 22-mm ID and head strap
- Adult patients

Size	Part number	QTY
S	10161775	10
L	10161774	10
XL	10161773	10



## BiTrac NIV circuit

- 72" smooth ID tubing, swivel connector with exhalation port
- Available with or without attached filter

Part number	QTY
10160057 with filter	10
10161575 without filter	10

\*Not for use with Hamilton Medical ventilators.



For more information and instructions for use, visit our e-Catalog at:  
[www.hamilton-medical.com/e-catalog](http://www.hamilton-medical.com/e-catalog)



**HAMILTON**  
**MEDICAL**

Distributed by

Hamilton Medical AG

Via Crusch 8, 7402 Bonaduz, Switzerland

+41 (0) 58 610 10 20

[info@hamilton-medical.com](mailto:info@hamilton-medical.com)

[www.hamilton-medical.com](http://www.hamilton-medical.com)

**Pulmodyne**<sup>®</sup>

Manufactured by

Pulmodyne, Inc

2055 Executive Dr, 46241 Indianapolis, US

+1 317 900 4093

[info@pulmodyne.com](mailto:info@pulmodyne.com)

[www.pulmodyne.com](http://www.pulmodyne.com)

ELO20250403N.00

The products shown here are not available for purchase by the general public. The information provided here is intended for healthcare professionals only. Always read the labels and follow the product's instructions for use. Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. All images are used for illustrative purposes only and may not accurately represent the product or its use. For all proprietary trademarks (®) and third-party trademarks used by Hamilton Medical AG see [www.hamilton-medical.com/trademarks](http://www.hamilton-medical.com/trademarks). ©2025 Hamilton Medical AG. All rights reserved.

NIV masks