

NEW

Supercharge your NIV therapy. Optimize patient synchrony with the HAMILTON-C1



IntelliSync®+ brings continuous, real-time patient-ventilator synchronization to the HAMILTON-C1

The IntelliSync+ technology mimics the expert's eye by continuously analyzing waveform shapes at least a hundred times per second. This allows IntelliSync+ to detect patient efforts immediately, and to initiate inspiration and expiration in real-time.

For maximum flexibility, IntelliSync+ can be activated to automate either the inspiratory trigger or the expiratory trigger, or both.

Features of IntelliSync+:

- ✓ For invasive or noninvasive ventilation modes
- ✓ During inspiration and expiration
- ✓ Noninvasive method
- ✓ No additional hardware or accessories required



[Learn more](#)

HAMILTON
MEDICAL