



Hamilton Connect

Digital solutions for respiratory care

We aim to support our customers and provide value for them throughout the entire ventilation workflow. The new Hamilton Connect Module and the associated Hamilton Connect App are the first steps in realizing this vision. They offer a simple solution to facilitate and streamline your daily ventilation workflow.

Hamilton Connect Module

The Hamilton Connect Module integrated in the ventilator provides the HAMILTON-C1/T1 with the capability of secure wireless communication with patient data management systems, monitoring systems, and the Hamilton Connect App.

Hamilton Connect App

The Hamilton Connect App connects your mobile device to your Hamilton Medical ventilator and allows you to access selected ventilation data.

Download it for free from the Apple App Store and Google Play store

The Hamilton Connect App runs on Android[®] or iOS[®] mobile devices. It is available for free on the Apple App Store[®] and Google Play[®] store. Download the app by scanning this QR-code:



or visit: <http://qr.ham.md/get-HCA>

Hamilton Connect App in detail

Mobile respiratory data and tools

Built to meet the needs of modern healthcare providers, the Hamilton Connect App gives you access to respiratory data. It provides you with a solution designed to facilitate and streamline your ventilation workflow. All you need is a mobile device and a Hamilton Medical ventilator with the Hamilton Connect Module.

Demo mode – No ventilator required

Test the app in demo mode, even without a ventilator. Explore all the features and functionalities for as long as you like.

Respiratory data – Wherever you are

The Hamilton Connect App transforms your mobile device into a portable secondary ventilator display. It streams near real-time data from the ventilator to your mobile device, so you can review your patient's respiratory data even when you are not within reach of the ventilator.

Review data

Review active control settings directly from the app and freeze waveforms for closer review. Change the time scale and display a full-screen waveform or an overview of all parameters. Review active alarms, as well as available waveforms and parameter values.


Customizable views

Choose how the monitoring parameters and waveforms are displayed on your mobile device, and create the layout that best suits your needs. Respiratory data can be displayed independently of what is currently being shown on the ventilator.

We will continue to enhance the app with new features and functionalities to support your respiratory care.

For more information visit:

www.hamilton-medical.com/Hamilton-Connect



Manufacturer:

Hamilton Medical AG

Via Crusch 8, 7402 Bonaduz, Switzerland

+41 (0)58 610 10 20

info@hamilton-medical.com

www.hamilton-medical.com

ELO20210214N.01

The Hamilton Connect App is not intended to replace the real-time display of data on the ventilator. DO NOT USE the app to supplement or replace any part of the hospital's device monitoring. Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. For all proprietary trademarks (®) and third-party trademarks (®) used by Hamilton Medical AG see www.hamilton-medical.com/trademarks. © 2021 Hamilton Medical AG. All rights reserved.

Hamilton Connect