

NEW



New SW 3.1.1 for HAMILTON-C1/T1/MR1

The software update v3.1.1 brings new features and improvements to the HAMILTON-C1/T1/MR1 ventilators.



HAMILTON-C1/T1:

- ✓ NIV-only option
- ✓ O2 assist
- ✓ IntelliSync®+

HAMILTON-MR1:

- ✓ IntelliSync®+


For detailed information about the update, refer to the respective brochure, software release notes v3.1.1, and operator's manual for your device.



O2 assist. Optimized oxygen, effortlessly controlled

Discover the convenience and simplicity of O2 assist, available at the touch of a button on our HAMILTON-C1/T1 ventilators. Take control of oxygen management by setting individualized SpO2 targets for your patients. The system continuously adjusts the oxygen supply to meet your patient's real-time needs, potentially protecting against both hyperoxemia and hypoxemia.



www.hamilton-medical.com/O2-assist



NIV-only option. The complete solution for noninvasive ventilation

Less can also be more! The HAMILTON-C1/T1 with the NIV-only option offers just the modes you need for high-performance noninvasive ventilation support. It enables you to transition seamlessly between noninvasive modes and high flow oxygen therapy (optional) for pediatrics and adults, without changing the device or even the breathing circuit.

www.hamilton-medical.com/NIV-only



IntelliSync@+. Keeping an eye on patient-ventilator synchrony

The well-trained eye of a ventilation expert is capable of detecting asynchronies by analyzing the shapes of the flow and pressure waveforms. IntelliSync@+ offers valuable workflow support by managing triggering and cycling automatically. This feature can help healthcare professionals spend less time adjusting these settings and free them up to focus more on other important bedside tasks.

www.hamilton-medical.com/intellisync

