



Hamilton Medical NIV mask portfolio

Ventilator application for adult patients



Full face masks

BiTrac™ NIV full face mask, single-use
with standard elbow 22 mm ID, OmniClip™ with silicone forehead pad and head strap.

For HAMILTON-G5/C1/C2/C3/T1/MR1 ventilators

PN	Size	QTY
313-9009W	S	10
313-9010W	M	10
313-9011W	L	10
313-9012W	XL	10



Nasal masks

BiTrac NIV nasal mask, single-use
with standard elbow 22 mm ID, OmniClip with silicone forehead pad and head strap.

For HAMILTON-G5/C1/C2/C3/T1/MR1 ventilators

PN	Size	QTY
313-9013W	S	10
313-9014W	M	10
313-9015W	L	10
313-9016W	XL	10



MaxShield™ masks

BiTrac MaxShield, single-use
with fixed standard elbow 22 mm ID and head strap.

For HAMILTON-G5/C1/C2/C3/T1/MR1 ventilators

PN	Size	QTY
313-9055W	S	10
313-9053W	L	10
313-9056W	XL	10

BiLEVEL or CPAP application for adult patients



Full face masks

BiTrac NIV full face mask, single-use
with anti-asphyxia elbow 22 mm ID, OmniClip with silicone forehead pad and head strap.

PN	Size	QTY
313-9001W	S	10
313-9002W	M	10
313-9003W	L	10
313-9004W	XL	10



Nasal masks

BiTrac NIV nasal mask, single-use
with standard elbow 22 mm ID, OmniClip with silicone forehead pad and head strap.

PN	Size	QTY
313-9013W	S	10
313-9014W	M	10
313-9015W	L	10
313-9016W	XL	10



MaxShield masks

BiTrac MaxShield, single-use
with anti-asphyxia elbow 22 mm ID and head strap.

PN	Size	QTY
313-9054W	S	10
313-9052W	L	10
313-9057W	XL	10



NIV circuits

BiTrac NIV circuit, single-use
72" Smooth ID tubing, swivel connector with exhalation port, cap plug; includes optional 24" tubing with „T“ and 22 mm ID x 22 mm ID adapter.

Available with or without attached filter.

PN	Size	QTY
313-9074 without filter	72"	10
313-9075 with filter	72"	10

For more information and instruction for use visit our e-catalog on www.hamilton-medical.com/e-catalog

Contact:

Hamilton Medical, Inc.
4990 Energy Way, Reno, NV 89502, United States
☎ +1 800 426 6331
✉ info@hamiltonmedical.net
www.hamilton-medical.com