

# Preassembled breathing circuit sets for HAMILTON-C2/C3/C6

## Single use, for adult and pediatric patients

- Coaxial breathing circuit available in two lengths
- Including expandable expiratory tube
- Proximal flow sensor and expiratory valve set (optional)



## Your benefits at a glance

- ✓ Increased efficiency and safety as supplied assembled ready to use
- ✓ Compact single limb design allows for clearer patient space
- ✓ Reduced risk of cross infection as circuits are single use
- ✓ Greater time and cost efficiency as reprocessing is eliminated
- ✓ Optimal positioning thanks to expandable expiratory tube
- ✓ Coaxial design provides passive warming of inspiratory gasses by heat exchange with expiratory flow

## Ordering information

Part number	With flow sensor	With expiratory valve set	Length	Packaging unit
260257	✓	✓	180 cm / 71 in	20/box
260207	✓		180 cm / 71 in	20/box
260206			180 cm / 71 in	20/box
260240	✓		240 cm / 94 in	10/box
260239			240 cm / 94 in	10/box