

Preassembled breathing circuit sets for HAMILTON-C1/T1/MR1

Single use, for adult and pediatric patients

- Coaxial breathing circuit available in four different lengths
- Proximal flow sensor and elbow adapters
- Expiratory valve set (optional)



PN 260128

Your benefits at a glance

- ✓ Increased efficiency and safety as supplied assembled ready to use
- ✓ Compact single limb design allows for clearer patient space
- ✓ Reduced risk of cross infection as circuits are single use
- ✓ Greater time and cost efficiency as reprocessing is eliminated
- ✓ Coaxial design provides passive warming of inspiratory gasses by heat exchange with expiratory flow

Ordering information

Part number	With flow sensor	With expiratory valve set	Length	Packaging unit
260128	✓	✓	180 cm / 71 in	20/box
260087	✓		180 cm / 71 in	20/box
260127	✓	✓	240 cm / 94 in	10/box
260094	✓		240 cm / 94 in	10/box
260167	✓	✓	300 cm / 118 in	10/box
260145	✓		300 cm / 118 in	10/box
260168	✓	✓	480 cm / 189 in	8/box
260144	✓		480 cm / 189 in	8/box